


Experiments in College Chemistry

AN OPENTEXTBOOK

Version 2.8 - November 22, 2024



Science Department
B.M.C.C.
The City University of New York



**THIS BOOK
NEEDS TO BE
SPIRAL BOUNDED
+
PRINTED DOUBLESIDED**

Contents

1	Isomers	5
2	Conductivity of electrolyte solutions	15
3	Softening of hard water	27
4	Freezing point depression	37
5	The iodine clock	55
6	Le Châtelier's Principle	69
7	Equilibrium constant	79
8	Titration curve	95
9	Bleach Strength	109
10	Relative strength of some acids	123
11	Entropy	135
12	A Solubility Product Constant: Solubility product calculation	149
13	Electrochemistry: galvanic cells	163
14	Radioactivity: half-life of K-40	177

1 Isomers

EXPERIMENT

Geometric isomers

Goal

The goal of this experiment is to understand the concept of **isomerism**. During this practice, the student will identify two different molecules based on their properties. Later, a simple reaction will be carried out in an attempt to transform each molecule into its corresponding isomer. The student should be able to identify the new products based on the first observations (i.e. the properties of each molecule) and to validate the success of the transformation. The experiment will also familiarize the students with the separation technique known as suction filtration.

Materials

- | | |
|---|---|
| <input type="checkbox"/> 4 small test tubes, with 2 suitable test tube stoppers | <input type="checkbox"/> Suction flask |
| <input type="checkbox"/> 2 test tube holders | <input type="checkbox"/> Round filter paper |
| <input type="checkbox"/> Bunsen burner | <input type="checkbox"/> plastic funnel |
| <input type="checkbox"/> 250 mL Erlenmeyer | <input type="checkbox"/> stirring rod |
| <input type="checkbox"/> Ring stand, Iron ring, Wire gauze, Clamps and Büchner funnel | <input type="checkbox"/> Fumaric Acid, and Maleic Acid, |
| <input type="checkbox"/> Rubber tubing | <input type="checkbox"/> 12 M Hydrochloric Acid |
| | <input type="checkbox"/> Ice |

Background

Consider a set of Lego® blocks with different sizes and colors. We can combine them in numerous ways to construct many different structures. Similarly, if those blocks were atoms, we would obtain different molecules from the same set of atoms. Isomers are molecules with the same number and nature of atoms, and therefore the same molecular formula, but with different structures. Diverse structures arise from variations in the bonds between atoms and lead to new molecules, with different names and properties.

Geometric isomers are a subclass of isomers where not only the number and nature of the atoms are the same, but also the bonds. In this case, it is the spatial distribution or 3D structure that differs. This little change can lead to considerable differences in basic properties such as melting points and solubilities. In molecules with four ligands distributed around a central atom or a double bond, the isomerism can arise from the distribution of those four ligands around the molecule, leading to *cis*- and *trans*- isomers. This is the case in the proposed experiment.





The process of **isomerization** is the transformation from one molecule into its corresponding isomer. This transformation does not take place spontaneously but typically requires external energy and, in some cases, specific conditions. Energy is necessary to break certain bonds to allow the molecule to restructure into the corresponding isomer.

Procedure

Part A. Identify the products by their melting point.

- ☐ *Step 1:* – Obtain the two unknowns. One unknown should be fumaric acid (*white powder*), and the other maleic acid (*white powder*). If both unknowns show the same properties in the following experiment, be prepared to select a different unknown.
- ☐ *Step 2:* – Label two small test tubes as A and B.
- ☐ *Step 3:* – Put a pea-size amount of one product in test tube A, and a similar amount of the other product in test tube B.
- ☐ *Step 4:* – Using a striker, light the Bunsen burner. Using a test tube holder for each sample, hold test tubes closely together, in the same hand. The samples must be heated in similar conditions. Make sure you never heat a test tube with a tube stopper on as they will explode.
- ☐ *Step 5:* – Write your results in the Results section, characterizing the melting point of each unknown by comparison as low or high.

Good Lab Practice

-  The burner should stand stable with no additional support. Coil the tubing at each end to stabilize the burner.
-  Work safely and avoid burns.
-  Always ensure easy access to the gas valve. Remove any items around the valve, and do not set the flame on your way.
-  Never leave the flame unattended. In case of accident or any risk, close the gas immediately.

- ☐ *Step 6:* – Start a timer. Swing the test tubes over the flame. Stop as soon as one of the two powders melts. Pay attention to the test tubes' label. Was it test tube A or test tube B? Write down the time it took to melt the isomer on the Results table.
- ☐ *Step 7:* – Refer to the known properties of maleic acid and fumaric acid (from the prelab) to determine which product is unknown A and which is unknown B.
- ☐ *Step 8:* – Allow the test tubes to cool down. Discard the contents in the indicated waste containers. Clean the test tubes.

Part B. Identify the products by their solubility in water.

- ☐ *Step 1:* – Label two small test tubes as A and B.
- ☐ *Step 2:* – Put a pea-size amount of one product in test tube A, and a similar amount of the other product in test tube B.
- ☐ *Step 3:* – Add similar amounts of distilled water to each test tube. Approximately 2 mL of water should be enough.
- ☐ *Step 4:* – Cover the tube using the stoppers and shake them until one of the 2 powders dissolves completely in the liquid.
- ☐ *Step 5:* – Refer to the known properties of maleic acid and fumaric acid (from the prelab) to determine which product is unknown A and which is unknown B.
- ☐ *Step 6:* – Your result should be consistent with Part A. If not, first repeat part B to double-check if the labeling was correct. If the mistake was not in part B, then repeat part A.
- ☐ *Step 7:* – Write your results in the Results section, characterizing the solubility of each unknown by comparison as low or high.

Part C. Isomerization reaction. Preparing the setup

- ☐ *Step 1:* – Obtain **about** 1.0 grams of the most soluble sample in Part B, using the scale and a weighing boat.
- ☐ *Step 2:* – Transfer the powder into a 250 mL Erlenmeyer flask and add **around** 10 mL of distilled water.
- ☐ *Step 3:* – Now, we will prepare a setup to boil the solution that will remove the vapors from the reaction.
- ☐ *Step 4:* – Set up a ring stand and the iron ring with the wire gauze on top. Do not put the Erlenmeyer with the solution yet. Make sure that the height of the platform is appropriate for the size of the flame.
- ☐ *Step 5:* – Use a second, smaller ring above the platform as protection.
- ☐ *Step 6:* – To remove the vapors improvise a little hood using a plastic funnel. Connect the funnel to the vacuum line using the rubber tubing and set it with a clamp, upside down, at a certain distance over the platform.
- ☐ *Step 7:* – Start the ventilation.
- ☐ *Step 8:* – Show the setup to your instructor, who will make sure it is sturdy enough. In particular, the instructor will sure all vapor is removed with the funnel, if not, adjust the setup right away. You want to avoid the release of any vapor into the air.
- ☐ *Step 9:* – Place the maleic acid solution in the Erlenmeyer, and light the burner to gently warming up the solution.

⚠ CAUTION!

- ⚠ Wear goggles at all times.
- ⚠ Handle concentrated acid with care to avoid chemical burns.
- ⚠ Use latex gloves when handling concentrated acid.
- ⚠ Ensure that the mouth of the test tube is pointed away from you and any other persons around you, including those on the opposite side of the bench.

Part C. Isomerization reaction. Adding the concentrated acid

- ☐ *Step 1:* – Make sure you wear some protective latex gloves.
- ☐ *Step 2:* – While the solution warms up obtain 10 mL (approx.) of concentrated hydrochloric acid, HCl, in a 25 mL graduated cylinder. Use a plastic dropper to transfer the acid from its container into the cylinder. Mind this is a very concentrated acid that can hurt you.
- ☐ *Step 3:* – Keep the cylinder all the time in the hood to minimize the release of acidic vapors.
- ☐ *Step 4:* – When the solution starts to boil, bring the acid from the hood, and add *carefully* the acid into the solution using a plastic dropper.
- ☐ *Step 5:* – There should be no acid left over. If so, get rid of the leftover acid right away.
- ☐ *Step 6:* – Wait until the boiling is re-established, and make sure no vapor escapes the ventilation system.
- ☐ *Step 7:* – Let the reaction boil for 1 minute and then turn off the gas. Continue preparing the next steps.

Part C. Isomerization reaction. Filtration

- ☐ *Step 1:* – Put 40 mL of distilled water in a clean 50 mL beaker. Cool down this water using an ice bath. Get ice in a large enough beaker and dip the 50 mL beaker in the ice.
- ☐ *Step 2:* – Go back to the Erlenmeyer with the solid dissolved and the acid. Allow the Erlenmeyer to cool down for 10 minutes. In the meantime, obtain ice in a 600 mL beaker and add some water to the ice.
- ☐ *Step 3:* – After the 10 minutes of cooling down, transfer the Erlenmeyer carefully and safely into the large ice bath.
- ☐ *Step 4:* – Get the suction filtration setup ready. Connect the suction flask to the vacuum and put the Büchner funnel on top.
- ☐ *Step 5:* – Obtain a flat, round filter paper with a diameter slightly smaller than the bottom of the funnel. Start the vacuum and add some distilled water with a wash bottle to test the suction.
- ☐ *Step 6:* – Filter the solution. Stir the mixture to help the precipitate to suspend and add the liquid carefully over the filter paper. Use the cold distilled water in the 50 mL beaker to remove any solid left in the Erlenmeyer and to wash the filtrate.
- ☐ *Step 7:* – Keep the suction on for another 5 minutes. This will help remove all the water and partially dry the filter.
- ☐ *Step 8:* – Turn off the vacuum and remove the filter from the funnel carefully.
- ☐ *Step 9:* – Keep the solid product until the end of the experiment.

Part C. Isomerization reaction. Testing the product

- ☐ *Step 1:* – Identify your product following the steps in Part B.
- ☐ *Step 2:* – Write your results in the Results section, characterizing the solubility of the isomerization product by comparison as low or high.

STUDENT INFO

Name:

Date:

Pre-lab Questions**Geometric isomers**

1. Complete the table for the molecules Fumaric Acid and Maleic Acid

	Fumaric Acid	Maleic Acid
Molecular formula		
Draw the structure		
Solubility in water		
Melting point		
What is it used for?		
Where is it found?		

2. Draw the setup for suction filtration (Büchner funnel, suction flask, tubing...)
3. What are the dangers of an improper installation of the ventilation system in this experiment? What chemical in this experiment can be dangerous?

4. Define:

(a) Isomers

(b) Geometric Isomers

(c) Isomerization

(d) Suction filtration

(e) Filtrate

5. When you mix hydrochloric acid and water, what do you add first the acid or water?

STUDENT INFO

Name:

Date:

**Results
EXPERIMENT****Geometric isomers**

This experiment is qualitative, which means that the results are based on observations. Record your observations for each sample during and after each part of the activity.

	Before isomerization		After isomerization
	Sample A	Sample B	
Part A. Time to melt (secs)			
Part A. Melting point (lower/higher)			
Part B. Solubility in water (smaller/larger)			
Isomer name (fumaric/maleic)			
Part C. Selected for isomerization? (yes/no)			

STUDENT INFO

Name:

Date:

Post-lab Questions**Geometric isomers**

1. A Similar amount of both unknowns are used for their characterization in Parts A and B. Why is it important to use similar amounts? What would change if the amounts used were very different?
2. How can you prove that the isomerization process was successful?
3. Why did you cool down the solution after the isomerization?
4. What would happen if hot water was used to wash the filtrate?

2 Conductivity of electrolyte solutions

EXPERIMENT

Conductivity of electrolyte solutions

Goal

The terms of **conductivity** and **electrolytes** will be introduced and related to the concept of dissociation. The use of the conductivity meter, very similar to the pH meter, will help to differentiate between strong and weak electrolytes. In the last part of the experiment, the presence of chemical reactions will be identified by measuring changes in conductivity.

Materials

- | | |
|--|---|
| <input type="checkbox"/> jumbo test tubes | <input type="checkbox"/> 0.1M solutions of KNO_3 , KCl , KOH , HCl , and $\text{Ca}(\text{NO}_3)_2$ |
| <input type="checkbox"/> conductivity meter | <input type="checkbox"/> an unknown solution |
| <input type="checkbox"/> 6M solution of NH_3 | <input type="checkbox"/> 5mL, 10mL, 25mL pipets (preferred Mohr's) |
| <input type="checkbox"/> 6M solution of CH_3COOH | <input type="checkbox"/> 250mL and 50mL beakers |
| <input type="checkbox"/> Ethyl alcohol | <input type="checkbox"/> rinsing water bottle |

Background

Playing with water near electrical appliances can be risky. Surprisingly, however, water does not conduct electricity. For electricity to propagate through media there must be free charges that move along the electric field. Pure water does not contain a significant number of charges and does not conduct electricity. However, when there is something dissolved in the water, therefore making an aqueous solution, it can be conductive.

Conductivity and aqueous solutions

Whether the aqueous solution is conductive or not depends on the nature of the solute. If the solute dissociates into ions when dissolved in water, the ions will conduct electricity. These types of chemicals are called electrolytes. If the solute dissolves as a molecule without dissociation, the solution will not conduct electricity and the solute will be called nonelectrolyte.

Conductivity units

The unit for conductivity is Siemens per centimeter, $\text{S}\cdot\text{cm}^{-1}$. Pay attention to the units during this experiment because the tool will often display prefixes like milliSiemens or microSiemens per centimeter. Temperature affects conductivity and must be recorded during the measurement of this property.

Types of electrolytes and factors affecting their conductivity.

Salt and sugar are soluble in water. When sugar dissolves in water its molecules remain intact because water can not break apart the strong covalent bonds between their atoms. As a neutral molecule, sugar in solution will not conduct electricity and it is an example of a **nonelectrolyte**. Salt, sodium chloride or NaCl , is an ionic compound that gives place to two ions, Na^+ and Cl^- when dissolved in water. These ions are responsible for the solution's conductivity. The number or more specifically

the concentration, and the mobility of these ions will determine how conductive the solution is. Based on these principles we can classify substances as:

- × **Non-electrolytes.** When the aqueous solution of that substance does not conduct electricity. In general terms, these electrolytes present low conductivities, lower than $10\mu\text{S}\cdot\text{cm}^{-1}$.
- × **Weak electrolyte.** When the aqueous solution of that substance conducts electricity poorly. In general terms, these electrolytes present medium conductivities, between $10\text{-}1000\mu\text{S}\cdot\text{cm}^{-1}$.
- × **Strong electrolyte.** When the conductivity of the aqueous solution of that substance is high. In general terms, these electrolytes present large conductivities, larger than $1000\mu\text{S}\cdot\text{cm}^{-1}$.

The **mobility** of the ions is related to their size. Small ions move fast which enhances conductivity, while large ions encounter more resistance to move and show lower conductivities. The table below displays a series of mobility values for positive and negative ions at infinite dilution and 298K and 1 atm.

Ion	H_3O^+	Li^+	Na^+	Mg^{2+}	OH^-	Cl^-	Br^-	NO_3^-
movility, u^∞ ($10^5\text{cm}^2\cdot\text{V}^{-1}\cdot\text{s}^{-1}$)	363	41	52	55	206	80	81	74

The effect of ions' mobility is easily appreciated in isomolar solutions, which are solutions with the same molarity. The more ions in the solution, the higher the conductivity will be. For the same substance, the number of ions is determined by the **concentration**; more concentration means more ions and more conductivity. It is important to recall that only ions contribute to conductivity. While some molecules dissociate completely, others dissociate only partially. The latter are weak electrolytes. Interestingly, when two solutions containing ions are mixed, the conductivity of the resulting solution will be equal to the sum of the two separate solutions with the same concentration as in the mixture. This is called the **additivity rule** and will hold as long as the ions do not react with each other. When two electrolytes undergo a chemical reaction, new products will form and the additivity rule will not apply. A large deviation from the sum of separate conductivities is an indication of a chemical reaction.

Example

The conductivity of a 0.050 M solutions of HNO_3 is measure to be 19.6 mS/cm. That of a 0.050 M solution of KCl is 9.6 mS/cm. A new solution is prepared mixing 10 mL of HNO_3 0.10 M and 10 mL of KCl 0.10 M. The conductivity of the resulting mixture is measure to be 29.6 mS/cm. Calculate the concentration of each spices in the mixture and determine if their compounds have reacted or not.

Answer: To calculate the new concentration use the formula for dilutions:

$$V_1 \times M_1 = V_2 \times M_2$$

where V_1 and M_1 are the initial volume and concentration and V_2 and M_2 are the values after the dilution. We look for M_2 , therefore

$$M_2 = V_1 \times \frac{M_1}{V_2}$$

for HNO_3 ,

$$M_2 = \frac{10\text{ mL} \times 0.10\text{ M}}{(10\text{ mL} + 10\text{ mL})} = 0.05\text{ M}$$

The same values apply for KCl .

We can see that the addition of the separate conductivities (κ) is very close to the total conductivity of the mixture with the same concentrations for each electrolyte.

$$\begin{aligned}\kappa(\text{HNO}_3\ 0.05\text{ M}) + \kappa(\text{KCl}\ 0.05\text{ M}) &= 19.6\text{ mS/cm} + 9.6\text{ mS/cm} = 29.2\text{ mS/cm} \\ \kappa(\text{HNO}_3\ 0.05\text{M} + \text{KCl}\ 0.05\text{M}) &= 29.6\text{ mS/cm}\end{aligned}$$

The additivity rules applies, meaning that there is no reaction between HNO_3 and KCl .

Procedure



Part A. Getting started: standardize the conductimeter

- ☐ *Step 1:* – Install the conductivity meter. The probe of the instrument is immersed in a clean solution that must be kept aside until the end of the experiment. Mind the probe needs to be hydrated at all times.
- ☐ *Step 2:* – To calibrate the conductivity meter, place 40mL of standard conductivity solution in a 50mL beaker while placing the probe in it.
- ☐ *Step 3:* – Make sure the conductivity cell is fully submerged in the liquid.
- ☐ *Step 4:* – Now press Setup (top right menu). Normal setups are: cell constant of 1.0, standard recognition with manual selected, reference temperature 25 degrees, temperature coefficient of 2.1%, and alarm limits off.
- ☐ *Step 5:* – Now press Standardize (top right menu) and input 1.000 mS/cm (or 1000 $\mu\text{S}/\text{cm}$) in the keypad and press confirm (this should be the conductivity value of the standard). Finally, press Confirm (top right menu) to save the standard value. Now you are ready to measure.

Part A. Getting started: measuring conductivity

- ☐ *Step 1:* – It is very important to follow some cleaning steps before you measure the conductivity of any solution. You will need a 100 mL beaker labeled as waste, a wash bottle, and a 250 mL beaker with distilled water.
 1. Rinse the probe using the wash bottle over the waste beaker.
 2. Dip the probe 3-4 times into the 250 mL beaker with clean distilled water.
 3. Dip the probe into the solution to be measured. Use a 50 mL beaker (or a jumbo test tube) with at least 25 mL of solution. Use the same volume for all conductivity measurements.
 4. Repeat cleaning steps 1 and 2.
- ☐ *Step 2:* – It is critical to properly read the conductivity units from the conductimeter. Every solution would give a different conductivity value and sometimes a different unit. The conductivity units are written next to the measurement in the conductivity. It is critical you copy the measurement with the unit given by the tool.

Good Lab Practice

-  To ensure the probe is well cleaned and the water in the 250 mL beaker remains clean, the conductivity meter should display a value close to that of the distilled water during the cleaning.
-  If the conductivity value of the cleaning water differs from the one of distilled water, replace the water in the 250 mL beaker with fresh distilled water, and thoroughly clean the probe.

Part B. Pure substances.

- ☐ *Step 1:* – Following the procedure outlined above, measure the conductivity of distilled water. Record the value in the Results section.
- ☐ *Step 2:* – Following the procedure outlined above, measure the conductivity of tap water. Record the value in the Results section.

- ☐ *Step 3:* – Following the procedure outlined above, measure the conductivity of Ethyl alcohol ($\text{C}_2\text{H}_5\text{OH}$). Record the value in the Results section.
- ☐ *Step 4:* – Classify the substance as strong, weak, or nonelectrolyte.

Part C. Series of electrolytes.

- ☐ *Step 5:* – In a clean 50 mL beaker, add 20 mL of 0.10 M HCl and 20 mL of distilled water (this is a $\frac{1}{2}$ dilution). Use a Mohrs pipet to transfer the liquids, making sure you rinse with distilled water between transfers. Mix well and measure the conductivity. Recalculate the concentration after the dilution and write it down in the Results section.
- ☐ *Step 6:* – In a clean 50 mL beaker, add 20 mL of 0.10 M KOH and 20 mL of distilled water (this is a $\frac{1}{2}$ dilution). Use a Mohrs pipet to transfer the liquids, making sure you rinse with distilled water between transfers. Mix well and measure the conductivity. Recalculate the concentration after the dilution and write it down in the Results section.
- ☐ *Step 7:* – In a clean 50 mL beaker, add 20 mL of 0.10 M KCl and 20 mL of distilled water (this is a $\frac{1}{2}$ dilution). Use a Mohrs pipet to transfer the liquids, making sure you rinse with distilled water between transfers. Mix well and measure the conductivity. Recalculate the concentration after the dilution and write it down in the Results section.
- ☐ *Step 8:* – In a clean 50 mL beaker, add 20 mL of 0.10 M KNO_3 and 20 mL of distilled water (this is a $\frac{1}{2}$ dilution). Use a Mohrs pipet to transfer the liquids, making sure you rinse with distilled water between transfers. Mix well and measure the conductivity. Recalculate the concentration after the dilution and write it down in the Results section.
- ☐ *Step 9:* – In a clean 50 mL beaker, add 20 mL of 0.10 M $\text{Ca}(\text{NO}_3)_2$ and 20 mL of distilled water (this is a $\frac{1}{2}$ dilution). Use a Mohrs pipet to transfer the liquids, making sure you rinse with distilled water between transfers. Mix well and measure the conductivity. Recalculate the concentration after the dilution and write it down in the Results section.
- ☐ *Step 10:* – Find an unknown solution. In a clean 50 mL beaker, add 20 mL of the unknown solution and 20 mL of distilled water (this is a $\frac{1}{2}$ dilution). Use a Mohrs pipet to transfer the liquids, making sure you rinse with distilled water between transfers. Mix well and measure the conductivity.
- ☐ *Step 11:* – Classify all substances as strong, weak, or nonelectrolytes and identify the unknown compound based on the conductivity values. Know that the unknown can also be water.

Part D. Series of Dilutions.

- ☐ *Step 12:* – Obtained the conductivity of a 0.1 M solution of HCl from previous section. Write down the result in the Results section.
- ☐ *Step 13:* – Now you are going to prepare two dilutions. Starting from the concentrated 0.1 M solution of HCl, you are going to prepare 40mL volumes of the following dilutions: 0.050 M and 0.020 M. The first is a $\frac{1}{2}$ dilution, that is 20ml of acid needs to be mixed with 20ml of water to get 40mL of 0.050M solution, whereas the second is a $\frac{1}{5}$ dilution, that is 8ml of acid needs to be mixed with 32ml of water. Use a Mohrs pipet to transfer the liquids, making sure you rinse with distilled water between transfers.
- ☐ *Step 14:* – Measure the conductivity of the two diluted HCl solutions (0.050 M, and 0.020 M) and record the values.

Part E. Mixtures of HCl and KNO_3 .

- ☐ *Step 15:* – Write down the conductivity of the 0.050 M solution of HCl from the previous section in the Results section.
- ☐ *Step 16:* – Write down the conductivity of the solution of KNO_3 from the previous section in the Results section.
- ☐ *Step 17:* – In a clean 50 mL beaker, add 20 mL of 0.10 M HCl and 20 mL of 0.10 M KNO_3 . Use a Mohrs pipet to transfer the liquids, making sure you rinse with distilled water between transfers. Mix well and measure the conductivity. Recalculate the concentration after the dilution and write it down in the Results section, minding that both solutions have now the same concentration and therefore only one number needs to be listed.

Part F. Mixtures of NH_3 and $\text{HC}_2\text{H}_3\text{O}_2$.

- ☐ *Step 18:* – Prepare the following dilutions: 100mL of acetic acid ($\text{HC}_2\text{H}_3\text{O}_2$) 0.12 M from 6M (add 2mL of concentrated solution and fill up with water until 100mL) and 40mL of ammonia (NH_3) 0.12 M from 6M (add 2mL of concentrated solution and fill up with water until 100mL). Carefully read the concentration of the initial concentrated solutions indicated in the labels. Work in the hood when handling highly concentrated solutions.
- ☐ *Step 19:* – In a clean 50 mL beaker, add 20 mL of the 0.12 M NH_3 dilution and 20 mL of distilled water (this is a $\frac{1}{2}$ dilution). Use a Mohrs pipet to transfer the liquids, making sure you rinse with distilled water between transfers. Mix well and measure the conductivity. Recalculate the concentration after the dilution and write it down in the Results section. Work in the hood when handling highly concentrated solutions.
- ☐ *Step 20:* – In a clean 50 mL beaker, add 20 mL of 0.12 M $\text{HC}_2\text{H}_3\text{O}_2$ and 20 mL of distilled water (this is a $\frac{1}{2}$ dilution). Use a Mohrs pipet to transfer the liquids, making sure you rinse with distilled water between transfers. Mix well and measure the conductivity. Recalculate the concentration after the dilution and write it down in the Results section.
- ☐ *Step 21:* – In a clean 50 mL beaker, add 20 mL of 0.12 M NH_3 and 20 mL of 0.12 M $\text{HC}_2\text{H}_3\text{O}_2$. Use a Mohrs pipet to transfer the liquids, making sure you rinse with distilled water between transfers. Mix well and measure the conductivity. Recalculate the concentration after the dilution and write it down in the Results section, minding that both solutions have now the same concentration and therefore only one number needs to be listed.

⚠ CAUTION!

- ⚠ Handle concentrated acid or bases with care to avoid chemical burns.
- ⚠ Are you wearing your goggles? Do you know where the eye-washer is?

STUDENT INFO

Name:

Date:

Pre-lab Questions

Conductivity of electrolyte solutions

1. Starting with a 6 M solution of NaCl, calculate the volume of solution and the volume of water necessary to prepare 10 mL of a 3 M solution of NaCl.
2. How can you dilute a 3 M solution to form a 1.5 M solution? (e.g. I would pick up X mL of the solution and YmL of distilled water.)
3. How can you prepare 40mL of acetic acid ($\text{HC}_2\text{H}_3\text{O}_2$) 0.12 M from 6M.
4. How can you prepare 40mL of ammonia (NH_3) 0.12 M from 6M.
5. Given the data below, classify the following electrolytes as strong, weak, or nonelectrolytes:

	Concentration (M)	Conductivity	Units	Type of electrolyte
CuSO_4	0.05M	4000	$\mu\text{S} \cdot \text{cm}^{-1}$	_____
FeCl_3	0.05M	10000	$\mu\text{S} \cdot \text{cm}^{-1}$	_____
$\text{C}_6\text{H}_{12}\text{O}_6$	0.05M	9	$\mu\text{S} \cdot \text{cm}^{-1}$	_____
CH_3COOH	0.05M	320	$\mu\text{S} \cdot \text{cm}^{-1}$	_____

STUDENT INFO

Name:

Date:

Results
EXPERIMENT

Conductivity of electrolyte solutions

	Conductivity	Units	Electrolyte type
Distilled water	_____	_____	_____
Tap water	_____	_____	_____
C ₂ H ₅ OH	_____	_____	_____

Part C: Electrolytes

	Concentration (M)	Conductivity	Units	Electrolyte type
HCl	_____	_____	_____	_____
KOH	_____	_____	_____	_____
KCl	_____	_____	_____	_____
KNO ₃	_____	_____	_____	_____
Ca(NO ₃) ₂	_____	_____	_____	_____
Unknown		_____	_____	
Unknown #	_____	Unknown identity	_____	

Part D: Concentration effect

	Concentration (M)	Conductivity	Units
HCl	_____	_____	_____
HCl	_____	_____	_____
HCl	_____	_____	_____

Part E: HCl+KNO₃ mixture

	Concentration (M)	Conductivity	Units
HCl	_____	_____	_____
KNO ₃	_____	_____	_____
HCl+KNO ₃	_____	_____	_____

Part F: HC₂H₃O₂+NH₃ mixture

	Concentration (M)	Conductivity	Units
HC ₂ H ₃ O ₂	_____	_____	_____
NH ₃	_____	_____	_____
HC ₂ H ₃ O ₂ +NH ₃	_____	_____	_____

STUDENT INFO

Name:

Date:

Post-lab Questions

Conductivity of electrolyte solutions

1. Based on your results on Part B, which conductivity was higher, the one for distilled water or for tap water? Why?
2. Based on your results on Part C, specifically the conductivity values of KCl and KOH, how does ion mobility affect conductivity (see table in the Background section)? Use your data to illustrate your answer.
3. Based on your results on Part D, how does concentration affect conductivity? Use your data to illustrate your answer.
4. Based on your results in Part E, determine if their compounds have reacted or not. Does the additivity rule hold for this case?
5. Based on your results in Part F, determine if their compounds have reacted or not. Does the additivity rule hold for this case?

3 Softening of hard water

EXPERIMENT

Softening of hard water

Goal

The goal of this experiment is to understand the concepts of hard and soft water and their implications in daily life. The student will also experiment with the *dissociation of ionic compounds* in water, from which the **hydrogen-donor** nature of the acids and the notion of *pH* are emphasized. Other principles, such as equilibrium and ion exchange, are also presented.

Materials

- | | |
|--|---|
| <input type="checkbox"/> 250mL and 50mL beakers | <input type="checkbox"/> $\text{NH}_3/\text{NH}_4\text{Cl}$ buffer |
| <input type="checkbox"/> cation-exchange resin | <input type="checkbox"/> 0.01M $\text{Na}_2\text{H}_2(\text{EDTA})$ |
| <input type="checkbox"/> 6M HCl | <input type="checkbox"/> indicator (Eriochrome black mixed with NaCl) |
| <input type="checkbox"/> distilled water (not tap water) | <input type="checkbox"/> dropper solutions (0.10 M $\text{Ca}(\text{NO}_3)_2$, 0.10 M Na_2CO_3 , and 0.10 M $\text{Mg}(\text{NO}_3)_2$) |
| <input type="checkbox"/> Litmus paper | <input type="checkbox"/> 0.10 M NaHCO_3 solution |
| <input type="checkbox"/> 50mL graduated cylinder | <input type="checkbox"/> solid CaCO_3 |
| <input type="checkbox"/> 125mL Erlenmeyer | |

Background


The term "**hard water**" refers to the concentration of certain minerals in the water. This topic concerns **aqueous solutions** in which water is the solvent and minerals are the solutes. These minerals are ionic compounds that dissociate in aqueous solutions producing ions. The hardness of water is an indication of the concentration of Ca^{2+} and Mg^{2+} ions.

Dissolving the insoluble


Ca^{2+} and Mg^{2+} ions originate in limestone (CaCO_3 , the material of the pyramids) and chalk (MgCO_3 , used to write on the chalkboard) deposits. However, chalk will not dissolve in water and neither the pyramids dissolve in the rain. How is dissolving these carbonates possible?

The answer is in the air


Water in contact with the air absorbs CO_2 in the atmosphere according to the equation:




Where H_2CO_3 is carbonic acid, which reacts with the carbonate to form hydrogen bicarbonate. For the case of CaCO_3



Calcium hydrogen carbonate is also called calcium bicarbonate. Bicarbonates are soluble in water (baking soda is sodium bicarbonate, NaCO_3) and, as ionic compounds, they are dissociated into its ions;




The overall reaction is the 3 previous equations added together:



and this is how carbonates dissolve in water.

Sodium hydrogen carbonate can decompose under different conditions. On one hand, in acidic conditions, it produces carbonic acid that quickly decomposes to give carbon dioxide:



On the other side, when heated it produces: $2\text{NaHCO}_3 (\text{aq}) \xrightarrow{\Delta} \text{Na}_2\text{CO}_3 + \text{CO}_2 (\text{g}) + \text{H}_2\text{O} (\text{l})$

Implications in daily life


The presence of these ions in water has several implications for routine actions in our lives. The problem arises when the latter reaction takes place in the opposite direction, releasing CO_2 gas in the process and depositing solid carbonate. This direction is favored when water is heated and explains the white residue in bathrooms and kettles. This is especially a problem for boilers, that suffer from scaling when operating hard water.

The experiment


The experiment consists of four parts: preparing the ion exchange resin, measuring the hardness of a sample of water, reducing the hardness also called softening, of the sample, and testing the hardness of the softened sample. Additionally, a series of simple reactions will be done to test the reactivity of Ca^{2+} and Mg^{2+} ions.

Measuring the hardness of the sample

The concentration of Ca^{2+} ions, or any other metallic cation, is usually measured by titration of the aqueous sample using a sodium salt as the titrant, a solution of known concentration added to the sample. First, an indicator will be added to the sample which will color the solution pink. The indicator chosen is Eriochrome Black T, which will turn blue when the Ca^{2+} ion is removed from the cation-indicator complex:



To achieve the color change $\text{H}_2\text{EDTA}^{2-}$ will be used. $\text{H}_2\text{EDTA}^{2-}$ (dihydrogen ethylenediaminetetraacetate) is an anion with a stronger affinity for the Ca^{2+} cations than the indicator. As $\text{H}_2\text{EDTA}^{2-}$ is added to the solution it will replace the indicator forming a new complex, leaving the indicator alone, which in turn will become blue.



Even when the example described above refers to calcium ions, the same chemistry applies to magnesium ions. It is important to note down the exact number of drops required to achieve the endpoint of the titration.

Softening the hard water



After measuring the hardness of the water sample, the next step is to soften the water. *Softening* is the process where ions, particularly calcium and magnesium ions, are removed from the water. When a hard water sample is softened we can refer to it as *softened water*. A cation-exchange resin will be used to soften the water. This resin contains many sulfonic groups, ($-\text{SO}_3\text{H}$). Every two hydrogens in these acidic groups will be replaced by a metallic cation, Ca^{2+} and Mg^{2+} , in a mechanism called *cation exchange*. Finally, it is necessary to evaluate whether the water sample has been softened. As such the titration performed in Part A will be repeated, but this time with the softened sample. The most relevant parameter is the amount of EDTA added to achieve the endpoint of the titration, to compare this result to the previous titration. This experiment aims to be a qualitative analysis of the water treatment. For that reason, and to simplify the process, the titration will be done using droppers instead of the more precise burette typically used in titrations.

Procedure

Part A. Preparing the resin.

- ☐ *Step 1:* – Obtain around 5g of the dry, ion exchange resin in a 250 mL beaker. Remember that this resin can be recycled. Do not discard it at the end of the experiment.
- ☐ *Step 2:* – Add enough 6 M HCl to cover the resin. By adding the acid, the resin is cleaned of metallic cations. To use the ion exchange resin, it must be clean, also known as "in H^+ form".

CAUTION!

-  Handle concentrated acid with care to avoid chemical burns.
-  Are you wearing your goggles? Do you know where the eye-washer is?

- ☐ *Step 3:* – Now you are going to wash the acid away. After 2 minutes, add 200 mL of distilled water to the resin-acid mixture. Allow the resin to settle in the bottom of the beaker.
- ☐ *Step 4:* – Carefully, pour the liquid off (decant) into another beaker. This liquid is an HCl solution and should be handled and discarded as such.
- ☐ *Step 5:* – To ensure that all the acid is gone, wash again the resin with 200 mL of distilled water. Allow the resin to settle in the bottom of the beaker and pour the liquid into another beaker.
- ☐ *Step 6:* – Before you discard the liquid you can test the pH. Dip the tip of a glass stirring rod in the liquid and poke the wet end of the rod into a piece of Blue Litmus paper. If the Litmus paper turns pink, there is still acid, repeat steps 5 and 6. If it does not turn pink, the resin is clean and ready to be used. Reserve the resin for Part B.

Part B. Measure the hardness of a sample of softened water

- ☐ *Step 1:* – Obtain a little bit more than 150 mL of hard water in a beaker. Hard water is already prepared in the lab.
- ☐ *Step 2:* – Using a graduated cylinder, add 20.0 mL of hard water to a 125 mL Erlenmeyer. Reserve this water for Part C. Use the rest of the water for this part of the experiment.
- ☐ *Step 3:* – Transfer the remaining hard water to the beaker with the resin and let it stand for 25 minutes. In the meantime, proceed with Parts C and D.
- ☐ *Step 4:* – Pour off about 40 mL of the water sample that has been in contact with the resin into a clean beaker. This sample would be referred to as the softened water.
- ☐ *Step 5:* – Using a graduated cylinder, add 20.0 mL of the softened water to a 125 mL Erlenmeyer.

- ☐ *Step 6:* – Add 5 mL of an $\text{NH}_3/\text{NH}_4\text{Cl}$ buffer solution (not to be confused with $\text{Na}_2\text{H}_2(\text{EDTA})$) to the Erlenmeyer and the tip of a spatula (less than a pea-sized amount) of the indicator (a solid mixture of eriochrome black T and sodium chloride). Adding too many indicators will ruin the experiment. The solution color will be rose pink after the addition of the indicator.
- ☐ *Step 7:* – Now, we are going to titrate the water sample. Obtain about 20 mL of 0.01 M $\text{Na}_2\text{H}_2(\text{EDTA})$ (the full name is Ethylenediaminetetraacetic acid) in a 50 mL beaker and a medicine dropper (or a plastic Pasteur pipette). Add the $\text{Na}_2\text{H}_2(\text{EDTA})$ drop by drop to the mixture while stirring. As the endpoint is approached the solution color will become lavender. Count the drops needed for the analyte to turn blue. When the addition of one drop turns the solution blue, all Ca_2^+ and Mg_2^+ ions have been removed from the sample
- ☐ *Step 8:* – Clean the resin as described in Part A, and return the resin to the instructor—there is a large beaker ready to collect wet resin.

Part C. Measure the hardness of a water sample

- ☐ *Step 1:* – While Part B takes place, and the water becomes softer, add 5 mL of a $\text{NH}_3/\text{NH}_4\text{Cl}$ buffer solution (not to be confused with $\text{Na}_2\text{H}_2(\text{EDTA})$) to the Erlenmeyer, and the tip of a spatula (less than a pea-sized amount) of the indicator. Remember that the indicator is just a solid mixture of eriochrome black T and sodium chloride. Adding too many indicators will ruin the experiment. The solution color will be rose pink after the addition of the indicator.
- ☐ *Step 2:* – Now, we are going to titrate the water sample. Obtain about 20 mL of 0.01 M $\text{Na}_2\text{H}_2(\text{EDTA})$ in a 50 mL beaker and a medicine dropper (or a plastic Pasteur pipette). Add the $\text{Na}_2\text{H}_2(\text{EDTA})$ drop by drop to the mixture while stirring. As the endpoint is approached the solution color will become lavender. Count the drops needed for the analyte to turn blue. When the addition of one drop turns the solution blue, all Ca_2^+ and Mg_2^+ ions have been removed from the sample

Part D. Observing reactions with Ca^{2+} and Mg^{2+} ions

- ☐ *Step 1:* – Obtain 4 test tubes and locate the following solutions: 0.10 M $\text{Ca}(\text{NO}_3)_2$, 0.10 M Na_2CO_3 , 0.10 M NaHCO_3 , and 0.10 M $\text{Mg}(\text{NO}_3)_2$.
- ☐ *Step 2:* – In test tube 1, add 10 drops of $\text{Ca}(\text{NO}_3)_2$. Add 10 drops of Na_2CO_3 and write down your observations.
- ☐ *Step 3:* – In test tube 2, add 10 drops of $\text{Ca}(\text{NO}_3)_2$. Add 2 drops of HCl and write down your observations.
- ☐ *Step 4:* – In test tube 3, add 10 drops of $\text{Ca}(\text{NO}_3)_2$. Add 10 drops of 0.1 M NaHCO_3 and write down your observations.
- ☐ *Step 5:* – Now, Heat test tube 3 in the flame of a Bunsen burner. Write down your observations.
- ☐ *Step 6:* – Now, cool down test tube 3. Add 2 drops of 6M HCl and write down your observations.
- ☐ *Step 7:* – In test tube 4, add a pea-sized amount of solid CaCO_3 . Add 2 drops of 6M HCl and write down your observations.
- ☐ *Step 8:* – Now, repeat the steps above using $\text{Mg}(\text{NO}_3)_2$ instead of $\text{Ca}(\text{NO}_3)_2$.
- ☐ *Step 9:* – At this point, complete the remaining steps of Part B.

STUDENT INFO

Name:

Date:

Pre-lab Questions

Softening of hard water

1. A student needs 10 drops of EDTA to reach the equivalency point when titrating a sample of hard water. After the students soften the water with the resin, it requires 12 drops of EDTA to reach the equivalency. Are these results correct? Explain.
2. Find the definition of the following concepts:
 - (a) Titration
 - (b) Titrant
 - (c) Litmus paper
 - (d) Ion exchange resin
3. Where is the term **hard water** derived from?
4. Compare the sizes used in the experiment. Circle the bigger size:
 - (a) Pea-size
 - (b) Tip of a spatula
5. In the experiment, find the chemical name of EDTA.

STUDENT INFO

Name:

Date:

**Results
EXPERIMENT**

Softening of hard water

Parts B and C: Comparing titrations

Indicate how many drops you used to achieve the endpoint for:

Hard water sample (Part C): _____ drops of $\text{Na}_2\text{H}_2(\text{EDTA})$ Softened water sample (Part B): _____ drops of $\text{Na}_2\text{H}_2(\text{EDTA})$ **Part D: Observing reactions with Ca^{2+} and Mg^{2+} ions**

Write down your observations such as gas evolution, the appearance of precipitate, dissolving of precipitate, and no observation.

	$\text{Ca}(\text{NO}_3)_2$	$\text{Mg}(\text{NO}_3)_2$
Na_2CO_3		
HCl		
NaHCO_3		
$\text{NaHCO}_3 + \text{heat}$		
$\text{NaHCO}_3 + \text{heat} + \text{HCl}$		
Solid $\text{CaCO}_3 + \text{HCl}$		

STUDENT INFO






Name:

Date:


Post-lab Questions**Softening of hard water**



1. Based on your measurements. Did the ion exchange resin soften the hard water sample? Explain.

2. Write the balanced equations for the chemical reactions below. Mind that calcium carbonate is an insoluble compound.



3. Write the balanced equations for the chemical reactions below. Mind that magnesium carbonate is an insoluble compound.





4. Using the equivalency between drops and mL (1mL=15drops), calculate the number of mL of $\text{Na}_2\text{H}_2(\text{EDTA})$ needed for the titration of the soft and hard water samples.

4 Freezing point depression

EXPERIMENT

Freezing point depression

Goal

The goal of this laboratory is to put into practice the theory seen for *properties of solutions* and to experimentally see the *freezing point depression* effect. The added value is to play detective by relating the nature of *colligative properties* to the molar mass, to identify the unknown powder.

Materials

- | | |
|--|---|
| <input type="checkbox"/> A ring stand and two iron rings with a wire gauze | <input type="checkbox"/> 20mL pipet |
| <input type="checkbox"/> two or three jumbo test tubes with stoppers | <input type="checkbox"/> Ice |
| <input type="checkbox"/> A series of 50mL, 100mL, 250mL and 600mL beaker | <input type="checkbox"/> Perforated stoppers with stirring wire and thermometer |
| <input type="checkbox"/> Cyclohexane | |

Background

During winter, salt is poured over the street to melt the ice on the floor. General knowledge says that the salt melts the ice, while some venture to claim that the salt lowers the freezing temperature of the water. Both assertions might be in the right direction but are, if not false, at least incorrect. When salt is mixed with water, there is no water or salt anymore. Instead, there is a *solution* of salt in water. Since the freezing point of the solution is lower than that of the pure water, the solution is no longer solid at outdoor temperatures. Notice that there is no need to state the nature of the solution or the solute because it is irrelevant. The depression of the freezing point is not a property specific to the salt, (i.e. the sodium chloride, NaCl), but rather a general effect of any solute. Salt is used because it is widely available, dissolves fast, and is inexpensive.

Freezing point depression

The formula governing the change (Δ) in the freezing point (T_f) is given by:

$$\Delta T_f = - K_f m$$

where m is the concentration in terms of *molality* and K_f is the *molal freezing-point depression constant*. K_f is characteristic of the solvent being experimentally determined. Remember that molality m refers to the moles of solute concerning the kilograms of solvent. ΔT_f is the difference between the freezing point of the pure solvent (T_f^{solvent}) and that of the solution (T_f^{solution}). The solvent used in this experiment will be *cyclohexane*, a very *volatile* organic solvent. Special precautions must be taken to avoid evaporation of the solvent. K_f for cyclohexane is 20.1 °C/m. The solute will be an unknown compound. Hence you will not know its molar mass and therefore the molarity will not be explicitly known here.

Example

The freezing point of pure benzene is 5.50 °C. A solution is prepared by dissolving 0.450 g of an unknown substance in 27.3 g of benzene. The new freezing point is determined to be 4.18 °C. What is the molar mass of the unknown substance? The freezing point constant, K_f , for benzene is 5.12 °C/m.

Answer: calculate the concentration first:

$$\Delta T_f = K_f \times m \quad m = \frac{\Delta T_f}{K_f} \quad m = \frac{5.50^\circ\text{C} - 4.18^\circ\text{C}}{5.12^\circ\text{C} \frac{\text{kg}}{\text{mol}}} = 0.258 \frac{\text{mol}}{\text{kg}}$$

then, calculate the molar mass:

$$0.258 \frac{\text{mol}}{\text{kg}} \times 27.3\text{g} \times \frac{1\text{kg}}{1000\text{g}} = 7.04 \times 10^{-3} \text{mol} \quad \text{molar mass} = \frac{0.450\text{g}}{7.04 \times 10^{-3} \text{mol}} = 63.9 \frac{\text{g}}{\text{mol}}$$

The experiment

The experiment today is divided into the following parts: finding the freezing point of the pure solvent, finding the freezing point of the solution, finding the freezing point after increasing the concentration, and solving the molar mass of the solute.


Cooling plots

During a phase change, all the energy transferred is used to reorganize the molecules and none is used to change the temperature of the substance. A liquid at a higher temperature than its surroundings gives up the heat and its temperature falls until it reaches the freezing point. The liquid continues to give up heat to the surroundings, but its temperature remains constant. Only when all the liquid is frozen will the temperature start to fall again.

If one plots temperature vs. time for the freezing of the pure liquid, a negative slope line will represent the cooling down, while the flat horizontal line (plateau) corresponds to the freezing process. These types of plots are called cooling plots. Notice that the transition from the cooling line to the horizontal line is not sharp. For that reason, the best way to determine the freezing point is to draw a straight line fitting most of the cooling points and another straight line fitting most of the cooling points, while finding the crossing point between two lines. You can measure the freezing point as the temperature measurement of this crossing point. For a pure solvent, the freezing point should be near the temperature corresponding to the plateau. This is the graphical method to find freezing points. In the case of a solution, this constant-temperature plateau does not exist anymore. In contrast, when a solution freezes its temperature decreases with time. However, the same graphical method can be used to obtain the freezing point of a solution. Mixing the solution well is critical to obtaining a good cooling plot. If you do not mix the solution well you will obtain a series of steps when the solution freezes that will not allow you to measure an accurate freezing point.

Finding the freezing point of the pure solvent.

In this part of the experiment, you will obtain the cooling plot of the pure solvent, cyclohexane. Using the graphical method you will obtain the freezing point as an average between two replications. You must plot the results carefully to obtain accurate results.



Finding the freezing point of a solution.

In this part of the experiment, you will obtain the cooling plot of a solution made of cyclohexane and an unknown solute. Using the graphical method you will obtain the freezing point. You will work with two different concentrations and will carry


out one replication per concentration. This part aims to calculate the molar mass of the solute.

Example

Use the data reported below to graph a cooling plot, and use the graphical method to calculate the freezing point of the substance. Is this a pure substance or a solution?

t (min)	T (°C)	t (min)	T (°C)
1	65.6	9	40
2	61.6	10	40
3	56.7	11	40
4	52.2	12	40
5	47.8	13	40
6	43.3	14	40
7	38.9	15	40
8	40	16	40

Answer: Based on the plot below, we have that at early times the liquid is being cooled down and its temperature decreases with time. At late times, the liquid freezes and its temperature reaches a plateau. The graphical method to obtain freezing points consists on plotting a line for early times and another line for late times. These two lines intersect in a point, and the temperature for this point is the freezing point. As the temperature reaches a plateau, we are dealing here with a pure substance.




Procedure

Preparing the setup.

- ☐ *Step 1:* – Prepare the set-up to cool down the solution. Set up a ring stand and an iron ring with wire gauze on top. The height of the platform will normally be very low.
- ☐ *Step 2:* – Use a second, larger ring above the platform as protection. Attach a small clamp above the second iron ring.
- ☐ *Step 3:* – Put ice with some water in a 500 mL beaker and set it on the platform. This will be your cool bath. Water helps create a homogeneous bath.

Good Lab Practice

 Never add the liquid or solid to be weighed to the container while on the scale. Take the container *off* the scale prior to adding the product. A spill on the scale will ruin both your measurement and the scale.



Preparing a cyclohexane sample.

- ☐ *Step 1:* – Find 100mL and 50mL beakers and a 20 mL pipet.
- ☐ *Step 2:* – Obtain a large test tube that is clean and dry, a suitable rubber stopper, and a 100 mL beaker. The beaker will be used to hold the test tube vertically. Get the mass of the test tube with stopper and the beaker in the scale and write it down in the Results section.
- ☐ *Step 3:* – Pour no more than 25 mL of cyclohexane into the 50 mL beaker. Using the pipet, transfer exactly 20.0 mL of cyclohexane from the 50 mL beaker to the test tube and close it with the stopper immediately. Mind that as the liquid is volatile the measured volume will not be accurate. Using the scale get the mass of the test tube with the cyclohexane, the stopper, and the beaker. Make sure you do not weigh the thermometer or the stirring system.
- ☐ *Step 4:* – Get the mass of the filled-test tube with stopper and the beaker in the scale and write it down in the Results section.

Part A. Freezing point of the pure solvent.

- ☐ *Step 1:* – Insert the test tube into the ice and hold it vertically using the clamp. Make sure all liquid is fully submerged in the cool bath. Find the perforated stoppers with the stirring wire and the thermometer. The round wire end must encircle the thermometer probe. Exchange the stopper on the test tube with this gadget.
- ☐ *Step 2:* – As soon as the temperature is 15°C (or lower) start stirring the solvent with the wire and record the temperature every 15 seconds. Use all the digits given by the thermometer down to the tenth of the degree Celsius.
- ☐ *Step 3:* – Start plotting the results on the graph. Plotting temperature (in °C) in the vertical axis and time in the horizontal axis.
- ☐ *Step 4:* – Stop the experiment when the temperature starts dropping after several minutes of having a constant temperature. Do not stop until the temperature reaches at least 7°C.
- ☐ *Step 5:* – Now you will repeat the measurement one more time.
- ☐ *Step 6:* – Lift the test tube with the clamp above the ice and allow the cyclohexane to melt. Avoid using your hands to accelerate the process. Before the cyclohexane gets to 15°C get a second measurement by repeating the step above.

Good Lab Practice

- ✍ Stir, and do not shake. It is very important to stir the cyclohexane constantly and gently until it is completely frozen.
- ✍ Make sure the wire stays around the thermometer and that the latter is centered in the test tube.
- ✍ The thermometer should not touch the walls of the test tube.
- ✍ A bad stirring method will lead to a partial freezing and the graph will not show straight lines.



Part B. Freezing point of the less concentrated solution.

- ☐ *Step 1:* – Pick up an unknown solute. You have to carry out Parts B and C at the same time.
- ☐ *Step 2:* – The same cyclohexane can be used in Part B. However, if you work in teams you can prepare another cyclohexane sample and carry out work in parallel.
- ☐ *Step 3:* – If you use the sample from part A, lift the test tube with the clamp above the ice and allow the cyclohexane to melt. Avoid using your hands to accelerate the process. During the next steps be careful to avoid the cyclohexane's temperature from getting above 15°C.
- ☐ *Step 4:* – If you prepare a new cyclohexane sample, see the steps above.
- ☐ *Step 5:* – Tare the scale with the weighing boat and obtain between 0.10 g and 0.11 g of the unknown solute. Remove the boat with the solute, tare the empty scale, and record the mass of the weighing boat together with the solute. Use the balance of maximum precision.
- ☐ *Step 6:* – Carefully add the solute into the test tube with the pure solvent. Make sure no solute is spilled. There is no need to add all the solute from the weighing boat, the amount added will be calculated by difference.
- ☐ *Step 7:* – Weight the weighing boat with the solute leftovers. Record the mass on the results page.
- ☐ *Step 8:* – Mix the solution well and insert the test tube in the ice while holding it vertically using the clamp.
- ☐ *Step 9:* – Start stirring the solution with the wire and record the temperature every 15 seconds. This time the temperature should not reach constant. Stop the experiment when the temperature reaches 1°C.

Part C. Freezing point of the more concentrated solution.

- ☐ *Step 1:* – Lift the test tube with the clamp above the ice and allow the cyclohexane to melt. Avoid using your hands to accelerate the process. During the next steps be careful to avoid the cyclohexane getting above 15°C.
- ☐ *Step 2:* – The same boat can be reused without cleaning. Tare the scale with the boat inside. Obtain between 0.24 g and 0.25 g of the unknown solute and record the mass of the weighing boat together with the solute using maximum precision.
- ☐ *Step 3:* – Add carefully the solute into the test tube with a fresh new 20mL sample of the pure solvent. Make sure no solute is spilled. There is no need to add all the solute from the weighing boat, the amount added will be calculated by difference.
- ☐ *Step 4:* – Weight the weighing boat with the solute leftovers. Record the mass on the results page.
- ☐ *Step 5:* – Mix the solution well and insert the test tube in the ice while holding it vertically using the clamp.
- ☐ *Step 6:* – Start stirring the solution with the wire and record the temperature every 15 seconds. This time the temperature should not reach a constant. Stop the experiment when the temperature reaches 0.5°C.

Calculations

① This is the mass of the empty test tube with a cork and placed in a beaker.

② This is the mass of the test tube filled with cyclohexane, with a cork, and placed in a beaker.

③ This is the mass of solvent, m_{liquid} :

$$m_{solvent} = \textcircled{2} - \textcircled{1}$$

④ This is the volume of cyclohexane added, where d is the density of cyclohexane:

$$v_{solvent} = \frac{m_{liquid}}{d} = \frac{\textcircled{3}}{0.779}$$

⑤ These are two replicates for the freezing point of the solvent obtained through the graphical method.

⑥ This is the freezing point of the solvent.

⑦ This is the mass of the weighing boat with solute before adding the solute to the solvent.

⑧ This is the mass of the weighing boat with solute after adding the solute to the solvent.

⑨ This is the mass of solute:

$$m_{solute} = \textcircled{7} - \textcircled{8}$$

⑩ These are the freezing points of two solutions.

⑪ These are the values of the freezing points depression:

$$\Delta T_f = \textcircled{10} - \textcircled{6}$$

⑫ These is the molar mass of the solute:

$$\text{MW} = -\frac{K_f \cdot m_{solute}}{\Delta T_f \cdot m_{solvent}} \times 1000 = -\frac{20.1 \cdot \textcircled{9}}{\textcircled{11} \cdot \textcircled{3}} \times 1000$$

STUDENT INFO

Name:


Date:

Pre-lab Questions**Freezing point depression**

1. What does volatile mean?
2. What are the colligative properties?
3. What is molality? What are the units?
4. Why is molality used in this experiment instead of Molarity?
5. Use the data reported below to graph a cooling plot and use the graphical method to calculate the freezing point of the substance. Is this a pure substance or a solution?

t (min)	T (°C)	t (min)	T (°C)
1	55.6	9	30
2	51.1	10	28.6
3	46.7	11	27.1
4	42.2	12	25.7
5	37.8	13	24.3
6	33.3	14	22.86
7	28.9	15	21.6
8	31.4	16	20

(the question continues on the back)



6. What is the freezing point you calculated in the previous question?

STUDENT INFO

Name:

Date:

Results
EXPERIMENT

Freezing point depression

Parts A. Freezing point of the pure solvent.

1

Mass of the empty test tube with cork and beaker:

g

2

Mass of the filled test tube with cork and beaker:

g

3

Mass of the cyclohexane:

g

4

Volume of the cyclohexane used:

mL

Trial 1				Trial 2			
Time (sec)	Temp. (°C)	Time (sec)	Temp. (°C)	Time (sec)	Temp. (°C)	Time (sec)	Temp. (°C)

_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____

⑤ Freezing point of the pure cyclohexane. Trial 1: _____ °C Trial 2: _____ °C

⑥ Mean: _____ °C

Parts B and C. Freezing point of solutions.

- 1

Mass of the empty test tube with cork and beaker:

_____ g
- 2

Mass of the filled test tube with cork and beaker:

_____ g
- 3

Mass of the cyclohexane:

_____ g

Low concentration				High concentration			
Unknown #:							
7	Mass of boat:	_____ g			Mass of boat:	_____ g	
8	Mass of boat (after adding):	_____ g			Mass of boat (after adding):	_____ g	
9	Mass of solute:	_____ g			Mass of solute:	_____ g	
Time (sec)	Temp. (°C)	Time (sec)	Temp. (°C)	Time (sec)	Temp. (°C)	Time (sec)	Temp. (°C)
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____

_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____	_____	_____

Low concentration

High concentration

⑩ Freezing point: _____ °C

Freezing point: _____ °C

⑪ ΔT_f : _____ °C ΔT_f : _____ °C

⑫ Molar mass: _____ g/mol

Molar mass: _____ g/mol

Mean molar mass: _____ g/mol




Figure 1: Temperature (Y axis) vs. time (X axis)




Figure 2: Temperature (Y axis) vs. time (X axis)




Figure 3: Temperature (Y axis) vs. time (X axis)




Figure 4: Temperature (Y axis) vs. time (X axis)

STUDENT INFO

Name:

Date:

Post-lab Questions

Freezing point depression

1. If the empirical formula of the solute used in this experiment is $\text{C}_3\text{H}_2\text{Cl}$, what is the molecular formula?
2. What would be the effect on the molar mass in the following situations? Explain why.
 - (a) A decrease in the mass measured for the solute.
 - (b) An increase in the mass measured for the solute.
 - (c) Dirt (soluble) from the stopper was added to the solution.
 - (d) A loss of cyclohexane.
3. How much would the molar mass change if 1.000 g of the cyclohexane evaporates during part B? What would be the induced error?

5 The iodine clock

EXPERIMENT

The iodine clock

Goal

The goal of this laboratory is to use the **method of initial rates** to calculate the rate of a chemical reaction. On top of that, the student will also visualize experimentally the effect of **temperature** and of a **catalyst** on the rate of a chemical reaction.

Background

The rate law expresses the progress of a chemical reaction in function of the concentration of the reactant

$$\text{Rate} = k[\text{A}]^a[\text{B}]^b,$$

where k is the proportionality constant called the **rate constant**, $[\text{A}]$ and $[\text{B}]$ are the respective concentrations for reactants A and B, and the exponents are the **orders of the reactants**, a way to establish their corresponding relevance in the reaction rate.

Initial rates method


The initial rates method is used to measure the rate for a chemical reaction and to find out the rate constant k and the order of the reactions. This method is based on the assumption that at the beginning of a chemical reaction, there is no reverse reaction, i.e. forming reactants from the products, can take place, and therefore the rate is solely affected by the reactants concentrations. The time for the reactant concentrations to change by the same amount is measured. Several experiments are carried out using all the same starting concentrations but for one reactant. The beauty of this experiment is the method used to make the reactant's concentration change by the same amount for all the experiments. This is probably the most often misunderstood part of this laboratory and it is key to understand it thoroughly.

The Iodine Clock

A watch system is a method used in the Navy to schedule regular periods aboard. A ship's bell is used to indicate the end of each period, allowing the crew to know when it is time to switch their task, take breaks, or have lunch.

A **chemical clock** is a chemical reaction with a sudden change in one of its properties, the tempo of which can be adjusted by adjusting the concentrations. The **iodine clock** is a chemical clock marked by a sudden, almost magical, color change. For the present experiment, the iodine clock will be connected to the reaction being studied and will indicate the moment our reaction reaches the change in concentration desired. It is very important to understand the mechanism of the two competing reactions involved in this experiment to fully understand this activity.


The reaction setting the time is:



This reaction consumes iodine, I_2 . Iodine reacts with starch turning the solution dark. The starch will be used in this experiment in a way similar to the indicator used in titration.

starch with iodine

starch, no I_2




black-dark blue

transparent


The iodine clock reaction will be running in parallel with another reaction that produces iodine. In the presence of starch, iodine becomes dark blue, almost black. Therefore, as long as there is thiosulfate ($\text{S}_2\text{O}_3^{2-}$) in the beaker the iodine will react with the latter and keep the solution transparent. We could say that the timer is adjusted by the concentration of thiosulfate ($\text{S}_2\text{O}_3^{2-}$).

The reaction to be studied

In this experiment, the student will determine the reaction rate of the reduction of potassium persulfate ($\text{K}_2\text{S}_2\text{O}_8$) with sodium iodide (NaI).



As ionic compounds dissociate in aqueous solutions it is more accurate to write the **net ionic equation**. These are the ionic net equations for the two reactions, the reduction of potassium persulfate with sodium iodide and the iodine clock:



To summarize: The first reaction will be running producing I_2 , which would turn dark blue in the presence of starch. However, the second reaction will avoid the color change by reducing the iodine into iodide. The solution will remain transparent as long as there is thiosulfate ($\text{S}_2\text{O}_3^{2-}$) available. Once the thiosulfate is finished, iodine will finally remain in the solution turning it dark blue. Since the amount of thiosulfate will be constant for all the experiments to be performed, the amount of iodine produced in the first reaction, and therefore the reactant's concentration change, will be the same each time the solution turns dark. The iodine clock is being used to schedule a constant change in the reactant's concentrations.

The experiment


The experiment today is divided into the following parts: obtaining the reaction rate constant, evaluating the influence of the temperature on the reaction rate, and evaluating the influence of a catalyst on the reaction rate.

Part A. Obtaining the reaction rate constant.

The reaction to be studied and its rate law are:





Using the **initial rates method** the reactants orders, a and b , can be solved. Find the suitable combinations of experiments where only one initial concentration of $\text{S}_2\text{O}_8^{2-}$ or I^- changes and compare the rates using a fraction. The initial concentrations must be recalculated using the volumes and their original concentrations. Once the reactant orders are known and the initial concentrations of $\text{S}_2\text{O}_8^{2-}$ and I^- for each reaction have been recalculated, the reaction constant, k , can be solved. See the example for more details:

	Volume					
	NaI	NaCl	Na ₂ S ₂ O ₃	Starch	K ₂ SO ₄	K ₂ S ₂ O ₈
	0.2M	0.2M	0.01M	2%	0.2M	0.2M
Exp. 1	2.0	2.0	2.0	1.0	2.0	2.0
Exp. 2	2.0	2.0	2.0	1.0	0	4.0
Exp. 3	4.0	0	2.0	1.0	2.0	2.0

Example

The reaction starts from a solution containing A (1 M) and B (0.5 M) and D (1 M). There are two reactions occurring in the solution. A slow reaction where C is produced and a fast reaction where C is reduced to E ions by D:



An indicator is added, such that the solution will turn dark in the presence of substance C.

Three experiments are performed using different initial volumes and the following times for the solution to turn dark are measured;

	Volumes (mL)				time
	A	B	D	Water	(s)
Exp. 1	4.0	2.0	2.0	2.0	20
Exp. 2	4.0	1.0	2.0	2.0	40
Exp. 3	2.0	1.0	2.0	4.0	160

Calculate the reactant orders and the reaction rate constant.

Answer: First write the reaction rate expression: $\text{Rate} = k[A]^a[B]^b$

Recalculate the concentrations for each reactant. Use the formula for dilutions $M_1 V_1 = M_2 V_2$, where V_2 is the total volume $V_A + V_B + V_D + V_{\text{water}}$.

For example, to calculate the concentration of A in Experiment 1:

$$M_1 V_1 = M_2 V_2 \quad M_2 = \frac{M_1 V_1}{V_2} \quad M_2 = \frac{1M \cdot 4mL}{10mL} = 0.4M$$

	[A]	[B]	[D]	rate (Ms ⁻¹)
Exp. 1	0.4 M	0.2 M	0.2 M	$5 \cdot 10^{-3}$
Exp. 2	0.4 M	0.1 M	0.2 M	$2.5 \cdot 10^{-3}$
Exp. 3	0.2 M	0.1 M	0.2 M	$6.25 \cdot 10^{-4}$

The reaction rate is the change in concentration per unit time. The change in concentration is determined by the second reaction (the clock). For the first experiment it takes 20 seconds to turn black. In those 20 seconds all of D has reacted, i.e. 0.2 M in 20 seconds.

Following the stoichiometry of the reaction, if 1 mole of D reacts with 1 mole of C and every 2 moles of C are produced when 1 mole of A reacts. The relation between D and A is 2:1.

$$1 \text{ mole of A} = 1 \text{ mole A} \frac{2 \text{ moles C}}{1 \text{ moles A}} \frac{1 \text{ mole D}}{1 \text{ mole C}} = 2 \text{ moles of D}$$

Therefore, 0.2 M of D corresponds to 0.1 M of A reacting in 20 s. The rate for experiment 1 is $5 \cdot 10^{-3}$ M/s. The order of the reactants are calculated by comparing the rate laws. Each time, choose two experiments where only one parameter is different.

$$\frac{\text{Rate 2}}{\text{Rate 3}} = \frac{k[A]^a[B]^b}{k[A]^a[B]^b} = \frac{k \cdot 0.4^a \cdot 0.1^b}{k \cdot 0.2^a \cdot 0.1^b} = \frac{2.5 \cdot 10^{-3}}{6.25 \cdot 10^{-4}}$$

$$2^a = 4 \quad a = 2$$

Similarly, using rate 2 over rate 1, $b=1$ is obtained.

The reactant order can be used to obtain the reaction rate.

$$\text{Rate} = k[A]^a[B]^b \quad k = \frac{\text{Rate}}{[A]^a[B]^b} = \frac{5 \cdot 10^{-3} \text{ M s}^{-1}}{(0.4 \text{ M})^2 (0.1 \text{ M})^1} = 0.15625 \text{ M}^{-2} \text{ s}^{-1}$$

Since experimentally it is more likely to obtain one slightly different k for each experiment, an average should be calculated.

Parts B. Evaluating the influence of the temperature on the reaction rate.

According to **Arrhenius equation**, the rate constant is affected by the temperature following the equation:

$$k = Ae^{-E_a/RT}$$

where A is a constant, E_a is the activation energy, R is the universal gas constant, and T is the temperature. A higher temperature would make the negative exponent smaller, and therefore the k would be greater. In other words, an increase in the temperature should accelerate the reaction. It is possible to calculate the activation energy by measuring the reaction rates at two different temperatures. After calculating the respective rate constants these are compared in their logarithmic forms:

$$\ln(k) = -\frac{E_a}{RT} + \ln(A)$$

After subtracting the expression above for both temperatures, T_1 and T_2 , the following equation is generated:

$$\ln\left(\frac{k_2}{k_1}\right) = \frac{E_a}{R} \left(\frac{1}{T_1} - \frac{1}{T_2}\right)$$

The values of k_1 and k_2 obtained respectively at T_1 and T_2 are used to calculate E_a . The temperatures should be given in Kelvins since $R = 8.3145 \text{ J/K mol}$.

Part C. Evaluating the influence of a catalyst on the reaction rate.

A catalyst is a substance that without being consumed, accelerates the reaction. The catalyst generates a new reaction path with a lower activation energy. Looking at the Arrhenius equation, the effect of lowering the activation energy works in the same way as an increase in temperature. The catalyst used in this experiment will be copper sulfate (CuSO_4). The Cu(II) cations are responsible for the catalytic effect.

Procedure

Part A. Obtaining the rate law. Preparing the mixtures

- ☐ *Step 1:* – Obtain three large test tubes with their corresponding stoppers, a suitable grid or rack, a Mohr pipet, and a thermometer.
- ☐ *Step 2:* – Label the test tubes with numbers 1, 2 and 3.
- ☐ *Step 3:* – Read the Good Lab Practice box

Good Lab Practice

- ✍ Pipettes are always used together with a suction bulb or a syringe. **Never** such a chemical using the mouth!
- ✍ **Never** leave a pipette laying on the table **while connected to the syringe**. Remaining liquid can leak into the syringe. The syringe will be damaged and the liquid will return contaminated to the pipette.
- ✍ Always wash pipettes with distilled water and rinse them using the same product to be measured. Such a small volume, disconnect the syringe and move the liquid about displacing the distilled water. Discard the rinsing liquid as waste. Rinse twice.

- ☐ **Step 4:** – Obtain about 25 mL **0.20 M NaI** in a small beaker.
- ☐ **Step 5:** – Clean and rinse the pipette with NaI. Pipette 2.0 mL in test tube 1, 2.0 mL in test tube 2 and 4.0 mL in test tube 3. Refer to the table below (Table of quantities) for the reagent quantities to add.
- ☐ **Step 6:** – Obtain no more than 20 mL **0.20 M NaCl** in a small beaker.
- ☐ **Step 7:** – Clean and rinse the pipette with NaCl. Pipette 2.0 mL in test tube 1 and 2.0 mL in test tube 2. Nothing in test tube 3. Refer to the table below (Table of quantities) for the reagent quantities to add.
- ☐ **Step 8:** – Obtain no more than 20 mL **0.010 M Na₂S₂O₃** in a small beaker.
- ☐ **Step 9:** – Clean and rinse the pipette with Na₂S₂O₃. Pipette 2.0 mL in each test tube. Refer to the table below (Table of quantities) for the reagent quantities to add.
- ☐ **Step 10:** – Obtain approximately 10 mL **2% starch solution** in a small beaker.
- ☐ **Step 11:** – Clean and rinse the pipette with starch. Pipette 1.0 mL in each test tube. Refer to the table below (Table of quantities) for the reagent quantities to add.
- ☐ **Step 12:** – Obtain no more than 20 mL **0.20 M K₂SO₄** in a small beaker.
- ☐ **Step 13:** – Clean and rinse the pipette with K₂SO₄. Pipette 2.0 mL in test tube 1, nothing in test tube 2 and 2.0 mL in test tube 3. Refer to the table below (Table of quantities) for the reagent quantities to add.
- ☐ **Step 14:** – Check that the volumes in the 3 test tubes make sense with the following table. Test tubes 1 and 3 should have the same volume. Test tube 2 should have slightly less volume.
- ☐ **Step 15:** – Using the thermometer, record the room temperature in Kelvins.

Table of quantities						
	Volume (mL)					⌚
	NaI	NaCl	Na ₂ S ₂ O ₃	Starch	K ₂ SO ₄	K ₂ S ₂ O ₈
	0.2M	0.2M	0.01M	2%	0.2M	0.2M
Exp. 1	2.0	2.0	2.0	1.0	2.0	2.0
Exp. 2	2.0	2.0	2.0	1.0	0	4.0
Exp. 3	4.0	0	2.0	1.0	2.0	2.0

Part A. Obtaining the rate law. Mixing the test tubes

- ☐ *Step 1:* – Now you are about to start up the reaction. It is critical to follow the steps below.
- ☐ *Step 2:* – Every reaction should be done in **the same way** to be comparable. The same steps, in the same order and at the same pace should be followed.
- ☐ *Step 3:* – You will be adding $K_2S_2O_8$. After adding the $K_2S_2O_8$, the test tube should be closed with the stopper, and the content should be mixed by inverting the test tube 16 times. Do not rush, turn the test tube *Mississippily*. Do the same for the following reactions.

⚠ CAUTION!

- ⚠ The stopper should always be placed before inverting the test tube.
- ⚠ Failure to place the stopper will cause the spill of the reaction.
- ⚠ There is a story behind every Caution Box.

Part A. Obtaining the rate law. Starting the reaction

- ☐ *Step 1:* – Obtain about 25 mL **0.20 M $K_2S_2O_8$** in a small beaker. Make sure you do not contaminate the bottle of reagent by always retrying out of the bottle into a beaker and never pipetting directly from the bottle. Clean and rinse the pipette with $K_2S_2O_8$.
- ☐ *Step 2:* – Get test tube 1 and a stopper. Set the time to zero.
- ☐ *Step 3:* – Start the stopwatch as you pipette 2.0 mL of $K_2S_2O_8$ in test tube 1. Time should start to be measured as soon as you add the first drop of $K_2S_2O_8$ in a test tube. Mix the reaction as explained above.
- ☐ *Step 4:* – Stop the stopwatch when the solution turns dark. The color will appear suddenly. If more than 5 minutes elapse without the solution turning dark, start doubting about what have you mixed.
- ☐ *Step 5:* – Repeat the previous steps now for test tube 2 and pipetting this time 4.0 mL of $K_2S_2O_8$ into the test tube. .
- ☐ *Step 6:* – Repeat one more time the previous steps now for test tube 3 and pipetting this time 2.0 mL of $K_2S_2O_8$ into the test tube.
- ☐ *Step 7:* – Now you need to repeat all steps above (mixture preparation and reaction) to do a second trial for the three test tubes.
- ☐ *Step 8:* – Discard the solutions in the indicated waste container and clean the test tubes.

Part B. Evaluating the influence of the temperature on the reaction rate.

- ☐ *Step 1:* – Label three large test tubes with the letters “T”, “C” and “K”. “T” stands for “Temperature”, “C” for “Catalyst” and “K” for “ $K_2S_2O_8$ ”.
- ☐ *Step 2:* – Add to test tubes “T” and “C” the same volumes of NaI, NaCl, $Na_2S_2O_3$, Starch and K_2SO_4 as for test tube 1. Do not add the $K_2S_2O_8$ yet. Reserve test tube “C” in the rack for Part C.
- ☐ *Step 3:* – Add about 5 mL of $K_2S_2O_8$ to test tube “K”.
- ☐ *Step 4:* – Obtain a 600 mL beaker and, with the help of the thermometer, add hot and cold water until you have about 400 mL of water between 40 and 45 degrees Celsius.

- ☐ *Step 5:* – Place the test tubes “T” and “K” in the hot water bath and allow them to temper for at least 5 minutes. Part C can be done during those 5 minutes. .
- ☐ *Step 6:* – Measure the temperature of the water after the tempering time. Write down this temperature in the results table. Start the reaction by pipetting 2 mL of $\text{K}_2\text{S}_2\text{O}_8$ from test tube “K” into test tube “T”. After the 16-Mississippi mixing, return the test tube to the bath and wait to record the time it takes for the solution to turn dark.

Part C. Evaluating the influence of a catalyst on the reaction rate.

- ☐ *Step 1:* – Add **one drop** of **0.2 M CuSO_4** to test tube “C”. Shake the test tube to mix the solution.
- ☐ *Step 2:* – Start the stopwatch as you pipette 2.0 mL of $\text{K}_2\text{S}_2\text{O}_8$ in the test tube. Mix the reaction as explained above.
- ☐ *Step 3:* – Record the time it takes for the solution to turn dark.

Calculations

① This is lab temperature. You will use this value for Part B as well.

① This is the volume of iodine added in each experiment.

② This is the volume of thiosulfate added in each experiment.

③ This is the volume of peroxydisulfate added in each experiment.

④ This is the total volume v_{total} of the reacting mixture.

⑤ This is the peroxy disulfate concentration in the mixture:

$$[\text{S}_2\text{O}_8^{-2}] = \frac{v_{\text{S}_2\text{O}_8^{-2}} \cdot c_{\text{S}_2\text{O}_8^{-2}}}{v_{total}} = \frac{\textcircled{3} \cdot 0.2}{\textcircled{4}}$$

⑥ This is the iodine concentration in the mixture:

$$[\text{I}^-] = \frac{v_{\text{I}^-} \cdot c_{\text{I}^-}}{v_{total}} = \frac{\textcircled{1} \cdot 0.2}{\textcircled{4}}$$

⑦ This is the change in thiosulfate concentration during the reaction:

$$\Delta[\text{S}_2\text{O}_3^{-2}] = \frac{v_{\text{S}_2\text{O}_3^{-2}} \cdot c_{\text{S}_2\text{O}_3^{-2}}}{v_{total}} = \frac{\textcircled{2} \cdot 0.01}{\textcircled{4}}$$

⑧ This is the average time Δt for each experiment.

⑨ This is the reaction rate:

$$r = \frac{\Delta[\text{S}_2\text{O}_3^{-2}]}{2 \cdot \Delta t} = \frac{\textcircled{7}}{2 \cdot \textcircled{8}}$$

⑩ This is the peroxy disulfate order in the rate law:

$$a = \frac{\ln\left(\frac{r_{Exp1}}{r_{Exp2}}\right)}{\ln\left(\frac{[S_2O_8^{-2}]_{Exp1}}{[S_2O_8^{-2}]_{Exp2}}\right)} = \frac{\ln\left(\frac{(9)_{Exp1}}{(9)_{Exp2}}\right)}{\ln\left(\frac{(5)_{Exp1}}{(5)_{Exp2}}\right)}$$

- (11) This is the iodide order in the rate law:

$$b = \frac{\ln\left(\frac{r_{Exp1}}{r_{Exp3}}\right)}{\ln\left(\frac{[I^-]_{Exp1}}{[I^-]_{Exp3}}\right)} = \frac{\ln\left(\frac{(9)_{Exp1}}{(9)_{Exp3}}\right)}{\ln\left(\frac{(6)_{Exp1}}{(6)_{Exp3}}\right)}$$

- (12) These are the rate constants calculated for each experiment:

$$k = \frac{r}{[S_2O_8^{-2}] \cdot [I^-]} = \frac{(9)}{(5) \cdot (6)}$$

- (13) This is the average rate constant.

- (14) The data on this row should be the same as in the first row in the Results table of Part B.

- (15) This is the activation energy. You need to use the two different reaction constants at different temperatures. R is 3.184 J/K·mol:

$$E_a = \frac{T_1 \cdot T_2}{T_2 - T_1} \cdot R \cdot \ln\left(\frac{k_2}{k_1}\right)$$

STUDENT INFO

Name:

Date:

Pre-lab Questions

The iodine clock

1. What is the difference between a test tube and a jumbo test tube?.
2. For the reaction $2A + 3B \rightarrow C$ what is the mathematical expression that relates the rate of A (r_a) and the rate of B (r_b)?
3. Explain the Iodine clock reaction used in this experiment.
4. What is starch and where can you find it in nature?
5. Low-quality boiled ham often contains potatoes to make it cheaper. How could you test if one piece of ham is pure meat or it contains potatoes?
6. Find out the difference between a **Mohr pipette** and a **transfer pipette**.

Results

Date:

K

K

Trial	Exp.	C_{NaI}	$C_{Na_2S_2O_3}$	C_{starch}	$C_{K_2SO_4}$	$C_{K_2S_2O_8}$	v_{Total}	time	$[S_2O_8^{-2}]$	$[I^-]$	$\Delta[S_2O_3^{-2}]$
		0.2M	0.01M	2%	0.2M	0.2M	(mL)	(s)	(M)	(M)	(M)
		①	②		③	④			⑤	⑥	⑦
1	1	2mL	2mL	1mL	2mL	2mL	11mL				
1	2	2mL	2mL	1mL	0mL	4mL	11mL				
1	3	4mL	2mL	1mL	2mL	2mL	11mL				
2	1	2mL	2mL	1mL	2mL	2mL	11mL				
2	2	2mL	2mL	1mL	0mL	4mL	11mL				
2	3	4mL	2mL	1mL	2mL	2mL	11mL				
⑧ Average time $\Delta t(s)$											

Part A. Obtaining the rate law.

Exp.	$\Delta[\text{S}_2\text{O}_3^{2-}]$ (M)	Δt (s)	r (M/s)	$[\text{S}_2\text{O}_8^{2-}]$ (M)	$[\text{I}^-]$ (M)	k
	(7)	(8)	(9)	(5)	(6)	(12)

(14)	1	_____	_____	_____	_____	_____	_____
	2	_____	_____	_____	_____	_____	_____
	3	_____	_____	_____	_____	_____	_____

$$r = k [\text{S}_2\text{O}_8^{2-}]^a [\text{I}^-]^b$$

(10) a= _____

(11) b= _____

(13) Average k (write down units)= _____

Parts B and C. Evaluating the influence of the temperature on the reaction rate and a catalyst

Exp.	Temperature (K)	$\Delta[\text{S}_2\text{O}_3^{2-}]$ (M)	Δt (s)	r (M/s)	$[\text{S}_2\text{O}_8^{2-}]$ (M)	$[\text{I}^-]$ (M)	k
		(7)	(8)	(9)	(5)	(6)	(12)

(14)	1	_____	_____	_____	_____	_____	_____
	1	_____	_____	_____	_____	_____	_____

(Results on the first row of Table A and Table B should be the same)

Catalyst 1 _____

(15) E_a = _____

J/K·mol

STUDENT INFO

Name:

Date:


Post-lab Questions

The iodine clock

1. Explain the effect of adding a catalyst to the reaction rate.
2. Using collision theory, explain how the concentration of the reactants affect the rate of the reaction.


Le Châtelier's Principle

In the case of Hexaaminenickel (II), you will also study the change in equilibrium position when adding an acid. The ligand will act as a base reacting with the acid and detaching from the metal.



Part C. Effect of pH on an indicator.


pH indicators are substances that change color depending on the medium pH. Typically they are weak bases or acids. As such, they dissociate slightly in water forming ions. Those ions might be colored species. The general expression for the dissociation of a divalent weak acid indicator (H_2In) is:



In the case of Methyl orange, the protonated form is red, while the conjugated base is yellow. Notice that hydronium is a product in this equilibrium. How will adding an acid affect this equilibrium?


Part D. Effect of pH on solubility.

Hydroxide salts of Group II elements (Ca, Sr, and Ba) are slightly soluble. Adding acid will neutralize the hydroxides in the solution while adding more hydroxide will increase the ion product. In both cases, the amount of precipitate will be affected.



Part E. Heat as a product

Lastly, you will experience Le Châtelier's Principle in an exothermic reaction, where the heat can be viewed as a product of the reaction.



Example

The following endothermic reaction is allowed to reach equilibrium:



Where A, the only colored compound, is blue. How will the following changes affect the color/equilibrium: (a) Adding more B reactant, (b) Adding more C product, (c) Heating the mixture.



Answer: (a) less blue/equilibrium shifts to the right. (b) more blue/ equilibrium shifts left. (c) Heat is a reactant; less blue/equilibrium shifts to the right.

Procedure

Part A. Colored complexed ions: FeSCN^{2+}

- ☐ *Step 1:* – Find 3 test tubes and a 100 mL beaker, clean them, and mark the test tubes with letters A, B, and R.
- ☐ *Step 2:* – Read the Good Lab Practice box

Good Lab Practice

-  Concentrated solutions of ammonia, sodium hydroxide and hydrochloric acid should be handled with care.
-  Do not shake open test tubes. Do not use your fingers as stoppers, even if you wear gloves.

- ☐ *Step 3:* – Add about 20 mL of distilled water to the beaker, 20 drops of 0.1 M $\text{Fe}(\text{NO}_3)_3$ and 20 drops of 0.1 M KSCN. Mix the solution until the color is homogeneous.
- ☐ *Step 4:* – Use a 10 mL graduated cylinder to add 3 mL of the mixture to each of the test tubes.
- ☐ *Step 5:* – Add 20 drops of 0.1 M $\text{Fe}(\text{NO}_3)_3$ to test tube A. Put a stopper and mix the solution.
- ☐ *Step 6:* – Add 20 drops of 0.1 M KSCN to test tube B. Put a stopper and mix the solution.
- ☐ *Step 7:* – Add 20 drops of distilled water to test tube R. Put a stopper and mix the solution.
- ☐ *Step 8:* – Compare the color of the test tubes A and B to R and write down your observations.

Part B. Colored complexed ions: $\text{Ni}(\text{NH}_3)_6^{2+}$.

- ☐ *Step 1:* – Obtain 1 test tube and clean it.
- ☐ *Step 2:* – Add 10 drops of 0.1 M $\text{Ni}(\text{NO}_3)_2$. Indicate the color on the results page.
- ☐ *Step 3:* – Add drops of 6 M NH_3 until the color changes.
- ☐ *Step 4:* – Add drops of 6 M HCl until the color changes.

Part C. Effect of PH on an indicator.

- ☐ *Step 1:* – Obtain two 50 mL beakers, clean them, and mark them with letters A (for Acid) and B (for Base). Find and clean a test tube.
- ☐ *Step 2:* – Add 10 mL of distilled water to each beaker and 1 mL of distilled water to a test tube.
- ☐ *Step 3:* – Add 4 drops of 6 M HCl to beaker A and stir it (this is the diluted acid).
- ☐ *Step 4:* – Add 4 drops of 6 M NH_3 to beaker B and stir it (this is the diluted base).
- ☐ *Step 5:* – Add 4 drops of the indicator (methyl orange) to the test tube.
- ☐ *Step 6:* – Add 2 drops from the diluted acid in beaker A to the test tube. Mix gently the solution. Record the color on the results page.

- ☐ *Step 7:* – Drop by drop, add the diluted ammonia solution from beaker B until the color changes. Homogenize the solution from time to time. Write down the number of drops.
- ☐ *Step 8:* – Attempt one more time for another color change by adding drop-by-drop diluted acid solution. Write down the number of drops.

Part D. Effect of pH on solubility.

- ☐ *Step 1:* – Add 5 mL of 6 M NaOH to a 50mL beaker using the graduated cylinder.
- ☐ *Step 2:* – Rinse a 10mL cylinder with water 3 times. Then, use it to add 5 mL of 1 M $\text{Ca}(\text{NO}_3)_2$ to the beaker.
- ☐ *Step 3:* – Stir the mixture.
- ☐ *Step 4:* – Make a cone with filter paper and place it in a funnel, on top of the Erlenmeyer. Filter the solution with the precipitate. Carry out a couple of washings with distilled water to make sure you collect all solids. Transfer with a spatula the white solid to a small clean beaker.
- ☐ *Step 5:* – Add 10 mL of distilled water to the beaker with the white solid. Stir the mixture. Do not expect the solid to completely dissolve. The solution is saturated.
- ☐ *Step 6:* – Add drops of 6 M HCl until a change is observed. Record your results.
- ☐ *Step 7:* – Add drops of 6 M NaOH until a change is observed. Record your results.

Part E. Heat as a product.

- ☐ *Step 1:* – Obtain a Bunsen burner, a stand, two iron rings, a wire gauze, two 250 mL beakers, and a test tube.
- ☐ *Step 2:* – Prepare the setup for the bunsen burner, using the second iron ring to protect the beaker from falling. Put about 100 mL of distilled water in the beaker and bring the water to a boil.
- ☐ *Step 3:* – Add 5 drops of 0.1 M $\text{Co}(\text{NO}_3)_2$ to the test tube.
- ☐ *Step 4:* – Record the color of the liquid solution after each step.
- ☐ *Step 5:* – Add drops of 12 M HCl until the color of the solution changes. You might need to stir the test tube to help mix the reagents.
- ☐ *Step 6:* – Add 5 drops of distilled water and mix.
- ☐ *Step 7:* – Place the test tube in the boiling water and wait for another color change.
- ☐ *Step 8:* – Attempt to reverse the reaction by placing the test tube in a beaker with ice or cold water.

STUDENT INFO

Name:

Date:

Pre-lab Questions

Le Châtelier's Principle

1. Write the equilibrium constant expressions for all reactions involved in this experiment.
2. How does pH affect the solubility of Ca(OH)_2 ?
3. Is heat being consumed (reactants) or produced (product) in an endothermic reaction?
4. Is heat being consumed (reactants) or produced (product) in an exothermic reaction?


STUDENT INFO

Name:


Date:

Results
EXPERIMENT


Le Châtelier's Principle

Part A. Colored complexed ions: FeSCN^{2+} 


After step #			Before color	After color
5	Add Fe^{3+}	Solution in test tube A	_____	_____
6	Add SCN^{-}	Solution in test tube B	_____	_____
7	Add H_2O	Solution in test tube R	_____	_____

Part B. Colored complexed ions: $\text{Ni}(\text{NH}_3)_6^{2+}$.


After step #	Color		
2	Add Ni^{2+}	_____	
3	Add NH_3	_____	# drops added _____
4	Add HCl	_____	# drops added _____

Part C. Effect of pH on an indicator.

After step #		Before color	After color	
5	Add indicator	_____	_____	
6	Add diluted acid	_____	_____	
7	Add diluted base	_____	_____	# drops added ____
8	Add diluted acid	_____	_____	# drops added ____

Part D. Effect of pH on solubility.

After step #		Indicate # drops added
6	Add acid	_____
7	Add base	_____

Part E. Heat as a product.

After step #		color	
3	Add Co^{2+}	_____	
5	Add acid	_____	# drops added ____
6	Add H_2O	_____	
7	Heat	_____	
8	Cool	_____	

STUDENT INFO

Name:

Date:

Post-lab Questions

Le Châtelier's Principle

1. In part A. Explain the different colors based on the equilibrium reaction and Le Châtelier's Principle.
2. In part A. Given that water is not involved in the equilibrium, explain the color change you observed after adding water to the mixture.
3. In part B. What did you observe after adding the acid and after adding the base? Explain the different colors based on the equilibrium reaction and Le Châtelier's Principle.
4. In Part C. What did you observe after adding the acid and after adding the base? Explain the changes based on the equilibrium reaction and Le Châtelier's Principle.
5. In Part D. What did you observe after adding the acid and after adding the base? Explain the changes based on the equilibrium reaction and Le Châtelier's Principle?

6. In Part E. What did you observe after heating and after cooling? Explain the changes based on the equilibrium reaction and Le Châtelier's Principle.

7 Equilibrium constant

EXPERIMENT


Equilibrium constant

Goal

The goal of this laboratory is to determine *the equilibrium constant* of a chemical reaction by using the *Lambert-Beer's law* to measure the concentration of the species in solution.

Background

Most chemical reactions do not proceed to completion, proceeding only to the point where both reactants and products have constant concentration. This is because most reactions are reversible and they can run in both forward and reverse directions. The concentrations of the species involved in the equilibrium will achieve an *equilibrium state* when the reaction rates in both directions equalize. For a certain chemical reaction:



The equilibrium can be characterized by the expression:

$$K_c = \frac{[C]^c [D]^d}{[A]^a [B]^b}$$

where K_c (capital letter) is the equilibrium constant in terms of molarity, the brackets represent the concentrations of the respective species and the exponents are their corresponding stoichiometric coefficients.

Spectrophotometry

Spectrophotometry is a technique that measures the amount of light absorbed by a chemical substance, typically in solution, using the use of a *spectrophotometer*. Spectrophotometers display the absorbance, the amount of light absorbed by a sample. To properly measure the absorbance of a chemical in a water-based solution you need to set the blank by using a sample containing only the solvent without a solute. The blank establishes the baseline of the measurement, eliminating any absorbance caused by the solvent.

The spectrometer allows the measurement of absorbance at a single wavelength. Interestingly, for each substance, there is a wavelength of *maximum absorbance*, which allows us to differentiate among chemical substances with different maximum absorbances. For example, in the current experiment, you will analyze mixtures containing $\text{Fe}_{(\text{aq})}^{3+}$, $\text{SCN}_{(\text{aq})}^-$, and $\text{Fe}(\text{SCN})_{(\text{aq})}^{2+}$. However, only $\text{Fe}(\text{SCN})_{(\text{aq})}^{2+}$ will absorb radiation at $\lambda=450$ nm, and hence we can track the absorbance of this chemical only in the mixture.

Lambert-Beer's law

Liquids *attenuated* the intensity of light passing through a colored solution. The effect can be comparable to a dirty window that attenuates the light passing through. The dirtier the window (higher concentration of dirt), the less one can see through it (because more light is absorbed, reflected, or refracted by the dirt). Lambert-Beer's law establishes the relationship between absorbance and concentration:

$$A = k_{\lambda} c$$

In this formula, A stands for the *absorbance* (light absorbed by the liquid), k_{λ} (lower case, not to be confused with the equilibrium constant, K) is a **proportionality constant**, which is specific to each compound and at the wavelength (λ) specified by the subindex, and c is the solution concentration. Lambert-Beer's law is a linear function of the form $y = mx + b$.


When graphing the absorbance (y -axes) versus the concentration (x -axes), the slope of the line corresponds to Lambert-Beer's proportionality constant. This line should cross the origin as there should be no absorbance when the concentration is null.

The experiment

The experiment today is divided into 2 parts: determining Lambert-Beer's constant k and determining the equilibrium constant K .

Part A. Determining Lambert-Beer's constant k .

This experiment studies the reaction between iron(III) and thiocyanate to give an iron-thiocyanate complex:



Interestingly, Fe^{3+} in solution is yellow and becomes red when reacting with thiocyanate to produce the $(\text{FeSCN})^{2+}$ complex. This intense red coloration allows us to use spectrophotometry to measure the $(\text{FeSCN})^{2+}$ concentration and hence, of all the species in equilibrium.

The first part of this experiment requires the determination of Lambert-Beer's constant, k_{λ} . However, here we deal with a chemical reaction. According to the equilibrium, $\text{Fe}(\text{SCN})^{2+}$ in solution will follow the reverse reaction forming iron(III) and thiocyanate ions. The shift due to the equilibrium makes the concentration of $\text{Fe}(\text{SCN})^{2+}$ unknown.

To solve this problem Fe^{3+} ion will be added in excess (0.0025 M) while adding only small amounts of SCN^{-} (10^{-4} M). This trick will displace the equilibrium to the right: according to Le Châtelier's, the equilibrium shifts to the products when one reactant is added. Hence, if one of the reactants is added in overwhelming excess, the other reactant will be consumed almost to exhaustion forming the product. Therefore, the FeSCN^{2+} concentration can be calculated using the concentration of SCN^{-} .

In this part, a set of solutions will be added to a clean cuvette, and absorbance at $\lambda=450$ nm will be measured using the spectrophotometer. With the 5 points obtained and the origin of the coordinates, a graph will be plotted to calculate the slope of the line. Notice again that in this linear representation, the line must pass through the origin by definition, since at concentration=0, absorbance must be 0.

Parts B. Determining the equilibrium constant K_c .

Using the absorbance measured and the value of $k_{450\text{nm}}$, in this second part of the experiment you will calculate the concentration in equilibrium of all species involved in the equilibrium. First, the concentration of the product will depend on the absorbance, and will be given by the formula below:

$$A = k_{450\text{nm}}[\text{Fe}(\text{SCN})^{2+}] \quad [\text{Fe}(\text{SCN})^{2+}] = \frac{A}{k_{450\text{nm}}}$$

An ICE chart given below will be used to calculate the equilibrium concentration of reactants, where (3) is calculated from the absorbance, and (2) and (1) are recalculated from the initial concentrations:

	$[\text{Fe}(\text{NO}_3)_3]$	$[\text{KSCN}]$	$[\text{Fe}(\text{SCN})^{2+}]$
Initial	(1)	(2)	0 M
Change	-(3)	-(3)	+(3)
Equilibrium	(1) - (3)	(2) - (3)	(3)

Procedure




Preparing a diluted KSCN solution



Step 1: – Clean 5 large test tubes and the matching stoppers and set them upside down on a rack. Number the test tubes from 1 to 5.

- ☐ *Step 2:* – Obtain a 100.0 mL volumetric flask and a Mohr pipette.
- ☐ *Step 3:* – Obtain about 15 mL of 0.0025 M KSCN in a 50 mL beaker. Rinse the Mohr pipette a couple of times with a small amount of this solution. Always discard the rinsing liquid. Then, transfer 4.00 mL into the volumetric flask.
- ☐ *Step 4:* – In a clean 100 mL beaker, obtain about 100 mL of distilled water. Carefully add water to the volumetric flask directly from the beaker and stop adding before you reach the 100.0 mL mark. Use a plastic dropper to level the water to the mark, drop by drop. This is the diluted KSCN solution.
- ☐ *Step 5:* – Rinse the Mohr pipette with the diluted solution of KSCN again.

Good Lab Practice

-  Pipettes are always used together with a suction bulb or a syringe. **Never** suck a chemical using the mouth!
-  **Never** leave a pipette lying on the table **while connected to the syringe**. Remaining liquid can leak into the syringe. The syringe will be damaged and the liquid will return to the pipette contaminated.
-  Always wash pipettes with distilled water and rinse them using the same product to be measured. Suck a small volume, disconnect the syringe and move the liquid about, displacing the distilled water. Discard the rinsing liquid as waste. Rinse twice.

Part A. Determining Lambert-Beer's constant k . Preparing the mixtures

- ☐ *Step 1:* – Following the Table in the description of Part A, add 1.00 mL of the diluted KSCN solution to jumbo test tube 1, 2.00 mL to test tube 2, 3.00 mL to test tube 3, 4.00 mL to test tube 4 and 5.00 mL to test tube 5. You will use jumbo test tubes throughout the experiment.
- ☐ *Step 2:* – Rinse the pipet with distilled water. Obtain about 50 mL of 0.25 M $\text{Fe}(\text{NO}_3)_3$ in a 100 mL beaker. Rinse the pipette with the solution of $\text{Fe}(\text{NO}_3)_3$ from the beaker.
- ☐ *Step 3:* – Add 5.0 mL of 0.25 M $\text{Fe}(\text{NO}_3)_3$ in each test tube.
- ☐ *Step 4:* – Rinse the pipette with distilled water. Obtain about 40 mL of 0.1 M HNO_3 in a 50 mL beaker. Rinse the pipette with the solution of HNO_3 from the beaker.
- ☐ *Step 5:* – Following the Table below, add 4.00 mL of the HNO_3 solution to test tube 1; 3.00 mL to test tube 2, 2.00 mL to test tube 3, and 1.00 mL to test tube 4.
- ☐ *Step 6:* – At this point, all your test tubes should have the same volume of liquid. If not, repeat the ones that diverge.
- ☐ *Step 7:* – Put the stoppers in the test tubes and mix the solutions.

Part A	$\text{Fe}(\text{NO}_3)_3$	KSCN	HNO_3	V_{total}
	0.25 M	0.0025 M	0.1 M	
Test Tube 1	5.0 mL	1.0 mL	4.0 mL	10.0 mL
Test Tube 2	5.0 mL	2.0 mL	3.0 mL	10.0 mL
Test Tube 3	5.0 mL	3.0 mL	2.0 mL	10.0 mL
Test Tube 4	5.0 mL	4.0 mL	1.0 mL	10.0 mL
Test Tube 5	5.0 mL	5.0 mL	0.0 mL	10.0 mL

Part A. Determining Lambert-Beer's constant k . Measuring absorbance

- ☐ *Step 1:* – Start the spectrophotometer. The bulb needs time to heat up to a stable temperature.
- ☐ *Step 2:* – Get a 250 mL beaker for waste. Refill the 100 mL beaker with distilled water.
- ☐ *Step 3:* – Using the plastic dropper, fill the cuvette 3/4 full with distilled water.
- ☐ *Step 4:* – Insert the cuvette into the spectrophotometer, set the wavelength to 450 nm, and press "blank" or "zero".
- ☐ *Step 5:* – Use a new cuvette. Fill the cuvette 3/4 full with the solution from test tube 1.
- ☐ *Step 6:* – Insert the cuvette into the spectrophotometer to measure the absorbance. Record the result.
- ☐ *Step 7:* – Repeat for test tubes 2, 3, 4, and 5, always using new cuvettes for each measurement.
- ☐ *Step 8:* – Plot absorbance (7) vs. concentration (6) in order to obtain the Lambert-Beer's constant k .
- ☐ *Step 9:* – Measure the absorbance of an unknown solution and use the linear regression (A/k) to obtain its concentration.
- ☐ *Step 10:* – Do not proceed unless you show the instructor your graph.

Good Lab Practice

- ✎ For best results, use the same cuvette, with the same orientation, for all the measurements.
- ✎ Keep the outer walls of the cuvette clean. Do not touch them. Wipe them with a piece of optical paper.
- ✎ Rinse the cuvette at the same time as you rinse the pipette.

Part B. Determining the equilibrium constant K . Preparing a different set of mixtures

- ☐ *Step 1:* – Prepare a new set of mixtures based on the proportions reported in the Table below.
- ☐ *Step 2:* – Make sure you use the 0.0025 M $\text{Fe}(\text{NO}_3)_3$ solution and not the 0.25M.
- ☐ *Step 3:* – The 0.0025 M KSCN solution to be used has already been prepared at the lab.
- ☐ *Step 4:* – Measure the absorbance for the new set of mixtures.

Part B	$\text{Fe}(\text{NO}_3)_3$ 0.0025 M	KSCN 0.0025 M	HNO_3 0.1 M	v_{total}
Test Tube 6	1.0 mL	5.0 mL	1.0 mL	7mL
Test Tube 7	1.0 mL	4.5 mL	1.5 mL	7mL
Test Tube 8	1.0 mL	4.0 mL	2.0 mL	7mL
Test Tube 9	1.0 mL	3.5 mL	2.5 mL	7mL
Test Tube 10	1.0 mL	3.0 mL	3.0 mL	7mL
Test Tube 11	2.0 mL	4.0 mL	1.0 mL	7mL
Test Tube 12	2.0 mL	3.5 mL	1.5 mL	7mL
Test Tube 13	2.0 mL	3.0 mL	2.0 mL	7mL
Test Tube 14	2.0 mL	2.5 mL	2.5 mL	7mL
Test Tube 15	2.0 mL	2.0 mL	3.0 mL	7mL

Calculations

① This is the volume of $\text{Fe}(\text{NO}_3)_3$ added.

② This is the volume of KSCN added.

③ This is the total volume of the mixtures.

④ This is the initial Iron(III) concentration in the mixtures:

$$[\text{Fe}^{+3}]_0 = \frac{v_{\text{Fe}^{+3}} \cdot c_{\text{Fe}^{+3}}}{v_{\text{total}}} = \frac{\textcircled{1} \cdot 0.25}{\textcircled{3}}$$

⑤ This is the initial thiocyanide concentration in the mixtures:

$$[\text{SCN}^-]_0 = \frac{v_{\text{SCN}^-} \cdot c_{\text{SCN}^-}}{v_{\text{total}}} = \frac{\textcircled{2} \cdot 0.0025}{\textcircled{3}}$$

⑥ This is the concentration of $\text{Fe}(\text{SCN})^{+2}$ in the mixture:

$$[\text{Fe}(\text{SCN})^{+2}] = [\text{SCN}^-] = \textcircled{5}$$

⑦ This is the measured absorbance of the mixture.

⑧ This step is needed to calculate numerically the Lambert-Beer's constant k :

$$[\text{Fe}(\text{SCN})^{+2}] \cdot A = \textcircled{6} \cdot \textcircled{7}$$

⑨ This step is needed to calculate numerically the Lambert-Beer's constant k :

$$[\text{Fe}(\text{SCN})^{+2}]^2 = \textcircled{6}^2$$

⑩ This step is needed to calculate numerically Lambert-Beer's constant k , and results from adding all values of ⑧ .

⑪ This step is needed to calculate numerically Lambert-Beer's constant k , and results from adding all values of ⑨ .

⑫ This step the Lambert-Beer's constant k :

$$k = \frac{\sum [\text{Fe}(\text{SCN})^{+2}] \cdot A}{\sum [\text{Fe}(\text{SCN})^{+2}]^2} = \frac{\textcircled{10}}{\textcircled{11}}$$

⑬ This is the volume of $\text{Fe}(\text{NO}_3)_3$ added.

⑭ This is the volume of KSCN added.

⑮ This is the total volume of the mixtures.

⑯ This is the initial Iron(III) concentration in the mixtures:

$$[\text{Fe}^{+3}]_0 = \frac{\nu_{\text{Fe}^{+3}} \cdot c_{\text{Fe}^{+3}}}{\nu_{\text{total}}} = \frac{(13) \cdot 0.0025}{(15)}$$

(17) This is the initial thiocyanide concentration in the mixtures:

$$[\text{SCN}^-]_0 = \frac{\nu_{\text{SCN}^-} \cdot c_{\text{SCN}^-}}{\nu_{\text{total}}} = \frac{(14) \cdot 0.0025}{(15)}$$

(18) This is the measured absorbance of the mixture.

(19) This is the concentration of $\text{Fe}(\text{SCN})^{+2}$ in equilibrium:

$$[\text{Fe}(\text{SCN})^{+2}]_{eq} = \frac{A}{k} = \frac{(18)}{(12)}$$

(20) This is the concentration of Fe^{+3} in equilibrium:

$$[\text{Fe}^{+3}]_{eq} = [\text{Fe}^{+3}]_0 - [\text{Fe}(\text{SCN})^{+2}]_{eq} = (16) - (19)$$

(21) This is the concentration of SCN^- in equilibrium:

$$[\text{SCN}^-]_{eq} = [\text{SCN}^-]_0 - [\text{Fe}(\text{SCN})^{+2}]_{eq} = (17) - (19)$$

(22) This is the equilibrium constant:

$$K = \frac{[\text{Fe}(\text{SCN})^{+2}]_{eq}}{[\text{Fe}^{+3}]_{eq} \cdot [\text{SCN}^-]_{eq}} = \frac{(19)}{(20) \cdot (21)}$$

STUDENT INFO

Name:

Date:

Pre-lab Questions

Equilibrium constant


1. You prepare a solution by following the procedure described next. You first obtain about 15 mL of 0.0025 M KSCN in a 50 mL beaker. Then you transfer 4.00 mL into the volumetric flask. In a clean 100 mL beaker, you obtain about 100 mL of distilled water, and carefully add water to the volumetric flask directly from the beaker stopping before reaching the 100.0 mL mark. Finally, you use a plastic dropper to level the water to the mark, drop by drop. Calculate the molarity of the solution you prepared.
2. You mix 1mL of $\text{Fe}(\text{NO}_3)_3$ 0.25 M with 5mL of KSCN 10^{-4} M and 1mL of HNO_3 0.1 M into a volumetric flask. Sketch a volumetric flask and calculate the initial concentration of Fe^{3+} in the mixture.
3. Write down the formula for the equilibrium constant in terms of concentration for the equilibrium below:
$$\text{Fe}^{3+}(\text{aq}) + \text{SCN}^{-}(\text{aq}) \rightleftharpoons \text{Fe}(\text{SCN})^{2+}(\text{aq})$$
4. A set of absorbances, A , for different concentrations, c , are given below.

(a) Plot A vs. c in the graph below.

(b) Compute the slope of the graph by using the formula:

$$k = \frac{\sum c \cdot A}{\sum c^2}$$

c (M)	A	$c \cdot A$ (M)	c^2 (M ²)
0.0120	0.681		
0.00960	0.540		
0.00720	0.389		
0.00480	0.270		
0.00240	0.133		
Sum			
k (M ⁻¹)			



(c) After calculating k , now write down the formula for absorbance in the form: $A = k \cdot c$

Results

Date:

Parts A. Determining Lambert-Beer's constant k .

[illegible]

Parts B. Determining the equilibrium constant K .

Test Tube	$\text{Fe}(\text{NO}_3)_3$	KSCN	HNO_3	ν_{tot}	$[\text{Fe}^{3+}]_0$	$[\text{SCN}^-]_0$	A	$[\text{Fe}(\text{SCN})^{2+}]_{\text{eq}}$	$[\text{Fe}^{3+}]_{\text{eq}}$	$[\text{SCN}^-]_{\text{eq}}$	K
	0.0025M	0.0025M	0.1M								
	(13)	(14)		(15)	(16)	(17)	(18)	(19)	(20)	(21)	(22)
6	1mL	5.0mL	1.0mL	7.0mL	_____	_____	_____	_____	_____	_____	_____
7	1mL	4.5mL	1.5mL	7.0mL	_____	_____	_____	_____	_____	_____	_____
8	1mL	4.0mL	2.0mL	7.0mL	_____	_____	_____	_____	_____	_____	_____
9	1mL	3.5mL	2.5mL	7.0mL	_____	_____	_____	_____	_____	_____	_____
10	1mL	3.0mL	3.0mL	7.0mL	_____	_____	_____	_____	_____	_____	_____
11	2mL	4.0mL	1.0mL	7.0mL	_____	_____	_____	_____	_____	_____	_____
12	2mL	3.5mL	1.5mL	7.0mL	_____	_____	_____	_____	_____	_____	_____
13	2mL	3.0mL	2.0mL	7.0mL	_____	_____	_____	_____	_____	_____	_____
14	2mL	2.5mL	2.5mL	7.0mL	_____	_____	_____	_____	_____	_____	_____
15	2mL	2.0mL	3.0mL	7.0mL	_____	_____	_____	_____	_____	_____	_____

Average K = _____




Figure 5: Column 6 $[\text{Fe}(\text{SCN})^{+2}]$ (X axis) vs. A Column 7 (Y axis)

STUDENT INFO

Name:

Date:


Post-lab Questions

Equilibrium constant

1. The absorbance of a colored substance in a colorless liquid is measured at each of a series of wavelengths, and the data is given below:

$\lambda(\text{nm})$	325	350	375	400	425	450	475	500	525
A	0.016	0.144	0.341	0.578	0.681	0.558	0.281	0.092	0.031

- (a) Plot A vs. λ in the graph below.
- (b) Calculate the λ value that gives a maximum A .



$\lambda(A_{max}) =$ _____

2. In parts A and B you use different iron and thiocyanate solutions: in part, A you use a 0.25M- $\text{Fe}(\text{NO}_3)_2$ solution and a 0.0025M-KSCN solution, whereas on part B you use a 0.0025M- $\text{Fe}(\text{NO}_3)_2$ solution and a 0.0025M-KSCN solution. Explain why you do this.
3. You mix 5mL of $\text{Fe}(\text{NO}_3)_3$ 0.25 M with 2mL of KSCN 0.0025 M and 3mL of HNO_3 0.1 M. Calculate the initial concentration of SCN^- in the mixture.
4. You mix 5mL of $\text{Fe}(\text{NO}_3)_3$ 0.25 M with 4mL of KSCN 0.0025 M and 1mL of HNO_3 0.1 M. Calculate the concentration of FeSCN^{2+} in the mixture.

8 Titration curve

EXPERIMENT

Titration Curve

Goal

The goal of this laboratory is to graph a titration curve and to use this curve to standardize an acid solution while determining the acidity constant of the acid.

Background

Titration curves were used in previous experiments to analyze the concentration of a solution, the so-called *analyte*. Using the *titrant*, a reactant of known concentration, it is possible to determine the concentration of the analyte if the volumes are measured carefully. The titrant is added from a buret until the solution changes color due to the presence of an indicator. The indicator denotes the moment when the *equivalence point* is achieved. The equivalence point is the moment at which equivalent molar amounts of both reactants have been added.

In an acid-base reaction, the acid reacts with the base, forming salt and water. This process is called *neutralization*; the protons released by the acid are neutralized by the hydroxide ions formed in the presence of the base. For strong acids and strong bases that dissociate completely in water, the equivalence point occurs at pH=7 and, therefore, an indicator with a color transition at pH=7 is used.

However, not all acids have the same strength. Weak acids do not release as many protons as strong acids do. The strength of a weak acid, and the consequent number of protons released, depend on its dissociation constant, K_a . Because of the partial dissociation of these species, the equivalence point in a weak acid titration does not occur at pH=7. The moment at which equivalent amounts of a weak acid and a strong base (titrant) have been added will take place at pH>7 (due to the non-complete dissociation). Determining the equivalence point is necessary to know the correct indicator to be used for titrations of weak acids.

The titration curve

During the neutralization process of an acid, the pH of the solution starts at a low value - lower than 7 since it is an acid - and increases with the addition of the base. If one plots the pH versus the volume of base added, the result is not a straight line, but an S-shape curve. This shape shows that the pH changes slowly at the beginning and the end of the process, but there is a region in the graph where the pH increases rapidly with small additions of the base. When that abrupt pH change takes place depends on the acid and the base used.

For the titration of a weak acid (analyte) using a strong base (titrant), the endpoint will be achieved at a pH higher than 7. For the titration of a weak base (analyte) using a strong acid (titrant), the endpoint will be achieved at a pH lower than 7. The titration curve for every base or acid is unique and it is useful for determining the dissociation constant (K_a or K_b) of the analyte.

The acid dissociation constant

The acid dissociation constant for a weak acid is a ratio of the amounts of products and reactants at the equilibrium position. The ratio is constant for a fixed temperature. The concentration of H^+ will always be among the products for the dissociation of an acid. Recalling that the concentration of H^+ defines the pH, the latter can be solved as follows:



$$K_a = \frac{[\text{H}^+][\text{A}^-]}{[\text{HA}]} \quad [\text{H}^+] = K_a \frac{[\text{HA}]}{[\text{A}^-]} \quad ; \quad \text{pH} = \text{p}K_a + \log \frac{[\text{A}^-]}{[\text{HA}]}$$

This expression is known as the *Henderson-Hasselbalch equation*, where [HA] is the concentration of the weak acid, and $[\text{A}^-]$ is the concentration of its conjugated base. It is possible to visualize the $\text{p}K_a$ from the titration curve by using the Henderson-Hasselbalch equation.

First, notice that at the equivalence point, the same amounts of acid and base have reacted and there will be no acid; $[\text{HA}] = 0$. But *halfway* to that point, when half of the acid has been converted to the same amount of its conjugated base, the concentrations of the two will be equal, $[\text{A}^-] = [\text{HA}]$. The fraction in the previous expression becomes equal to 1 and, since $\log 1 = 0$, then $\text{pH} = \text{p}K_a$.

$$\text{pH} = \text{p}K_a + \log \frac{[\text{A}^-]}{[\text{HA}]} \quad ; \quad \text{pH} = \text{p}K_a + \log 1 \quad ; \quad \text{pH} = \text{p}K_a$$

The equivalence point would correspond to the inflection point of the titration curve. By reading the volume of base added at that point, it is possible to calculate the volume needed to reach the halfway point ($V_{eq}/2$). The pH when that volume of base was added is taken from the titration curve and corresponds to the weak acid's $\text{p}K_a$.

The location of the equivalence point can be estimated by fitting two straight lines to the straight sections of the curve, as depicted in Figure 1, and a third line equidistant to the previous two. The point at which the third line and the titration curve intersect is the equivalence point.

The experiment

The experimental part consists of the titration only. The data will be used later to plot the curve. From the plotted curve the equivalence point and the acid $\text{p}K_a$ should be calculated.

The weak acid studied in this experiment is potassium hydrogen phthalate, $\text{KHC}_8\text{H}_4\text{O}_4$ (KHP). The base used as a titrant is sodium hydroxide (NaOH), a strong base. It is very important to obtain precise measurements of the volumes and masses used during this experiment.

The analyte (KHP) will be dissolved in a beaker. The titrant (NaOH) will be added from a buret. A pH meter will be used to monitor the pH and magnetic stirring will be used during the measurement to homogenize the solution pH.

Example

A solution of hydrogen cyanide (HCN), a weak acid, is titrated using NaOH, a strong base. The pH is monitored during the titration using a pH-meter. The following plot is obtained:




Figure 6: Titration curve of a HCN solution using NaOH. The numbers indicates the steps to follow as indicated in the answer.

Use the graphical method to identify the equivalence point for the titration, the volume of base added at the equivalence point, the midpoint, and the HCN K_a .

Answer: The equivalence point is obtained from the graph. (1) Draw the lines for the 2 straight sectors. (2) Draw the equidistant line. (3) Identify the equivalence point where the equidistant line intercepts the curve. (4) Draw a vertical line to find the volume; 50.0 mL. (5) The half-way point would be: $50.0 \text{ mL}/2 = 25.0 \text{ mL}$. (6) Draw a vertical line up to the curve. (7) From the curve draw an horizontal line to find the pH; 9.21. The $pK_a = 9.21$ (7), and therefore $K_a = 10^{-9.21} = 6.2 \times 10^{-10}$

Procedure

Preparing the setup

- ☐ *Step 1:* – Obtain a 50 mL buret, rinse it with distilled water, and verify that it operates correctly (you should be able to add drops but also larger quantities).
- ☐ *Step 2:* – Obtain a hot plate with magnetic stirring, a magnetic stirring rod, a stand, two 200 mL beakers, and a pH meter.
- ☐ *Step 3:* – Arrange the setup so that the beaker is on the hot plate, the pH meter is to its left, with the probe inserted in one of the beakers, and the tip of the buret is also slightly inserted into the beaker. The number markings on the buret should be facing you, and the valve should be accessible and operative.
- ☐ *Step 4:* – Rinse the buret with a small amount of 0.1 M NaOH solution. Discard the rinsing liquid. Fill the buret with fresh NaOH solution. Your instructor will either give you an unknown or instruct you to use the already prepared NaOH solution. Discard the first 3-4 mL (the volume of solution in the tip of the buret used for the rinsing).

Good Lab Practice

- 🔍 Burets are read from top to bottom.
- 🔍 The volume added is calculated by difference. It is not necessary to start from 0 mL.
- 🔍 Use a funnel to fill the buret but remove it always before the first reading. Drops from the funnel will affect the measurements.
- 🔍 Do **never** refill the buret with the Erlenmeyer/beaker of the analyte below it. A spill would ruin the titration.

Preparing the KHP solutions

- ☐ *Step 1:* – Put one of the beakers on the scale and set it to zero (T/O). The beaker must be dry. Remove the beaker from the scale to add a small amount of KHP. You need an amount between 0.45 and 0.55 grams. Record the value with all the decimals.
- ☐ *Step 2:* – Obtain 50.0 mL of distilled water in a graduated cylinder.
- ☐ *Step 3:* – Carefully add the 50.0 mL of distilled water from the graduated cylinder into the beaker containing the KHP.
- ☐ *Step 4:* – Repeat the procedure to prepare a second solution on a second beaker.

Measuring

- ☐ *Step 1:* – Add the magnetic stirrer to one of the beakers, put the beaker on the hotplate, and set a gentle stirring speed. Do not start the heat on the plate!
- ☐ *Step 2:* – Clean the pH-meter probe before immersing it in the KHP solution.

- ☐ *Step 3:* – Either you are given an unknown or you will use the NaOH already prepared from the lab. Read and record the initial pH before adding any NaOH.
- ☐ *Step 4:* – Start adding the titrant (the base). After every addition (see next step) record the **buret reading** (not the volume added) and the corresponding pH measurement.
- ☐ *Step 5:* – NaOH should be added initially in 1 mL increments, but at a certain moment, additions should be reduced to 0.2 mL portions. Mind that the pH will change very slowly at the beginning and the end of the titration. However, around the equivalence point, the pH will change extremely fast with small additions of the titrant.
- ☐ *Step 6:* – Read the buret and record the volume with all the digits given by the buret plus the estimated digit.
- ☐ *Step 7:* – Add the titrant in 1 mL portions making sure the buret readings are precise.
- ☐ *Step 8:* – After every addition, estimate how much the pH changed from the previous pH measurement. If the pH has changed by less than 0.3 pH units, add 1 mL of titrant in the next step. If the pH has changed by 0.3 pH-units or more, begin adding the titrant in 0.2 mL portions.
- ☐ *Step 9:* – When adding 0.2 mL portions, if the pH change starts to slow down and increases by less than 0.3 pH-units between measurements **for at least 3 measurements in a row**, then you can start adding 1 mL portions again.
- ☐ *Step 10:* – Stop the titration when the pH exceeds 11.5.
- ☐ *Step 11:* – Do a second titration trial. Clean and reuse the magnetic stirred; in the results section, plot by hand the titration plot.

Calculations

① This is the molarity of the base solution, written on the label.

② This is the mass of your sample.

③ These are the moles on your sample:

$$n_{\text{KHP}} = \frac{\textcircled{1}}{204.22}$$

④ These are the molarity on your acid sample:

$$c_a = \frac{\textcircled{2}}{50} \times 10^3$$

⑤ This is the burette reading in mL.

⑥ This is the PH as given by the OH meter.

⑦ This is the cumulative volume of base added.

STUDENT INFO

Name:


Date:

Pre-lab Questions**Titration Curve**

1. The data below represents the titration of a weak base with a strong acid. Use the data below to answer the following questions:

V_{HCl}^{Added} mL	PH	V_{HCl}^{Added} mL	PH	V_{HCl}^{Added} mL	PH	V_{HCl}^{Added} mL	PH
0	11.12	45.01	8.3	49.97	6.1	50.32	3.5
0.71	10.9	45.95	8.2	49.98	5.9	50.4	3.4
1.48	10.7	47.89	7.9	49.99	5.7	50.63	3.2
1.96	10.6	48.31	7.8	50.00	5.4	50.8	3.1
4.03	10.3	48.65	7.7	50.01	5	51.01	3
6.18	10.1	49.31	7.4	50.02	4.7	51.27	2.9
11.06	9.8	49.45	7.3	50.03	4.5	51.61	2.8
13.18	9.7	49.56	7.2	50.04	4.4	52.04	2.7
18.12	9.5	49.72	7.0	50.06	4.2	52.58	2.6
20.85	9.4	49.78	6.9	50.08	4.1	53.27	2.5
26.58	9.2	49.83	6.8	50.1	4	54.15	2.4
29.41	9.1	49.89	6.6	50.13	3.9	55.28	2.3
37.02	8.8	49.91	6.5	50.16	3.8	56.73	2.2
40.94	8.6	49.93	6.4	50.2	3.7	58.63	2.1
42.53	8.5	49.95	6.3	50.25	3.6	61.11	2

- 1.1. Plot PH (vertical axis) vs. volume of titrant added (horizontal axis), identifying the initial, mid, and endpoints.



1.2. Fill the table below based on the plot above:

Experimental Data						
	Initial Point		Midpoint		Endpoint	
	$V_{\text{HCl}}^{\text{Added}}$	PH	$V_{\text{HCl}}^{\text{Added}}$	PH	$V_{\text{HCl}}^{\text{Added}}$	PH
	0					

1.3. Use the midpoint to compute the basicity constant of the titrate given that at that point $\text{PH} = \text{PK}_b$. Write down the K_b value under the graph.

2. Calculate the PH of a 0.5M solution of an acid with $K_a = 1.4 \times 10^{-5}$.

STUDENT INFO	
Name:	Date:

Date:

Titration Curve

① $c_{\text{NaOH}} \text{ (M)} =$ _____

[illegible][illegible]

[illegible][illegible]

Experimental Data						
	Initial Point		Midpoint		Endpoint	
	$V_{\text{NaOH}}^{\text{Added}}$	PH	$V_{\text{NaOH}}^{\text{Added}}$	PH	$V_{\text{NaOH}}^{\text{Added}}$	PH
Trial 1	0					




Figure 7: PH (Y axis) vs. $V_{\text{NaOH}}^{\text{Added}}$ (X axis)

Experimental Data						
	Initial Point		Midpoint		Endpoint	
	$V_{\text{NaOH}}^{\text{Added}}$	PH	$V_{\text{NaOH}}^{\text{Added}}$	PH	$V_{\text{NaOH}}^{\text{Added}}$	PH
Trial 2	0					




Figure 8: PH (Y axis) vs. $V_{\text{NaOH}}^{\text{Added}}$ (X axis)

STUDENT INFO

Name:

Date:

Post-lab Questions

Titration Curve

1. Preparing a NaOH solution, starting from a solid NaOH sample, is difficult because NaOH is very hygroscopic (absorbs water from the air). This affects the mass measured and therefore the real concentration expected. The NaOH solution used was estimated to be 0.10 M but it might not be correct. Using the information from the titration (mass of the phthalate), and knowing that the acid and the base react 1:1, calculate the real concentration from each trial, and the mean value.

	Trial 1	Trial 2
KHP mass	_____	_____
KHP moles	_____	_____
V_{NaOH} @ equiv.	_____	_____
Real c_b , [NaOH]	_____	_____
Mean c_b , [NaOH]	_____	

2. Calculate the volume of NaOH for the halfway point (half of the volume added at the equivalence point) and locate in the graph the pH for the solution after that volume was added. Use the value of the pH to calculate the KHP pK_a .

V_{NaOH} @ halfway	_____	_____
pH @ halfway	_____	_____
pK_a	_____	_____
Mean pK_a	_____	

3. Using the pK_a mean value, calculate KHP's K_a . Write the whole expression for K_a .

4. Using the KHP's K_a and its initial concentration, calculate the theoretical pH before any NaOH was added, and compare the calculated values to the pHs measured experimentally.

	Trial 1	Trial 2
KHP moles	_____	_____
c_a , [KHP]	_____	_____
Theoretical pH	_____	_____
Experimental pH	_____	_____

5. Using the K_a for KHP, calculate K_b for its conjugated base, $\text{C}_8\text{H}_4\text{O}_4^{2-}$, given that $K_w = K_a \times K_b$.

9 Bleach Strength

EXPERIMENT

Strength of Laundry Bleach

Goal

The goal of this laboratory is to perform an **oxidation-reduction titration** to determine the concentration of sodium hypochlorite (NaOCl) in commercial bleach samples.

Materials

- | | |
|--|---|
| <input type="checkbox"/> A buret, a stand and a clamp | <input type="checkbox"/> 10% diluted bleach |
| <input type="checkbox"/> 100mL, 50mL beakers | <input type="checkbox"/> 10mL transfer pipet |
| <input type="checkbox"/> 125mL Erlenmeyer | <input type="checkbox"/> 2M HCl solution |
| <input type="checkbox"/> 0.025M $\text{Na}_2\text{S}_2\text{O}_3$ solution | <input type="checkbox"/> 0.2% starch solution |
| <input type="checkbox"/> solid KI | |

Background

The properties of chlorine for bleaching fabrics were discovered in the 18th century by the French scientist Claude Berthollet. He was able to make an aqueous solution of sodium hypochlorite which was named after the quarter in Paris where it was produced, Javel. *Eau de Javel* (Javel water) is still used in French to name bleach.

At the beginning of the 18th century, the pharmacist Antoine Germain Labarraque discovered the disinfecting ability of another hypochlorite solution, this time calcium hypochlorite. The so-called "Eau de Labarraque" was employed in the disinfection of animal guts to make musical instrument strings. Hypochlorite was the first antiseptic product, helping to treat gangrene and putrescent wounds in people in the 1820s, and ending with the propagation of "cadaveric particles" in morgues.

Nowadays bleach is a common product at home, used to whiten clothes, remove stains, and as a disinfectant. Because of its strong bactericidal properties, hypochlorite is also used for disinfecting and to prevent the proliferation of algae in swimming pools. Bleach can be consumed. The CDC recommends adding eight drops of bleach per gallon of clear water. For cloudy water, the bleach quantity should be doubled to 16 drops per gallon. Bleach can be found in powder form as well, in the form of Bleaching powder. This is indeed, a solid combination of chlorine and slaked lime, introduced in 1799 by Scottish chemist Charles Tennant.

Chlorine


Chlorine forms four different polyatomic anions with oxygen:



These are called **oxoanions** and illustrate 4 possible oxidation states of chlorine atoms; +1, +3, +5, and +7 respectively. When the oxidation state of the chlorine decreases, for instance going from ClO_4^- to ClO_3^- (+7 \rightarrow +5), it reduces. For reduction to happen another species must be oxidized. Chlorine is then acting as the **oxidizing agent**. Perchlorate is the strongest oxidizing agent of the 4 oxoanions, but even ClO^- can still be reduced to Cl^- (Cl oxidation state -1).

The first oxidation


To determine the concentration of ClO^- , two consecutive oxidation-reduction reactions will be performed. In the first reaction, all the ClO^- present in the sample will be oxidized into Cl^- . The **oxidizing agent** will be I^- , which will be in turn reduced into I_2 . The reaction, carried out in an acidic medium, will read as:



The reaction above is not balanced. You can determine the number of electrons being transferred using the oxidation states to help you balance the reaction. An excess amount of I^- will be used to ensure that all the ClO^- is oxidized. The iodine (I_2) formed in this reaction will be therefore limited by the amount of hypochlorite present in the sample. Keep in mind that the amount of hypochlorite is what we want to determine, and it can be stoichiometrically calculated from the amount of iodine produced. In the following step, the amount of I_2 will be determined using an oxidation-reduction titration. The solution might turn slightly brownish due to the presence of a Lewis complex between iodine and water.

The oxidation-reduction titration


The titration to be performed accounts for the following redox reaction:



Again, this reaction is not balanced. Iodine (I_2) is reduced back to iodide I^- . Thiosulfate $\text{S}_2\text{O}_3^{2-}$ is the reducing agent, being itself oxidized into tetrathionate ions ($\text{S}_4\text{O}_6^{2-}$). The titration endpoint will be determined by adding starch as an indicator at the beginning of the reaction. Remember from previous experiments that starch and I_2 react turning the solution into a dark blue color, therefore, the titration should be stopped when all the iodine has been reduced and the solution becomes transparent.

Redox reactions

Redox reaction are different than acid-base or precipitation reaction. They contain the same chemical element in two different states resulting from the loss or win of electrons. For example:



We have that neither $\text{Al}_{(\text{s})}$ or $\text{Cu}^{2+}_{(\text{aq})}$ are an acid or a base. Also there is no product precipitate. Hence, this reaction is not an acid-base reaction or a precipitation reaction. This is a redox reaction, as we have Al in two different states: as metallic $\text{Al}_{(\text{s})}$ and as ionic $\text{Al}^{3+}_{(\text{aq})}$, which result from the loss of an electron. In redox reaction there is always elements in the chemicals that lose electrons and chemicals winning electrons. The redox number, also called oxidation number or oxidation state helps compare elements that have lost (or gained) different charge. When comparing two elements with different redox numbers, the larger this number the more electrons the element has lost. Similarly, the smaller this number the smaller this number the more electrons the element has won.

We indicate redox numbers with roman number on top of the element. For example the redox number of manganese in this compound is +7: $\text{Mn}^{\text{VII}}\text{O}_4^-$. The redox number can be a non integer number. There is five rules to identify the redox number of an element.

- 🔗 **Rule 1** Single atoms or elements have zero redox number. Examples are Na or H_2 , both with redox zero.
- 🔗 **Rule 2** Monoatomic ions have redox number equal to their charge. Examples are Na^+ or Cl^- with redox +1 and -1, respectively.
- 🔗 **Rule 3** The redox number of fluorine is -1
- 🔗 **Rule 4** The redox number of hydrogen on its covalent (e.g. H_2O) compounds is +1.
- 🔗 **Rule 5** The redox number of oxygen in normal oxides (e.g. MgO) is normally -2, with the exception of peroxides (e.g. H_2O_2) in which is -1.

How do we calculate the redox number for example of manganese in this chemical: MnO_4^- , permanganate. In order to do this, we need to set up a formula so that the redox numbers of all elements in the molecule-taking into account the number of atoms in the molecule-equals to the charge. In the case of permanganate, let us call x to the redox number of manganese. We know the redox of oxygen is -2 and we have four oxygens in the molecule. We also know the charge of the ion is -1. Therefore we have:

$$x + 4 \cdot (-2) = -1$$

If we solve for x we obtain a redox number of manganese of VII.

Example

Calculate the redox number of the elements underlined in the following molecules: (a) $K_2\text{C}\underline{\text{O}}_3$ and (b) $\text{H}_2\underline{\text{C}}\text{O}$.

Answer: Let us set up the redox equation for the first compound, knowing that the redox of oxygen is -2 and potassium $+1$. The unknown variable x represents the redox number of the underlined element. We have:


$$2 \cdot (+1) + x + 3 \cdot (-2) = 0$$

Mind we have two potassium and three oxygens hence we need to time the redox of K by two and the redox of O by three. If we solve for x we obtain a redox number for carbon of IV. The redox equation for the second example is:


$$2 \cdot (+1) + x + (-2) = 0$$

Mind that according to the redox rules, the redox number of oxygen is $+1$. Solving for x we have a redox number of zero.


By comparing the redox number of the same element in two different compounds we can figure out in what compound the element has lost or gained electrons. Look for example the case of $\text{Cr}^{\text{VI}}_2\text{O}_7^{2-}$ and $\text{Cr}^{\text{III}}_2\text{O}_3$. The same element in two different molecules has two different redox numbers. In the case of dichromate, the redox of Cr is VI, whereas in the case of chromium(III) oxide the redox of Cr is III. The larger the redox number the more oxidized is an element, and that means the element has lost electrons. The smaller the redox number the more reduced is the element and that means it has gained electrons. If we compare both case, we have that Cr in dichromate is oxidized—it lost electrons—and Cr in chromium(III) oxide is reduced—it gained electrons. The goal is to identify the element that undergoes oxidation and reduction in a chemical reaction. We can reach this goal by using the half-reaction method. Every redox reaction is composed of two process, a reduction and the oxidation. These two processes can be separated into two half-reactions so that the combination of both half-reactions lead to the balanced redox. Let us work on a simple unbalanced redox reaction:




Solid $\text{Al}_{(\text{s})}$ on the reactant side has zero redox number, whereas ionic $\text{Al}_{(\text{aq})}^{+3}$ on the product side has redox number equal to III. Al has undergone oxidation as its redox number increases from zero to three. Al has lost three electrons. We can write the oxidation half-reaction:




Mind that electrons have negative charge and we add electrons to compensate the charge of $\text{Al}_{(\text{aq})}^{+3}$. Now let us compare the redox number of Cu. In the reactant side we have $\text{Cu}_{(\text{aq})}^{+2}$ with redox of II. In the product side we have $\text{Cu}_{(\text{s})}$ with zero redox number. Cu has undergone reduction as its redox number has decreases. This means it has gained electrons, in particular two electrons:




The goal here is to balance a redox chemical reaction by combining two half-reactions. In the example above the oxidation and reaction involve different number of electrons. Hence in order to be able to add both redox we need to time each half-reaction by a number so that the number of electrons cancel out. As the first reaction involved three electrons and the second two, we will do:




The overall balanced redox equation is:




Redox reactions happen in either basic or acidic medium. Here we will go over how to balance redox reactions in acidic medium. In order to do this we will first separate the reaction in two half-reactions. In each semi-reaction we will balance all elements but hydrogen and oxygen. When all elements are balanced, we will proceed to balance O by adding H_2O molecules and we will balance H by adding H^+ . Finally, we all add electrons to compensate the charge of the reaction. Let us work on an example:




One of the semi-reactions involve Manganese whereas the other involves Iron. The redox number of Mn in permanganate is VII hence Manganese is being reduced, as its redox number decreases from VII to II, whereas Iron is being oxidized as its redox number increases from II to III. The oxidation half-reaction does not contain hydrogen or oxygen hence we will only have to balance the charge with one electron:




The reduction half-reaction contains oxygen. Hence, we will have to add H_2O molecules to balance oxygen and H^+ to balance hydrogen. In particular, we will need two water molecules—as MnO_4^- has four oxygens and we will have to add eight protons as we are adding four molecules of water. Finally, we need to add three electrons to equalize the charge:




As the oxidation involves one electron and the reduction three, we need to time the oxidation by three:




In order to balance a redox in basic medium we need first to balance the reaction in acidic medium. After, we will compensate all H^+ with OH^- in both sides of the reaction. Mind that when combining H^+ with OH^- we obtain H_2O . For example, in order to balance the following reaction in basic medium:




we will add four OH^- in both sides:



And after cancelling the four protons with the four hydroxyls, we have:



Now we have four water molecules in the left and two in the right. We will cancel two molecules from each side in order not to report water molecules twice:




Example

Balance the following redox in acidic medium:




Answer:


We locate the same element in both sides of the reaction with different redox number. We found Chromium in the form of dichromate ($\text{Cr}_2\text{O}_7^{2-}$) with redox number VI and Chromium in the product side with redox III. Therefore Chromium is being reduced, as its redox number decreases. Differently, we found sulphur in the reactants side in the form of sulfite (SO_3^{2-}) with redox number of IV and in the product side with redox number of VI. Therefore Sulphur is being oxidized. We will first set up the oxidation half-reaction knowing that we have different amounts of Cr in both side and that we will have to add water molecules to balance O and protons to balance H. As we have seven oxygens in $\text{Cr}_2\text{O}_7^{2-}$ we will have to add seven water molecules. Also, as we add seven water molecules we will have to add fourteen protons. We will need six electrons to compensate the charge:




For the reduction half-reaction, we have a difference of one oxygen atoms and hence we will need one water molecule and two protons; we will need two electrons to compensate the charge:



In order to add both half-reactions as the reduction involves two electrons and the oxidation six, we will have to multiply the reduction half-reaction by three:



After we cancel the electrons and eliminate protons and water molecules, the balanced reaction is:



The experiment

The experimental part consists of a trial titration and an exact titration. The trial titration will be used to estimate the volume added at the endpoint of the titration. It is conducted quickly, adding portions of 1 mL until the endpoint is exceeded. The final result is not used for calculations. In the exact titration, an initial amount is added rapidly. This amount is calculated by subtracting 5 mL from the volume used in the trial titration. Like this it is possible to stop before the endpoint of the titration and to continue carefully, adding the titrant on a drop-by-drop basis until the precise endpoint. If time allows, the exact titration should be repeated and the mean volume used in calculations.

Procedure

Part A. Trial titration.

- ☐ *Step 1:* – Obtain a solution of bleach (the instructor might give you an unknown bleach sample, if so write down the Unknown # in the results section).
- ☐ *Step 2:* – Obtain a buret, a stand, and a clamp. Rinse the buret with distilled water and make sure the valve works properly.
- ☐ *Step 3:* – In a 100 mL beaker, obtain 70 mL of 0.025 M solution of $\text{Na}_2\text{S}_2\text{O}_3$. Use a small portion of the solution to rinse the buret. Discard the rinsing liquid and fill up the buret above the zero mark. Bring the volume in the buret to zero or below, discarding the excess liquid through the tip of the buret.
- ☐ *Step 4:* – In a 125 mL Erlenmeyer, obtain about one inch of a spatula of KI.

⚠ CAUTION!

- ⚠ Bleach is a base.
- ⚠ Bases and acids can cause burns and clothes color loss.
- ⚠ Always wear eye protection and handle carefully.

- ☐ *Step 5:* – In a 50 mL beaker, obtain about 30 mL of a 10% bleach solution.
- ☐ *Step 6:* – Obtain a 10 mL transfer pipet. Pipet 10.0 mL of the bleach into the KI containing Erlenmeyer. Add 20 mL of distilled water from a graduated column and 20 drops of HCl 2 M. Homogenize the mixture by stirring the flask.
- ☐ *Step 7:* – Record the buret initial reading to the closest hundredth of mL.
- ☐ *Step 8:* – Place the Erlenmeyer under the buret and start the titration by adding 1 mL portions until the solution turns from brown to yellow.
- ☐ *Step 9:* – Add 40 drops of a 0.2% starch solution to the yellow solution in the Erlenmeyer.
- ☐ *Step 10:* – Proceed with the titration adding the $\text{Na}_2\text{S}_2\text{O}_3$ in 1 mL portions until it turns from dark color to colorless.
- ☐ *Step 11:* – Record the buret final reading and calculate the volume added.
- ☐ *Step 12:* – Discard the solution in the Erlenmeyer in the corresponding waste container.

Part B. Exact titration.

- ☐ *Step 1:* – Refill the buret with $\text{Na}_2\text{S}_2\text{O}_3$ 0.025 M. At this point you know how much volume you will need from the trial titration. You don not need to fill the buret up to zero but add a few milliliters more than the expected volume.
- ☐ *Step 2:* – In a 125 mL Erlenmeyer, obtain about one inch of a spatula of KI.
- ☐ *Step 3:* – Pipet 10.0 mL of the diluted bleach into the KI containing Erlenmeyer. Add 20 mL of distilled water from a graduated column and 20 drops of HCl 2 M. Homogenize the mixture by stirring the flask.
- ☐ *Step 4:* – Record the buret initial reading to the closest hundredth of mL.
- ☐ *Step 5:* – Place the Erlenmeyer under the buret and add 3 mL less than the volume used in the trial titration.
- ☐ *Step 6:* – Add 40 drops of a 0.2% starch solution to the solution in the Erlenmeyer.
- ☐ *Step 7:* – Continue the titration adding the titrant drop by drop until the solution turns transparent.
- ☐ *Step 8:* – Record the buret final reading in your notes.
- ☐ *Step 9:* – Discard the solution in the Erlenmeyer. If the time allows, do two exact titration replicates.
- ☐ *Step 10:* – When finished empty the buret and rinse it with plenty of distilled water. The buret should be stored empty, open, and upside down.

Calculations

- ① This is the initial reading of the buret (often is zero if the buret is all filled up).
- ② This is the final reading of the buret after the endpoint has been reached.
- ③ This is the volume of titrant used, $\nu_{\text{S}_2\text{O}_3^{-2}}$:

$$\nu_{\text{S}_2\text{O}_3^{-2}} = \textcircled{2} - \textcircled{1}$$

- ④ This is the moles of thiosulfate used in the titration:

$$n_{\text{S}_2\text{O}_3^{-2}} = \textcircled{3} \times 0.025 \times 10^{-3}$$

5 These are the moles of bleach in the (diluted) sample:

$$n_{\text{OCl}^{-1}} = \frac{n_{\text{S}_2\text{O}_3^{2-}}}{2} = \frac{4}{2}$$

6 This is the molarity of bleach in the (diluted) sample:

$$c_{\text{OCl}^{-1}}^{\text{diluted}} = \frac{n_{\text{OCl}^{-1}}}{1 \times 10^{-2}} = \frac{5}{1 \times 10^{-2}}$$

7 This is the molarity of bleach in the (original, concentrated) sample:

$$c_{\text{OCl}^{-1}}^{\text{concentrated}} = c_{\text{OCl}^{-1}}^{\text{diluted}} \times 10 = 6 \times 10$$

STUDENT INFO


Name:

Date:


Pre-lab Questions

Strength of Laundry Bleach

1. Balance the following redox reaction, indicating the oxidizing and reducing agents:



2. Balance the following redox reaction, indicating the oxidizing and reducing agents:



3. What is the proportion of bleach to water necessary to make water safe to drink?
4. What is the active ingredient in liquid bleach? And the active ingredient in bleach powder?
5. Indicate the oxidation state of all elements in the following species:

× H^{+}

× H_2O

× Cl^{-}

× ClO^{-}

× I^{-}

× I_2

× $\text{S}_2\text{O}_3^{2-}$

× $\text{S}_4\text{O}_6^{2-}$

STUDENT INFO

Name:

Date:

**Results
EXPERIMENT****Strength of Laundry Bleach**

		Trial		Exact	
		1	2	3	
①	Initial buret reading (mL)	_____	_____	_____	_____
②	Final buret reading (mL)	_____	_____	_____	_____
③	$\text{Na}_2\text{S}_2\text{O}_3$ added (mL)	_____	_____	_____	_____
④	Moles $\text{Na}_2\text{S}_2\text{O}_3$	_____	_____	_____	_____
⑤	Moles NaOCl	_____	_____	_____	_____
⑥	[NaOCl] (diluted) (M)	_____	_____	_____	_____
⑦	[NaOCl] (original) (M)	_____	_____	_____	_____
Average Molarity, [NaOCl] (M)		_____			
Unknown number (if given)		_____			

STUDENT INFO

Name:

Date:

Post-lab Questions

Strength of Laundry Bleach

1. Commercial bleach claims a 5.25% NaClO by weight. Assuming a density of 1.0 g/mL for NaClO obtain the molarity of the solution. Consider your results in the laboratory, is the claim true?
2. Considering the bleach-water ratio of eight bleach drops per water gallon, what is the concentration, in molarity, of NaClO in the dilution?

10 Relative strength of some acids

EXPERIMENT

Acids Strength

Goal

The goal of this laboratory is to learn three different techniques to measure pH.

Materials

- | | |
|--|---|
| <input type="checkbox"/> 5 jumbo testtube and rack | <input type="checkbox"/> 0.1 M H_3PO_4 , 0.1 M CH_3COOH , 0.1 M NaH_2PO_4 , 0.1 M $\text{Zn}(\text{NO}_3)_2$, and 0.1 M NH_4NO_3 |
| <input type="checkbox"/> 100mL, 50mL beaker | <input type="checkbox"/> Distilled water and tap water |
| <input type="checkbox"/> Indicator solutions | <input type="checkbox"/> Unknown solution |

Background

Acids are considered dangerous substances and while it is true that special care must be taken when handling some of these products, many acids are used in daily life with great simplicity. For example, aspirin's active principle ingredient is acetylsalicylic acid, vinegar is a 5% solution of acetic acid, and vitamin C is ascorbic acid.

Acids are defined by Brønsted-Lowry as proton (H^+) donors because when dissolved in water, acids dissociate by releasing protons. Strong acids dissociate completely into two ions: the anion is the acid conjugate base and the cation is a proton. For stoichiometric reasons, the concentration of each ion in the solution equals the initial concentration of the acid. Weak acids do not dissociate completely but, instead, only do so to a certain extent. The strength of weak acids depends on how much they dissociate, a ratio characterized by the acid dissociation constant, K_a .

For a certain weak acid



The equilibrium can be characterized by the expression:

$$K_a = \frac{[\text{H}^+][\text{A}^-]}{[\text{HA}]}$$

where K_a (capital letter) is the dissociation constant and the brackets represent the concentrations of the respective species.

pH

The acidity depends on the concentration of H^+ in the solution, which will always be a value between 10^{-14} and 1. To avoid using such a large and complex range, a logarithmic scale is preferred, ranging from 0 to 14.

This scale is called the pH, defined by Danish biochemist Sørensen in 1909 as the "power of Hydrogen".


$$\text{pH} = -\log[\text{H}^+]$$

Considering that the strength of acids is their power to generate H^+ , and to compare the strength of different acids, pH should be measured.

In this experiment, the pH for different solutions will be measured using three different approaches; indicators, pH paper, and a pH meter. Each of these techniques has a different accuracy. While using **pH indicators** it is only possible to determine a range for the pH, **pH paper** will provide an estimated unit, and the **pH-meter** will read down to the decimals.

Indicators

There are several techniques to locate the endpoint during a titration. On one hand, we can use a PH meter and measure PH during the titration. The endpoint locates between the initial and final plateau. We can also use an automatic titrator, a machine that automatically mixes the acid and the base while displaying the PH. Another technique uses an indicator that changes color when the PH changes from acidic to basic or the opposite. Acid-base indicators are water-soluble organic dyes that have a different color depending on PH. Halochromism is the effect of changing color due to changes on a external force such as PH. Many flowers or plants exhibit this quality. For example, the french hydrangea is a plant that changes color depending on the acidity of the soil. Turmeric is a specie that turn from yellow to brown when the PH increases. indicator can be extracted from many plants and fruits such as red cabbage, geranium, poppy, rose petals or blueberries. Indeed, indicators are weak acids that exist in a conjugate base form as well. Both the acid (the indicator form with the hydrogen attached) and the conjugate base have a different color due to differences in their electronic structure that transform into differences absorbing light. For example, phenolphthalein is a common acid-base indicator that is transparent when the PH is acidic and pink when the PH is basic. The equivalency point of a titration (when acid and base are neutralized) is not necessarily the end point (when the indicator changes color). However, if both points are close enough chemical indicators are an accurate way to visually locate the equivalency point during a titration. As indicators are weak acids, each particular indicator is characterized by a constant of acidity referred to as K_{In} so that








The change in color is considered to happen when the difference between the acidic and basic form of the indicator is an order or magniture—ten times—different. This is the ten times rule. For this difference, we have that

$$\text{PH} = \text{PK}_{In} \pm 1$$

We can always use a different value for the difference (Δ) between the acid and basic form of an indicator so that the PH range will be given by

$$\text{PH} = \text{PK}_{In} \pm \log(\Delta)$$

This expression helps with the selection of the indicator for a given titration. For example, if we want to select an indicator that changes color when the difference between acid and base form is 5 times ($\Delta = 5$), then the OH range would be: $\text{PH} = \text{PK}_{In} \pm 0.7$. A good indicator will turn into a different color within the transition between acid and basic PH during a titration. For example, for the titration of hydrochloric acid with sodium hydroxide, the PH changes from 5 before the equivalency point to 8 after. Any indicator with PK_{In} between 5 and 8 would be able to accurately indicate the stoichiometric point. Examples are bromothymol blue, m-nitrophenol, and even phenolphthalein that transition color at 8. For the titration of a weak acid and a strong base, the selection of an indicator is more difficult as the difference between the acidic and basic plateaus during the titration is less pronounced. We must choose an indicator that changes color more closely to the end point. Litmus paper is an indicator heavily used for quick PH testing in labs. It contains a mixture of dyes that were extracted originally from lichens. There are two main types of litmus paper, red and blue litmus. Light blue litmus paper turns red when the medium is acidic, whereas red litmus paper turns blue in basic mediums. By wetting a piece of litmus paper one can also test for the acidity of gases, and for example, ammonia vapor would turn blue a piece of wet red litmus paper.

Acid-base indicators showing ample representations of their colors				
Name	PH Range	Color acid form		Color basic form
Thymol blue (acid range)	1.2-2.8	red		yellow
Thymol blue (basic range)	8.0-9.6	yellow		blue
Methyl orange	3.1-4.4	red		orange
Methyl red	4.4-6.2	red		yellow
Bromomethyl blue	6.0-7.6	yellow		blue

The experiment

The experiment today is divided into the following parts:

- Determining pH range using indicators.
- Determining pH using pH paper.
- Determining pH using a pH meter.

Part A. Determining pH range using indicators.

A pH indicator is a chemical substance that, when added to a solution, changes its color depending on the solution's pH. The pH at which color changes is called the **equivalence point**. Some indicators have one equivalence point while others have more than one equivalence points. By using an indicator it will be possible to determine a pH range only. Combining different indicators the range can be narrowed down.

Example

We test a sample with three different indicators by dividing the sample into three small samples. Each small sample is tested using an indicator from the following list: i) Indicator A is colorless for pH values lower than 10 and turns pink for a pH greater than 10. ii) Indicator B is blue up to pH 3, then it turns colorless. iii) Indicator C is colorless for a pH lower than 6, when it turns red. The red color changes to yellow between pH 9 and pH 10, and stays yellow for higher pH values. Estimate the pH range for the solution if the three samples are colorless.

Answer: Indicator A tells that the solution's pH is lower than 10. Therefore, the range must be $0 < \text{pH} < 10$. Indicator B shows that the solution's pH is greater than 3 (3 to 14). This information, combined with that of indicator A, narrows the pH range to $3 < \text{pH} < 10$. Indicator C tells that the solution's pH is lower than 5. The final range is $3 < \text{pH} < 5$.

Parts B. Determining pH using pH-paper.

pH paper, or **universal pH paper** is a strip of paper impregnated with a mixture of indicators. The paper changes color smoothly for each pH value and the pH is determined by matching the color of the paper with that in a chart attached to the pH paper. In the best case, by using pH paper it will be possible to select one value for the pH, although in some situations when the color is not one in the chart, it will be possible to indicate a range of two numbers only.

Parts C. Determining pH using a pH-meter.

pH meters are scientific instruments that measure the pH of solutions. The tool has a probe that is inserted in the liquid and measures the concentration of protons based on the voltage between two electrodes. Their use is similar to that of a conductivity-meter, and some instruments can perform both tasks. Probes need to be cleaned after each measurement to avoid contamination and should never be touched with the hands or paper towels. Calibration of the tool with solutions of known pH ensures the best accuracy. After use, probes are best kept moist, typically in distilled water. A pH meter can offer a precision down to the hundredth of the unit of pH.

Procedure

Part A. Determining pH range using indicators.

- ☐ **Step 1:** – Obtain 5 jumbo test tubes and a rack. Clean the test tubes and let them drain upside down in the rack.
- ☐ **Step 2:** – Obtain a 100 mL beaker with plenty of distilled water and transfer 4.0 mL in a graduated cylinder.
- ☐ **Step 3:** – Label one test tube with the letter R, and pour the 4.0 mL of distilled water into it. This test tube will be used as a volume reference for the following test tubes.
- ☐ **Step 4:** – In a 50 mL beaker get 40 mL of one of the following list of solutions (you will use the others further in the experiment). These solutions are going to be referred to here as the acid or base solutions.

- | | | | |
|-----------------------------------|-------------------------------------|-----------------------------------|---------------------|
| 1. 0.1 M H_3PO_4 | 3. 0.1 M NaH_2PO_4 | 5. 0.1 M NH_4NO_3 | 7. Tap water |
| 2. 0.1 M CH_3COOH | 4. 0.1 M $\text{Zn}(\text{NO}_3)_2$ | 6. Distilled water | 8. Unknown solution |

⚠ CAUTION!

- ⚠ Handle concentrated acids with care to avoid chemical burns.
- ⚠ Are you wearing your goggles? Do you know where the eye-washer is?

- ☐ *Step 5:* – Obtain a 250 mL beaker for waste, a medicine dropper or a plastic pipette
- ☐ *Step 6:* – Label the rest of the test tubes from 1 to 4. Use test tube R as a reference for the volume to add in the following test tubes. Rinse the dropper with the distilled water (discarding it in the waste beaker) between different solutions. Add 4.0 mL of the solution chosen to each test tube.
- ☐ *Step 7:* – Add two drops of Thymol Blue to test tube 1. Using the color table provided in the background section, indicate the PH range of the solution.
- ☐ *Step 8:* – Add two drops of Methyl Orange to test tube 2. Using the color table provided in the background section, indicate the PH range of the solution.
- ☐ *Step 9:* – Add two drops of Methyl Red to test tube 3. Using the color table provided in the background section, indicate the PH range of the solution.
- ☐ *Step 10:* – Add two drops of Bromothymol Blue to test tube 4. Using the color table provided in the background section, indicate the PH range of the solution.
- ☐ *Step 11:* – Record the colors on the results page.
- ☐ *Step 12:* – Reserve the acid or base solution for Parts B and C (the beaker you used to fill the test tubes). Set the beaker on a paper towel and label it.
- ☐ *Step 13:* – Discard the colored solutions tested with the indicator and clean the test tubes.
- ☐ *Step 14:* – Repeat steps 4 to 12 all other solutions given in Step 4. Make sure to rinse the dropper.

Part B. Determining pH using pH-paper.

- ☐ *Step 1:* – Obtain a stirring rod and the universal indicator paper. Be careful not to contaminate the paper with dirty gloves or bare hands. Always set the (small) pieces of paper on a clean paper towel.
- ☐ *Step 2:* – Rinse the stirring rod with the wash bottle over the waste beaker. Rinse thoroughly between different solutions.
- ☐ *Step 3:* – Immerse the tip of the stirring rod in one solution and touch a piece of pH paper with the rod.
- ☐ *Step 4:* – Match the color of the drop on the paper with the color chart on the plastic case. Record the pH value in the results section.
- ☐ *Step 5:* – Repeat the same process; rinsing, dipping, and touching, for the rest of the solutions.
- ☐ *Step 6:* – Reserve the solutions in the beakers for part C.

Part C. Determining pH using a pH-meter.

- ☐ *Step 1:* – Install the pH meter. The probe of the instrument is immersed in a buffer solution that must be kept aside until the end of the experiment.
- ☐ *Step 2:* – Press Calibrate (top right menu) and make sure the PH is set to the buffer solution's value (i.e. 4.00).
- ☐ *Step 3:* – Cleaning and measuring. It is very important to follow the cleaning steps before each measurement. You will need a 100 mL beaker labeled waste a wash bottle and a 100 mL beaker with distilled water.
 - 1. Rinse the probe using the wash bottle over the waste beaker.
 - 2. Dip the probe 3-4 times into the 250 mL beaker with clean distilled water.
 - 3. Dip the probe into the solution to be measured. Use a 50 mL beaker with at least 25 mL of solution.
 - 4. Repeat cleaning steps 1 and 2.
- ☐ *Step 4:* – Following the procedure outlined in the step above, measure the eight solutions.
- ☐ *Step 5:* – Record all the values.
- ☐ *Step 6:* – Make sure you leave the PH meter's probe immersed in the buffer solution, as you originally found it.

STUDENT INFO

Name:

Date:

Pre-lab Questions

Acids Strength

- Find the definition of the following concepts:
 - Acids (according to Brønsted-Lowry definition).
 - pH.
- You add a few drops of thymol blue to a clear solution and its color quickly changes to red. What is the pH range of the solution?
- You add a few drops of thymol blue to a clear solution and its color quickly changes to orange. What is the pH range of the solution?
- Name the three methods to determine pH to be used in this laboratory.
- Order the three methods from the most accurate to the least accurate, according to your expectations.
- We test a sample with three different indicators by dividing the sample into three small samples. Each small sample is tested using an indicator from the following list: i) Indicator A is colorless for pH values lower than 10 and turns pink for a pH greater than 10. ii) Indicator B is blue up to pH 3, then it turns colorless. iii) Indicator C is colorless for a pH lower than 6 when it turns red. The red color changes to yellow between pH 9 and pH 10, and stays yellow for higher pH values. Estimate the pH range for the solution if the samples tested with indicators A and B become colorless and the one tested with indicator C turns reddish.

STUDENT INFO

Name:

Date:

Results
EXPERIMENT

Acids Strength

Solution	Indicators						PH Paper	PH Meter
	Thymol Blue		Methyl Orange		Methyl Red		PH	PH
	Color	PH Range	Color	PH Range	Color	PH Range		
H ₃ PO ₄	_____	_____	_____	_____	_____	_____	_____	_____
CH ₃ COOH	_____	_____	_____	_____	_____	_____	_____	_____
NaH ₂ PO ₄	_____	_____	_____	_____	_____	_____	_____	_____
Zn(NO ₃) ₂	_____	_____	_____	_____	_____	_____	_____	_____
NH ₄ NO ₃	_____	_____	_____	_____	_____	_____	_____	_____
Distilled water	_____	_____	_____	_____	_____	_____	_____	_____
Tap water	_____	_____	_____	_____	_____	_____	_____	_____
Unknown	_____	_____	_____	_____	_____	_____	_____	_____
Unknown # _____							Unknown Identity _____	

STUDENT INFO

Name:

Date:

Post-lab Questions

Acids Strength

1. Based on your observation, identify the unknown.
2. Based on your observation, which of the three techniques is the most accurate?
3. Calculate the concentration of H_3O^+ ions in the most acidic solution, using the pH-meter value.
4. Calculate the concentration of H_3O^+ ions in the less acidic solution, using the pH-meter value.

11 Entropy

EXPERIMENT

Entropy and spontaneity

Goal

The goal of this laboratory experiment is to calculate the minimum entropy change for a chemical reaction involving the dissolution of sodium nitrate.

Materials

- ☐ A coffee-cup calorimeter
- ☐ $\text{NaNO}_3(\text{s})$
- ☐ A thermometer and a lit for the calorimeter

Background

Standard molar entropy

The standard molar entropy of a substance is the absolute entropy of one mole of substance at 1 atm. Entropies are absolute properties in contrast to other relative thermodynamic functions like the enthalpy in which a set of molecular states act as a reference. The tabulated standard entropy values are listed at the end of the chapter at one bar and 25 °C, the conventional temperature for reporting thermodynamic data. The units of molar entropy are J/molK. Molar entropies tend to be small values—as the unit joule is a small unit of energy. Finally, entropies can be positive or negative values.

Factors affecting entropy

We can qualitatively rationalize the trends between the standard molar entropy values for different substances and different conditions (temperature, volume, pressure) by using some simple rules described below:

🔗 **The state of matter:** The standard entropy of gases is larger than the standard entropy of liquids as gases present a larger degree of freedom (possible configurations) than liquids. The standard entropy of liquids is larger than the standard entropy of solids as liquids present a large degree of freedom than solids. The molecules of a liquid can jiggle more in comparison with the molecules of a solid. Overall, we can assume that liquids and solids have almost null entropy whereas gases have very large entropy. For example, the entropies of ice, water, and steam are respectively 41, 69.95, and 69.95 J/molK.

🔗 **Molar mass:** For monoatomic substances (e.g. Ne, Ar, etc.) the larger the atomic weight the larger entropy. This is because the energy levels of larger molecules are more packed, and hence are easier to access at a given temperature.

🔗 **Molecular complexity:** For substances with comparable molar mass (O_3 and F_2), the more complex the molecule the larger entropy, as the number of degrees of freedom is directly correlated with the number of microstates.

🔗 **Temperature:** Temperature increases entropy as the system has more accessible microstates. The plot below represents the entropy change with temperature for a substance.

🔗 **Number of particles:** The larger the number of particles of a system, the larger the entropy, as the more particles the more microstates or possible configurations.

🔗 **Volume:** The larger the volume of a system, the larger the entropy, as the larger volume the more microstates of possible configurations.

Calculating entropy changes in reactions

We can calculate the standard entropy of a reaction in a similar way as we calculate the standard enthalpy of a reaction:

$$\Delta S_R^\circ = \Delta S_{products}^\circ - \Delta S_{reactants}^\circ \quad \text{Entropy change} \quad (1)$$

where:

ΔS_R° is the standard entropy change of the reaction

$\Delta S_{products}^\circ$ is the standard entropy of all products

$\Delta S_{reactants}^\circ$ is the standard entropy of all reactants

It is important to take into account the stoichiometric coefficients. For example, for the reaction:



We have the entropy values of: $S^\circ(\text{CH}_3\text{OH}_{(g)})=239.7 \text{ J/K} \cdot \text{mol}$, $S^\circ(\text{O}_{2(g)})=161.1 \text{ J/K} \cdot \text{mol}$, $S^\circ(\text{CO}_{2(g)})=213.79 \text{ J/K} \cdot \text{mol}$, and $S^\circ(\text{H}_2\text{O}_{(l)})=69.95 \text{ J/K} \cdot \text{mol}$. We can calculate ΔS_R° :

$$\begin{aligned} \Delta S_R^\circ &= \Delta S_{products}^\circ - \Delta S_{reactants}^\circ = \\ &= \left(2 \cdot S^\circ(\text{CO}_{2(g)}) + 4 \cdot S^\circ(\text{H}_2\text{O}_{(l)}) \right) - \left(2 \cdot S^\circ(\text{CH}_3\text{OH}_{(g)}) + 3 \cdot S^\circ(\text{O}_{2(g)}) \right) \\ &= \left(2 \cdot 213.79 + 4 \cdot 69.95 \right) - \left(2 \cdot 239.7 + 3 \cdot 161.1 \right) = -255.32 \text{ J/K} \end{aligned}$$

Estimating the sign of the entropy changes in reactions

Oftentimes times we are more interested in predicting the sign of the entropy change of a reaction than in computing the exact value. This is because the sign can be used to estimate whether a reaction proceeds spontaneously. Two basic rules are used to estimate the entropy change sign:

¶ **The state of matter:** solids and liquids have very low entropy in comparison to gases. For example, in the case below, we have the production of liquid water from ice produces entropy as liquids have more entropy than solids:



On the other hand, the condensation of water vapor to produce a liquid consumes entropy, as liquids have less entropy than gases:



¶ **The number of molecules:** the larger the number of molecules of gas the larger entropy. For example, in the reaction below we have that we produce three molecules from two molecules. Hence, the entropy increases.



However, this rule only works if we only take into account the number of gas molecules. For example, in the reaction below we produce two gas molecules from three gas molecules and hence we lose entropy:



Gibbs free energy

The Gibbs free energy is just a combination of enthalpy and entropy for a given temperature:

$$G = H - T \cdot S \quad (2)$$

Working at constant temperature, we can compute the change in Gibbs free energy as:

$$\Delta G^T = \Delta H - T \Delta S \quad \text{Gibbs free-energy} \quad (3)$$

where:

ΔG^T is the Gibbs free-energy change at constant temperature

ΔH is the enthalpy change

ΔS is the entropy change

T is the temperature

Gibbs free energy is a state function that only depends on the final and initial state and not the path followed. At the same time, Gibbs free energy depends on temperature and pressure—we will discuss more about this at the end of the section. More importantly, the change in Gibbs free energy is associated with the spontaneity of the process. Gibbs free-energy decreases with temperature and the decrease is sharper for gases in comparison to liquids, and at the same time, the decrease is also sharper for liquids in comparison to solids.

Gibbs free energy and spontaneity

The Gibbs free-energy change of a reaction is associated with the spontaneity of the process or with its state of equilibrium. In particular, reactions that produce Gibbs free energy are nonspontaneous. Differently, reactions consuming free energy are indeed spontaneous. Finally, reactions without a change in free energy are in equilibrium.

⚡ $\Delta G < 0$: The reaction is spontaneous

⚡ $\Delta G = 0$: The reaction is in equilibrium

⚡ $\Delta G > 0$: The reaction is nonspontaneous

This is because the expression of Gibbs free energy is related to the total entropy change of the universe given by the second law of thermodynamics:

$$\Delta G^T = \Delta H - T\Delta S \quad \text{and} \quad T\Delta S_{univ}^{T,P} = -(\Delta H - T\Delta S) \quad (4)$$

Hence, we have that

$$\Delta G^T = T\Delta S_{univ}^{T,P} \quad (5)$$

This last expression means that an increase in total entropy is accompanied by a decrease in Gibbs free energy at constant pressure and temperature. At the same time, the expression indicates that the direction of spontaneous change corresponds to a decrease of Gibbs free energy, at constant temperature and pressure.

How to find extrapolated temperature


A calorimeter is ideally an isolated system. That means the heat produced does not leave the cup. If a reaction is exothermic—if it produces heat—the temperature inside the calorimeter will increase. Differently, for an endothermic reaction—a reaction that consumes heat—the temperature inside the calorimeter will decrease. When the reaction stops, as there is no more heat being produced or consumed, the temperature inside the calorimeter will change. To estimate the temperature at that point, we need to extrapolate the temperature inside the calorimeter. If you plot temperature versus time you will see that for early times, the temperature inside a calorimeter in which a reaction happens changes abruptly. After that, the temperature slowly decreases with time. To calculate the temperature when the reaction stops you need to connect the dots of the later stages, when the reaction already has stopped. We call this interpolation. All these points will be part of a line. If you extend the line, it will cross the Y-axis, and that crossing point is the temperature at which the reaction has stopped. We call this the extrapolation method.

Example

The following data reports the temperature inside a cup calorimeter. Indicate whether the reaction is exothermic or endothermic and calculate the temperature when the reaction stops employing the extrapolation method.

t (min)	0	1	2	3	4	5	6	7	8	9	10
T (°C)	10	10	20	28	29	30	29.5	29	28.5	27.5	27

Answer: if we plot temperature versus time we will see the two different regions of the plot. In early times, the temperature increased and hence the reaction in the calorimeter was exothermic. At later times, the reaction stops and temperature decreases due to the lack of heat release. We will use the last points and connect them with a line. The crossing point with the Y axis is 33 °C and that is the final temperature of the reaction.





Procedure

- ☐ *Step 1:* – Obtain a coffee cup calorimeter with a thermometer and a plastic lid.
- ☐ *Step 2:* – Weight about 8.5g of NaNO_3 . Write down the mass you weighed in the results section.
- ☐ *Step 3:* – Use a 100mL graduated cylinder to add 100mL of distilled water in the calorimeter.
- ☐ *Step 4:* – Measure the temperature of water in the calorimeter. Record this temperatures as your initial temperature t_i .
- ☐ *Step 5:* – Get the timer ready. Time should be taken as soon as the solid reactant is added.
- ☐ *Step 6:* – Add the NaNO_3 to the calorimeter. If you are using a magnetic stirrer set the speed to a medium speed and make sure the heating option is off. Then add the solid. Start the timer immediately and place the lid on the calorimeter.
- ☐ *Step 7:* – Record the temperature for 4 min. Do not stop the timer during the experiment. Record the results in the results table.
- ☐ *Step 8:* – Replicate the experiment one more time.
- ☐ *Step 9:* – Plot the temperature of the calorimeter against time and use a straight line to extrapolate your results against the starting time (time = 0 s). Plot the results for both replicates in the same plot. Record the extrapolated temperature as t_f .

- ☐ *Step 10:* – Do not carry out any calculation before finishing the plot from the previous step. Calculate the change of enthalpy, ΔH .
- ☐ *Step 11:* – Calculate the average enthalpy change using the two replicated experiments.
- ☐ *Step 12:* – If your instructor assigns you an unknown solute, repeat the previous procedure with this new solute.

Good Lab Practice

 Make sure you always use the same thermometer all through the experiment so that the results are consistent. 

Calculations

- ① This is the initial temperature before you add the solid.
- ② This is the final temperature after the dissolution reaction stops. You need to plot the data and extrapolate to get this value.
- ③ This is the maximum temperature during the experiment. Comparing the initial and final temperatures you can figure out if the reaction is endothermic or exothermic.
- ④ This is the mass of salt weighed.
- ⑤ This is the number of moles of salt:

$$\text{④} / 84.9947$$

- ⑥ You can calculate the enthalpy change (ΔH in KJ/mol) by means of the following formula:

$$\Delta H = - \frac{1 \times 10^{-3}}{\text{⑤}} \cdot (4.184 \cdot 1.0 \cdot 100.0 \cdot (\text{②} - \text{①}) + 10 \cdot (\text{②} - \text{①}))$$

- ⑦ This is the average enthalpy.
- ⑧ This is the average maximum temperature.
- ⑨ Calculate the minimum entropy change for equilibrium using the maximum temperature and the average enthalpy (remember enthalpy is expressed in KJ/mol and entropy in J/molK)

$$\Delta S_{min} = \frac{\Delta H}{T_m} = \frac{\text{⑦} \times 1000}{\text{⑧}}$$

STUDENT INFO

Name:

Date:


Pre-lab Questions

Entropy and spontaneity

1. Ammonium nitrate (NH_4NO_3) is dissolved in a coffee-cup calorimeter. We recorded the temperatures with time and the data is given below. The initial temperature, t_i was 25.0°C .

t (s)	0	30	60	90	120	150	180	210	240
T ($^\circ\text{C}$)	25	40.8	40.7	40.6	40.5	40.4	40.3	40.2	40.1

- (a) Plot T (vertical axis) vs. t (horizontal axis).



- (b) Obtain t_f , the final temperature after the reaction starts, by extrapolating to the time of mixing (time = 0 s) with a straight line.
- (c) Indicate whether this is an exothermic or an endothermic reaction.

(d) Explain the change of temperature inside the calorimeter.

This area in the pdf is writable only with Acrobat Reader...

2. Predict the sign of the entropy change for the dissolution of NH_4NO_3 .

This area in the pdf is writable only with Acrobat Reader...

3. Predict the spontaneity of the dissolution of NH_4NO_3 using the signs of the entropy and enthalpy estimated above.

This area in the pdf is writable only with Acrobat Reader...

4. Calculate the NaNO_3 (MW=84.9947 $\text{g} \cdot \text{mol}^{-1}$) mass needed to prepare 100mL of a 1M solution of this salt.

Results

Date:

NaNO ₃ Trial 1															
t (s)	0	30	60	90	120	150	180	210	240	T_i (°C) <div>(1)</div>	T_f (°C) <div>(2)</div>	T_m (°C) <div>(3)</div>	m (g) <div>(4)</div>	n (moles) <div>(5)</div>	$\Delta H(kJ/mol)$ <div>(6)</div>
T (°C)															
NaNO ₃ Trial 2															
t (s)	0	30	60	90	120	150	180	210	240	T_i (°C)	T_f (°C)	T_m (°C)	m (g)	n (moles)	$\Delta H(kJ/mol)$
T (°C)															
										Average ΔH (kJ/mol) =					
										Average T_m (K) =					
										ΔS_{min} (J/molK) =					

Unknown Trial 1																	
t (s)	0	30	60	90	120	150	180	210	240		T_i (°C)	T_f (°C)	T_m (°C)	m (g)	n (moles)	$\Delta H(kJ/mol)$	
											①	②	③	④	⑤	⑥	
T (°C)																	
Unknown Trial 2																	
t (s)	0	30	60	90	120	150	180	210	240		T_i (°C)	T_f (°C)	T_m (°C)	m (g)	n (moles)	$\Delta H(kJ/mol)$	
T (°C)																	
											⑦	Average ΔH (kJ/mol) =					
											⑧	Average T_m (K) =					
											⑨	ΔS_{min} (J/molK) =					




Figure 9: time (X axis) vs. Temperature (Y axis)

STUDENT INFO

Name:

Date:

Post-lab Questions**Entropy and spontaneity**

1. Predict the sign of the entropy change for the dissolution of NaNO_3
2. Use the data below to compute the entropy change in the dissolution process of sodium nitrate at 298K. Compare your experimental result with the calculated giving reasons for the difference.

	$\text{NaNO}_3(\text{s})$	\longrightarrow	$\text{Na}^+_{(\text{aq})}$	+	$\text{NO}_3^-_{(\text{aq})}$
$S^\circ (\text{J/molK})$	116.32		59.0		146.4
n					
$S^\circ \cdot n$					
$\Delta S^\circ (\text{J/molK})$					

3. Use the data below to compute the enthalpy change in the dissolution process of sodium nitrate at 298K. Compare your experimental result with the calculated giving reasons for the difference.

	$\text{NaNO}_3(\text{s})$	\longrightarrow	$\text{Na}^+_{(\text{aq})}$	+	$\text{NO}_3^-_{(\text{aq})}$
$H^\circ (\text{KJ/mol})$	-466.68		-240.1		-207.36
n					
$H^\circ \cdot n$					
$\Delta H^\circ (\text{KJ/mol})$					

4. Was the dissolution of sodium nitrate spontaneous in the conditions of the experiment?

5. Calculate the value of Gibbs free energy at 298K using the entropy and enthalpy calculated above.

12 A Solubility Product Constant: Solubility product calculation

EXPERIMENT

A Solubility Product Constant: Solubility product of an iodate salt

Goal

The goal of this laboratory experiment is to calculate the solubility product of a iodate salt.

Materials

- | | |
|--|--|
| <input type="checkbox"/> A 10mL pipet, a 150mL beaker, stirring rod | <input type="checkbox"/> solid KI _(s) |
| <input type="checkbox"/> Filter paper, funnel, stand, ring | <input type="checkbox"/> A 0.2M KIO _{3(aq)} solution |
| <input type="checkbox"/> A burette, stand and fisher clamp | <input type="checkbox"/> 2M HCl _(aq) and 2% starch solution |
| <input type="checkbox"/> A 1M Ca(NO ₃) _{2(aq)} solution | <input type="checkbox"/> A series of 125mL (or 250mL) erlenmeyers |


Background

When insoluble compounds dissolve in water an equilibrium between the solid and the dissolved ions establishes. Based on the stoichiometry of the compound we will obtain a different number of moles of ions in solution.


This experiment can be carried out in two different forms. The instructor can either give you an unknown saturated solution, or you can prepare it yourself. At the same time, the standardization of sodium thiosulfate can be carried out optionally, based on the instructor's indications.

Solubility equilibrium

Insoluble compounds dissolve in small amounts in water. Solubility equilibrium is the equilibrium that describes the dissolution of an insoluble compound to produce an aqueous solution of ions. For example, in the case of silver chloride, we have:



The solubility equilibrium above represents the dissolution of the insoluble salt—a solid—to produce ions, silver, and chloride in solution. In general terms, solubility equilibria start with a solid and generate ions in solution. To break down the insoluble compound into ions, we just need to take into account the stoichiometry of the compound. For example, for the case of Cobalt(II) phosphate, we have:



Solubility product in terms of molarities, K_{sp}

Silver chloride, AgCl is an insoluble compound in water which means this compound will not fully dissolve in water. Still, small quantities of the salt will indeed dissolve leading to a small amount of silver and chlorine ions. As this compound is insoluble the ion concentration will be very small. The dissociation equilibrium involved in the dissolution indicated below, is characterized by an equilibrium constant K_{sp} called *solubility product constant* or simply *solubility product*:



As pure solids are not included in any equilibrium constant, the formula above does not include $\text{AgCl}_{(s)}$. Hence, all solubility products simply result from the product of the molarities of the ions involved in the dissociation with the corresponding powers based on the stoichiometry of the salt. This is why the term *product* is included in its name. K_{sp} will have a different explicit expression depending on the formula of the insoluble compound—depending on its stoichiometry. For example, K_{sp} for calcium fluoride—a 1:2 compound—would be:



Solubility products are related to the solubility of chemicals. However, the relationship is not one-on-one. In other words, a larger K_{sp} does not necessarily imply a larger solubility. For example, K_{sp} for PbBr_2 is $6.6 \times 10^{-3}\text{M}$ whereas K_{sp} for MgCO_3 is $4.0 \times 10^{-5}\text{M}$. We have that K_{sp} for PbBr_2 is smaller than K_{sp} for MgCO_3 . However, the solubility of PbBr_2 is indeed larger than the solubility of MgCO_3 .

Relating solubility to the concentration of ions

We have that insoluble compounds dissolve to produce small amounts of ions in solutions. The concentration of these ions is directly related to solubility s using the stoichiometry of the compound. For example, for silver chloride—a salt with a 1:1 stoichiometry—we have that solubility has a 1 and 1 relationship with the concentration of ions:

$$[\text{Ag}^{+}] = s \text{ and } [\text{Cl}^{-}] = s$$

Similarly, for barium sulfate we have that solubility has a 1:1 relationship with the concentration of ions:

$$[\text{Ba}^{+2}] = s \text{ and } [\text{SO}_4^{-2}] = s$$

When the stoichiometry of the compound is not 1:1 we need to include the stoichiometry coefficients in the relationship between solubility and ion concentration. For example, for Ag_2SO_4 —a compound with 2:1 stoichiometry, we have that

$$[\text{Ag}^{+}] = 2 \cdot s \text{ and } [\text{SO}_4^{-2}] = s$$

This is because, for every mole of silver sulfate, we produce two moles of silver—and hence the factor two—and one mole of sulfate in solution. Similarly, for $\text{Nd}_2(\text{CO}_3)_3$ we have

$$[\text{Nd}^{+3}] = 2 \cdot s \text{ and } [\text{CO}_3^{-2}] = 3 \cdot s$$

Relating solubility to the concentration of ions

The solubility product is directly related to molar solubility s . We will demonstrate how to obtain this relationship using three examples. First, the solubility product of AgCl is:



As the concentration of each ion, Ag^{+} and Cl^{-} , is related to the molar solubility of the salt, we have

$$K_{sp} = [\text{Ag}^{+}] \cdot [\text{Cl}^{-}] = (s) \cdot (s) = s^2$$

Second, for silver sulfide, an insoluble compound with a more complex stoichiometry, we have:



Third, the solubility equilibrium and K_{sp} expression for $\text{Nd}_2(\text{CO}_3)_3$ is



As the ion concentrations are related to solubility, we have

$$K_{sp} = (2s)^2 \cdot (3s)^3 = 108s^5$$

K_{sp} in terms of molar solubility: general formula

For any insoluble salt A_xB_y , we have that K_{sp} is related to s by means of a general formula:

$$K_{sp} = a \cdot s^b \quad (6)$$

where:

$$a \text{ is } x^x \cdot y^y$$

$$b \text{ is } x + y$$

For example, for BaF_2 the constant a would be $1^1 \cdot 2^2$, that is four, whereas the constant b will be $1 + 2$ which equals three. As such, the expression of K_{sp} in terms of s would be: $K_{sp} = 4s^3$. This approach is useful when we need to compute the solubility product constant given the molar solubility.

Molar solubility in terms of K_{sp} : general formula

We previously explored the relationship between K_{sp} and molar solubility. Here we will explore the relationship between molar solubility and K_{sp} , simply solving for s in Equation 6. Again, for any insoluble salt A_xB_y , we have:

$$s = \left(\frac{K_{sp}}{a} \right)^{\frac{1}{b}} \quad (7)$$

where:

$$a \text{ is } x^x \cdot y^y$$

$$b \text{ is } x + y$$

For example, for BaF_2 the constant a would be $1^1 \cdot 2^2$, that is four, whereas the constant b will be $1 + 2$ which equals three. As such, the expression of s in terms of K_{sp} would be: $s = \left(\frac{K_{sp}}{4} \right)^{\frac{1}{3}}$. This approach is useful when the solubility product constant is given and we need to calculate the molar solubility. As the relation between molar solubility and K_{sp} is not a one-to-one relationship, K_{sp} is not directly related to solubility. For salts with a similar stoichiometry, it would be safe to compare solubilities in terms of K_{sp} . For example:

CuS	$K_{sp} = 8 \times 10^{-37}$	$s = 6 \times 10^{-19}$
PbS	$K_{sp} = 3 \times 10^{-28}$	$s = 1 \times 10^{-14}$
LiF	$K_{sp} = 2 \times 10^{-3}$	$s = 3 \times 10^{-2}$

and we have that the lower K_{sp} the lower solubility. When the salt stoichiometry differs

MgF ₂	$K_{sp} = 5 \times 10^{-11}$	$s = 2 \times 10^{-4}$
Li ₃ PO ₄	$K_{sp} = 2 \times 10^{-4}$	$s = 7 \times 10^{-2}$
Li ₂ CO ₃	$K_{sp} = 8 \times 10^{-4}$	$s = 6 \times 10^{-2}$

then an increase in K_{sp} does not necessarily follow an increase in solubility.

Procedure

Part A: Preparation of a saturated solution

- ☐ Step 1: – If your instructor assign you an unknow you can skip Part A.
- ☐ Step 2: – In a 150mL beaker mix 50mL of a 0.2M KIO₃ solution and 20mL of a 1M Ca(NO₃)₂ solution. Make sure you pick the reagents with the indicated concentration as in the lab there could be more diluted potassium iodate. Stir the resulting mix with a stirring rod. A precipitate will appear.
- ☐ Step 3: – Using gravity filtration (use a filter funnel and not a Büchner funnel), filter the precipitate. Make sure you rinse the beaker with water, to collect all precipitate leftovers.

- ☐ *Step 4:* – Wash the precipitate on the filter paper with three different volumes of water.
- ☐ *Step 5:* – Use a spatula to separate the precipitate into two halves. Dispose of each part in a clean labeled 100mL beaker. Put one of the beakers aside to be used in case of an accident. Use the other beaker in the rest of the experiment. If you do not use the first half of the precipitate, dispose of it at the end of the experiment following your instructor's indications.
- ☐ *Step 6:* – Add 40mL of distilled water to the beaker containing the precipitate. Use a stirring rod to stir the mixtures.
- ☐ *Step 7:* – This is the saturated solution. Let this solution rest for 30 minutes, occasionally stirring.
- ☐ *Step 8:* – It is important not to add additional water to this solution.

Part B: Standardization of thiosulfate: mixture preparation

- ☐ *Step 1:* – Locate the KIO_3 solution and write down its molarity in the results section.
- ☐ *Step 2:* – Prepare a diluted KIO_3 solution by pipetting 5mL of the KIO_3 solution into a 100mL volumetric flask and fill up with distilled water.
- ☐ *Step 3:* – Pipet 10mL of the diluted KIO_3 solution prepared in the previous step into a 125mL (or 250mL) Erlenmeyer flask. First, add 20mL of distilled water from a cylinder. Add 20 drops of 2M HCl and the equivalent of 1mL of solid KI, measured with a 10mL cylinder. Mind you can approximately measure the volume of a solid with a cylinder. Finally, add 40 drops of a 0.2% starch solution. Swirl the mixture to obtain a homogeneous mixture.
- ☐ *Step 4:* – Repeat the previous step two more times, preparing a total of three Erlenmeyers containing the iodate mixture.

Part B: Standardization of thiosulfate

- ☐ *Step 1:* – Ask your instructor of the thiosulfate solution is already standardized. If so, you can skip Part B.
- ☐ *Step 2:* – Rinse a buret with a small portion of $\text{Na}_2\text{S}_2\text{O}_3$ solution, discarding the solution used in the rinses.
- ☐ *Step 3:* – Fill a buret with $\text{Na}_2\text{S}_2\text{O}_3$ solution.
- ☐ *Step 4:* – Titrate one of the iodate samples adding increments of 1mL of sodium thiosulfate at a time. Swirl the mixture after each addition. You will have reached the equivalency point when the solution becomes colorless. Try placing the flask containing the analyte on a piece of paper. The white background will help see the color change when it occurs.
- ☐ *Step 5:* – Write down the volume of thiosulfate added to reach the equivalency point to the nearest 0.01mL. This will give an estimate of the volume needed to reach the endpoint.
- ☐ *Step 6:* – After having estimated the volume of titrate needed to reach the endpoint, we will now speed up the titration by rapidly titrating the second thiosulfate sample.
- ☐ *Step 7:* – Add thiosulfate ($\text{Na}_2\text{S}_2\text{O}_3$) to a second iodate sample, stopping 3mL before the estimated endpoint.
- ☐ *Step 8:* – Add $\text{Na}_2\text{S}_2\text{O}_3$ drop by drop until you reach the endpoint and the solution becomes colorless. Record the buret reading to the nearest 0.01mL.
- ☐ *Step 9:* – Repeat the titration with the last iodate sample.

Part C: Preparing three mixtures of the saturated solution

- ☐ *Step 1:* – The goal now is to prepare three samples of an acidified, saturated solution containing iodide. To do this we will first filter the saturated solution to eliminate the solid. After that, we will add acid and iodine to an exact volume of this solution.
- ☐ *Step 2:* – Filter the saturated solution obtained earlier in the experiment. Make sure the funnel and the paper are dry. Do not wash the precipitate with water. The filtrate is the solution to be titrated.
- ☐ *Step 3:* – Rinse a 10mL pipet with distilled water making sure the pipet is free of water at the end. Now rinse the same pipet two times with a small amount of the filtrate collected in the previous step (2mL), discarding the solution from the pipet.
- ☐ *Step 4:* – Pipet 10mL from the filtrate into a 125mL (or 250mL) Erlenmeyer flask. First, add 20mL of distilled water from a cylinder. Add 20 drops of 2M HCl and the equivalent of 1mL of solid KI, measured with a 10mL cylinder. Mind you can approximately measure the volume of a solid with a cylinder. Finally, add 40 drops of a 0.2% starch solution. Swirl the mixture to obtain a homogeneous mixture.
- ☐ *Step 5:* – Repeat the previous step two more times to prepare a total of three mixtures.

Part D: Titration of the saturated solution mixtures: trial titration

- ☐ *Step 1:* – The final goal here is to calculate the molarity of the saturated solution, which corresponds to the solubility of the insoluble compound. Here we will estimate the equivalence point of the titration of sodium thiosulfate and sodium iodate.
- ☐ *Step 2:* – Rinse a buret with a small portion of $\text{Na}_2\text{S}_2\text{O}_3$ solution, discarding the solution used in the rinses.
- ☐ *Step 3:* – Fill a buret with $\text{Na}_2\text{S}_2\text{O}_3$ solution.
- ☐ *Step 4:* – Titrate one of the samples adding increments of 1mL of sodium thiosulfate at a time. Swirl the mixture after each addition. You will have reached the equivalency point when the solution becomes colorless. Try placing the flask containing the analyte on a piece of paper. The white background will help see the color change when it occurs.
- ☐ *Step 5:* – Write down the volume of thiosulfate added to reach the equivalency point to the nearest 0.01mL. This will give an estimate of the volume needed to reach the endpoint.

Part D: Titration of the saturated solution mixtures: exact titration

- ☐ *Step 1:* – After having estimated the volume of titrate needed to reach the endpoint, we will now speed up the titration by rapidly titrating the second thiosulfate sample.
- ☐ *Step 2:* – Add thiosulfate to a second sample of saturated solution, stopping 3mL before the estimated endpoint.
- ☐ *Step 3:* – Add $\text{Na}_2\text{S}_2\text{O}_3$ drop by drop until you reach the endpoint and the solution becomes colorless. Record the buret reading to the nearest 0.01mL.
- ☐ *Step 4:* – Repeat the procedure with the last saturated sample.

Calculations

0 This is the molarity of the potassium iodate solution, c_{KIO_3} .

1 This is the initial volume on the buret.

2 This is the final volume in the buret.

3 This is the volume of sodium thiosulfate used in the standardization, $v_{\text{Na}_2\text{S}_3\text{O}_3}$:

$$v_{\text{Na}_2\text{S}_3\text{O}_3} = 2 - 1$$

4 This is the molarity of thiosulfate, $c_{\text{Na}_2\text{S}_3\text{O}_3}$:

$$c_{\text{Na}_2\text{S}_3\text{O}_3} = 10 \times c_{\text{KIO}_3}^{\text{diluted}} / v_{\text{Na}_2\text{S}_3\text{O}_3} \times 6 = 10 \times 0.05 \times 0 / 3 \times 6$$

5 This is the average sodium thiosulfate concentration. If you do not standardize the reagent this value should be written in the bottle (0.025M).

6 This is the initial volume in the buret.

7 This is the final volume in the buret.

8 This is the volume of thiosulfate used, $v_{\text{Na}_2\text{S}_3\text{O}_3}$:

$$v_{\text{Na}_2\text{S}_3\text{O}_3} = 7 - 6$$

9 This is the concentration of iodate IO_3^- in the saturated solution, $c_{\text{IO}_3^-}$:

$$c_{\text{IO}_3^-} = v_{\text{Na}_2\text{S}_3\text{O}_3} \times c_{\text{Na}_2\text{S}_3\text{O}_3} / 60 = 8 \times 5 / 60$$

10 This is the solubility of the iodate salt, s .

$$s = c_{\text{IO}_3^-} / 2 = 9 / 2$$

11 This is the average solubility.

12 You obtain the solubility product by doing:

$$s = 4 \times 11^3$$

STUDENT INFO


Name:

Date:

Pre-lab Questions


A Solubility Product Constant: Solubility product of an iodate salt

1. In the experiment, you will react potassium iodate with calcium nitrate to produce a precipitate. Indicate the name of the precipitate formed:




This area in the pdf is writable only with Acrobat Reader...

2. In this experiment the iodate (IO_3^-) and iodide ions (I^-) react to produce molecular iodine (I_2) in acidic medium. Hence, in this redox reaction, a reactant is converted into products by simultaneous oxidation and reduction reactions. This type of redox reaction is referred to as disproportionation. Balance the corresponding disproportionation reaction.




This area in the pdf is writable only with Acrobat Reader...

3. In this experiment thiosulfate ($\text{S}_2\text{O}_3^{2-}$) is used in a redox titration. This compound reacts with molecular iodine (I_2), which has a reddish-brown color, to produce the colorless iodide (I^-) and tetrathionate. Hence, when all thiosulfate is finished along with the titration the reaction becomes colored. Balance the corresponding redox reaction in acidic medium.




This area in the pdf is writable only with Acrobat Reader...

4. Both reactions in the previous questions happen consecutively during the experiment. In the redox titration, first iodate reacts with iodine to produce iodine, which also reacts with thiosulfate to produce iodine. Combine both reactions to generate the overall reaction happening in the titration.



This area in the pdf is writable only with Acrobat Reader...

5. You need 25mL of sodium thiosulfate 0.01M to reach the equivalency point in the titration of a 10mL sample of iodate. Using the balanced reaction from the previous question, carry stoichiometric calculations to calculate the molarity of iodate in the solution.



This area in the pdf is writable only with Acrobat Reader...

STUDENT INFO

Name:

Date:

**Results
EXPERIMENT**

A Solubility Product Constant: Solubility product of an iodate salt

Part B: $\text{Na}_2\text{S}_2\text{O}_3$ standarization

This table in the pdf is writable only with Acrobat Reader...

① Molarity of KIO_3 , c_{KIO_3} =

	1	2	3	4
① Initial Buret Volume (mL)				
② Final Buret Volume (mL)				
③ $\text{Na}_2\text{S}_2\text{O}_3$ Volume used (mL)				
④ $\text{Na}_2\text{S}_2\text{O}_3$ Concentration (M)				

⑤ Mean $\text{Na}_2\text{S}_2\text{O}_3$ Concentration (M) =

Part D: Titration of saturated solution

This table in the pdf is writable only with Acrobat Reader...

⑤ Molarity of $\text{Na}_2\text{S}_2\text{O}_3$, $c_{\text{Na}_2\text{S}_2\text{O}_3} =$

		1	2	3	4
⑥	Initial Buret Volume (mL)				
⑦	Final Buret Volume (mL)				
⑧	$\text{Na}_2\text{S}_2\text{O}_3$ Volume used (mL)				
⑨	IO_3^- Concentration (M)				
⑩	Solubility (M)				

Unknown Number = _____

⑪ Mean solubility (M) =

⑫ Mean $K_{sp} =$

STUDENT INFO

Name:

Date:

Post-lab Questions

A Solubility Product Constant: Solubility product of an iodate salt

1. Write down the solubility equilibrium of calcium iodate.

This area in the pdf is writable only with Acrobat Reader...

2. Write down the expression of the solubility constant of calcium iodate in terms of ion molarities.

This area in the pdf is writable only with Acrobat Reader...

3. Write down the expression of the solubility constant of calcium iodate in terms of solubilities.

This area in the pdf is writable only with Acrobat Reader...

4. Write down the expression of the solubility of calcium iodate in terms of the solubility product constant.

This area in the pdf is writable only with Acrobat Reader...

13 Electrochemistry: galvanic cells

EXPERIMENT

Electrochemistry: galvanic cells

Goal

The goal of this laboratory experiment is to build up a series of galvanic cell as well as a concentration cell.

Materials

- ☐ A series of 150mL beaker
- ☐ A salt bridge
- ☐ A series of electrodes: Zn, 2 Cu electrodes, Pb, Fe, 2 graphite electrodes, Pt
- ☐ Two sets of banana plug to alligator clip wires and a voltmeter
- ☐ A series of solutions: 0.1M $\text{Cu}(\text{NO}_3)_2$, 0.1M $\text{Zn}(\text{NO}_3)_2$, 6M HNO_3 , 0.1M $\text{Pb}(\text{NO}_3)_2$, 0.05M $\text{Fe}_2(\text{SO}_4)_3$ mixed with 0.05M FeSO_4 in 1M H_2SO_4

Background

Introduction to galvanic cells

Galvanic cells—also known as voltaic cells or piles—are electrochemical cells that generate electricity from spontaneous redox reactions. They contain two different metals immersed in electrolyte solutions. Each half-cell contains a metal in contact with a liquid solution of the same metal in ionic form. Both half-cells are either connected by a salt bridge or separated by a porous membrane. Galvanic cells differ from batteries. Batteries are composed of multiple single cells working out of equilibrium while producing electricity. Galvanic cells are textbook batteries. They are reversible devices unable to produce electricity.

Components of a galvanic cell

Galvanic cells are composed of two different electrodes, an anode, and a cathode, connected by means of a salt bridge or a membrane. The role of the salt bridge or membrane is to complete the electrical circuit. Anodes are sources of electrons, whereas cathodes are electron sinks. At the same time, anodes generate positive cations, whereas cathodes generate negative anions. The role of the salt bridge or membrane is also to allow the charges generated in the cathode to be compensated by the charges generated in the anode. Electrodes contain two different redox states of the same element in contact with each other. An example of an electrode would be a piece of metallic copper Cu in contact with a solution of $\text{Cu}^{2+}_{(\text{aq})}$ ions. However, electrodes are not always made of metals. For example, electrodes can contain gas in contact (e.g. $\text{H}_{2(\text{g})}$) with an electrolyte solution (e.g. $\text{H}^{+}_{(\text{aq})}$). Electrodes without a metal being directly involved in the redox reaction need to include an external metal to support the charge transfer. Metals such as Pt are normally used for this purpose.

The electrodes: anode and cathode

Every galvanic cell is composed of two electrodes: an anode and a cathode. Electrodes produce ionic and electronic charges. The oxidation occurs on the anode which is indicated with a negative (−) sign. Electrons are being produced in the anode resulting from an oxidation reaction. The reduction occurs on the cathode, indicated with a positive (+) sign. Electrons are being consumed in the cathode resulting from a reduction reaction. These electrodes also produce ions, in particular cations and anions which have a tendency to migrate inside the cell. Anodes generate cations (and consume anions, depending on

the chemical reaction involved), whereas cathodes generate anions (and consume cations, again, depending on the chemical reaction involved). As such, anions have a tendency to migrate to the anode, as their concentration is lower there, whereas cations have a tendency to migrate to the cathode. The excess ionic charge is compensated in the interface between the electrode and the salt bridge or the membrane. Mind that in a galvanic cell only electrons flow through the circuit, by means of the wire connecting both electrodes. The ions involved have a tendency to migrate inside each of the electrodes but do not leave the electrode. The name of the electrodes—anode and cathode—results from the ionic flow involved in the galvanic cell: anions have a tendency to migrate towards the anode and cations to the cathode.

Cell potential

Water flows down a waterfall due to the difference of potential energy between the high and low parts of the waterfall. Similarly, heat flows between a hot and a cold reservoir due to the difference in temperature between both locations. The force that drives the flow is heat is temperature. Electricity flows through a galvanic cell resulting from the difference of cell potential $\Delta\mathcal{E}$ between both electrodes, the anode, and the cathode. The cell potential—also referred to as cell voltage, cell electromotive force, or cell *emf*—is the force that drives the flow of electrons. Anodes and cathodes have a characteristic cell potential associated with the electrochemical half-reaction happening in the electrode. The voltage of the anode ($\mathcal{E}_{\text{anode}}$) is always lower than the one from the cathode ($\mathcal{E}_{\text{cathode}}$). The combination of the anodic and cathodic voltage gives the overall cell potential measured in a galvanic cell. In particular, the overall voltage results from the voltage of the cathode with respect to the anode, so that the overall voltage of a galvanic cell is always positive.

A galvanic cell example

Below we display a representation of the Daniell cell, a classical galvanic cell in which copper is oxidized by zinc. The name of the cell is in honor of John Daniell, a British chemist from the nineteenth century who was trying to develop an electric power supply to sustain telegraphy. He connected with a metallic wire a zinc electrode in contact with a zinc sulfate solution to a copper electrode in contact with a copper (II) sulfate solution. In this cell, Copper(II) ions are converted into metallic copper in the cathode by means of the reaction $\text{Cu}_{(\text{aq})}^{2+} + 2\text{e}^- \longrightarrow \text{Cu}_{(\text{s})}$, whereas metallic zinc is converted into zinc ions in the anode following the reaction $\text{Zn}_{(\text{s})} \longrightarrow \text{Zn}_{(\text{aq})}^{2+} + 2\text{e}^-$. Copper(II) ions are being reduced whereas zinc is being oxidized. The cathode solution becomes negatively charged whereas the anode solution becomes positively charged. A porous membrane permeable to ions was used to avoid the charge buildup.




Figure 10 The Daniell galvanic cell

Anodes and cathodes

The standard reduction potential of an electrode predicts the tendency of an electrode to act as an anode or cathode. Imagine we have two different electrodes involving the following reaction with the standard potentials indicated on the side, and we need to set up a galvanic cell:



How to determine which electrode will act as anode and which will act as cathode? The rule is the smaller the electrode potential the more tendency of the electrode to act as an anode. If we compare the copper and zinc electrodes, as the electrode potential of zinc is smaller—more negative—than the electrode potential of copper. Therefore, zinc will act as an

anode and copper will act as a cathode. The standard reduction potential is not affected by the stoichiometry of the reaction. For example, the standard reduction potential of cesium is -3.03V :



At the same time, the standard reduction potential of two moles of cesium is still -3.03V .



Differently, inverting the reduction reaction switches the sign of the reduction potential. For example, the standard reduction potential of Titanium(II) is -1.63V




while the standard potential of oxidation of Titanium into Titanium(II) is $+1.63\text{V}$



The Nernst equation applied to a galvanic cell


We can also apply the Nernst's equation to a galvanic cell, for example the Daniell cell:



The redox reaction behind this cell oxidizes zinc and reduces copper, with a redox standard potential of 1.10V , with all concentrations fixed at 1M :



For nonstandard conditions we have that as Zinc is being oxidized, its ionic concentration will increase, whereas as Copper is being reduced its concentration will decrease. One can envision the following galvanic cell,



in which the following redox reaction occurs



Applying Nernst's equation at 298.15K we have:

$$\Delta\mathcal{E}_{\text{cell}} = \Delta\mathcal{E}_{\text{cell}}^\circ - \frac{R \cdot T}{n_e \cdot F} \ln Q$$


$$\Delta\mathcal{E}_{\text{cell}} = 1.10 - \frac{8.314 \cdot 298.15}{2 \cdot 96485} \ln \frac{1.2}{0.8} = 1.09\text{V}$$

We have that as the galvanic cell is being consumed, the overall cell potential has decreased. For overall redox reactions behind galvanic cells involving a single ionic species in the anode and cathode, the anode molarity will be on the top side of Q whereas the cathode molarity will be on the bottom:

$$Q = \frac{\text{anode}}{\text{cathode}}$$


Example

For the galvanic cell below:



(a) Calculate the standard voltage of the cell, given $\mathcal{E}^\circ (\text{Au}^+/\text{Au}) = 1.83\text{V}$ and $\mathcal{E}^\circ (\text{Mg}^{2+}/\text{Mg}) = -2.37\text{V}$ (b) calculate the cell potential at 298.15K

Answer: The redox reaction being this galvanic cell is indicated below with the standard cell voltage obtained by doing $\mathcal{E}_{\text{cathode}}^\circ - \mathcal{E}_{\text{anode}}^\circ$:




We have that two electrons flow through the cell. To calculate the cell potential at 298.15K we apply Nernst equation. While calculating Q we have to place the anodic concentration on top and the cathodic on the bottom:

$$\Delta \mathcal{E}_{cell} = \Delta \mathcal{E}_{cell}^{\circ} - \frac{R \cdot T}{n_e \cdot F} \ln Q = 4.20 - \frac{8.314 \cdot 298.15}{2 \cdot 96485} \ln \frac{10^{-4}}{(10^{-3})^2} = 4.14 \text{ V}$$

Concentration cells

What caused a voltage difference in a galvanic cell is the differences in reduction potentials of the anode and cathode. By selecting cathodes with larger positive \mathcal{E}° and anodes with large negative \mathcal{E}° we can achieve large cell potentials. At the same time, the Nernst equations tell us that \mathcal{E} depends on molarity, hence we could potentially select the same anode and cathode and use the molarity terms in the Nernst equation to drive the potential difference. This is the principle behind concentration cells. These types of galvanic cells contain the same material in the anode and cathode but with different concentrations. The anode is the less concentrated electrode, whereas the cathode is the most concentrated electrode. The concentration difference is what drives the cell potential. Concentration cells in general produce very low voltages, even with very high concentration differences. For example, for the cell:




we have that $\mathcal{E}^{\circ} = +0.29\text{V}$. At the same time, good concentration cells have very low anodic concentrations and very large cathodic concentrations. For concentration cells we have that as both electrodes are made of the same material $\Delta \mathcal{E}^{\circ} = 0$. Hence, Nernst equation for a concentration cell is simplified to:

$$\Delta \mathcal{E}_{cell} = - \frac{R \cdot T}{n_e \cdot F} \ln Q \quad \text{where } Q = \frac{[\text{anode}]}{[\text{cathode}]}$$

Example

For the concentration cell below:



(a) Identify the anode and the cathode (b) calculate the cell potential at 298.15K

Answer: The electrode potential depends on concentration. The smaller the concentration the lower the cell potential. On the other hand, we have that electrodes with lower cell potential have a stronger tendency to act as anodes, whereas electrodes with higher cell potential have a stronger tendency to act as cathodes. The less concentrated electrode on the left is the anode, whereas the more concentrated electrode on the right is the cathode. To calculate the cell voltage we use Nernst equation at 298.15K and given that two electrons flow through the cell:

$$\Delta \mathcal{E}_{cell} = - \frac{R \cdot T}{n_e \cdot F} \ln Q = - \frac{8.314 \cdot 298.15}{2 \cdot 96485} \ln \frac{10^{-8}}{10^8} = 0.47 \text{ V}$$

Table 1 Standard reduction potentials at 298K

Element	Reaction	\mathcal{E}° (V)	Element	Reaction	\mathcal{E}° (V)
↑ Sr	$\text{Sr}^+ + \text{e}^- \rightleftharpoons \text{Sr}_{(\text{s})}$	-4.10	H	$2\text{H}^+ + 2\text{e}^- \rightleftharpoons \text{H}_{2(\text{g})}$	0.00
Ca	$\text{Ca}^+ + \text{e}^- \rightleftharpoons \text{Ca}_{(\text{s})}$	-3.80	Ag	$\text{AgBr}_{(\text{s})} + \text{e}^- \rightleftharpoons \text{Ag}_{(\text{s})} + \text{Br}^-$	+0.07
Li	$\text{Li}^+ + \text{e}^- \rightleftharpoons \text{Li}_{(\text{s})}$	-3.04	S	$\text{S}_4\text{O}_2^{6-} + 2\text{e}^- \rightleftharpoons 2\text{S}_2\text{O}_2 - 3$	+0.08
Cs	$\text{Cs}^+ + \text{e}^- \rightleftharpoons \text{Cs}_{(\text{s})}$	-3.03	N	$\text{N}_{2(\text{g})} + 2\text{H}_2\text{O} + 6\text{H}^+ + 6\text{e}^- \rightleftharpoons 2\text{NH}_4\text{OH}_{(\text{aq})}$	+0.09
Ca	$\text{Ca}(\text{OH})_2 + 2\text{e}^- \rightleftharpoons \text{Ca}_{(\text{s})} + 2\text{OH}^-$	-3.02	Hg	$\text{HgO}_{(\text{s})} + \text{H}_2\text{O} + 2\text{e}^- \rightleftharpoons \text{Hg}_{(\text{l})} + 2\text{OH}^-$	+0.10
Ba	$\text{Ba}(\text{OH})_2 + 2\text{e}^- \rightleftharpoons \text{Ba}_{(\text{s})} + 2\text{OH}^-$	-2.99	C	$\text{C}_{(\text{s})} + 4\text{H}^+ + 4\text{e}^- \rightleftharpoons \text{CH}_{4(\text{g})}$	+0.13
Rb	$\text{Rb}^+ + \text{e}^- \rightleftharpoons \text{Rb}_{(\text{s})}$	-2.98	Sn	$\text{Sn}^{4+} + 2\text{e}^- \rightleftharpoons \text{Sn}_2^{2+}$	+0.15
K	$\text{K}^+ + \text{e}^- \rightleftharpoons \text{K}_{(\text{s})}$	-2.93	Cu	$\text{Cu}^{2+} + \text{e}^- \rightleftharpoons \text{Cu}^+$	+0.159
Ba	$\text{Ba}^{2+} + 2\text{e}^- \rightleftharpoons \text{Ba}_{(\text{s})}$	-2.91	Fe	$3\text{Fe}_2\text{O}_3_{(\text{s})} + 2\text{H}^+ + 2\text{e}^- \rightleftharpoons 2\text{Fe}_3\text{O}_4_{(\text{s})} + \text{H}_2\text{O}$	+0.22
Sr	$\text{Sr}^{2+} + 2\text{e}^- \rightleftharpoons \text{Sr}_{(\text{s})}$	-2.90	Ag	$\text{AgCl}_{(\text{s})} + \text{e}^- \rightleftharpoons \text{Ag}_{(\text{s})} + \text{Cl}^-$	+0.22
Sr	$\text{Sr}(\text{OH})_2 + 2\text{e}^- \rightleftharpoons \text{Sr}_{(\text{s})} + 2\text{OH}^-$	-2.88	Cu	$\text{Cu}^{2+} + 2\text{e}^- \rightleftharpoons \text{Cu}_{(\text{s})}$	+0.34
Ca	$\text{Ca}^{2+} + 2\text{e}^- \rightleftharpoons \text{Ca}_{(\text{s})}$	-2.87	Fe	$\text{Fe}^+ + \text{e}^- \rightleftharpoons \text{Fe}_{(\text{s})}$	+0.40
Li	$\text{Li}^+ + \text{C}_6_{(\text{s})} + \text{e}^- \rightleftharpoons \text{LiC}_6_{(\text{s})}$	-2.84	O	$\text{O}_{2(\text{g})} + 2\text{H}_2\text{O} + 4\text{e}^- \rightleftharpoons 4\text{OH}^-_{(\text{aq})}$	+0.40
Na	$\text{Na}^+ + \text{e}^- \rightleftharpoons \text{Na}_{(\text{s})}$	-2.71	Cu	$\text{Cu}^+ + \text{e}^- \rightleftharpoons \text{Cu}_{(\text{s})}$	+0.520
Mg	$\text{Mg}(\text{OH})_2 + 2\text{e}^- \rightleftharpoons \text{Mg}_{(\text{s})} + 2\text{OH}^-$	-2.69	C	$\text{CO}_{(\text{g})} + 2\text{H}^+ + 2\text{e}^- \rightleftharpoons \text{C}_{(\text{s})} + \text{H}_2\text{O}$	+0.52
Mg	$\text{Mg}^{2+} + 2\text{e}^- \rightleftharpoons \text{Mg}_{(\text{s})}$	-2.37	I	$\text{I}_{2(\text{s})} + 2\text{e}^- \rightleftharpoons 2\text{I}^-$	+0.54
H	$\text{H}_{2(\text{g})} + 2\text{e}^- \rightleftharpoons 2\text{H}^-$	-2.23	Mn	$\text{MnO}_4^- + 2\text{H}_2\text{O} + 3\text{e}^- \rightleftharpoons \text{MnO}_{2(\text{s})} + 4\text{OH}^-$	+0.595
Sr	$\text{Sr}^{2+} + 2\text{e}^- \rightleftharpoons \text{Sr}_{(\text{Hg})}$	-1.79	O	$\text{O}_{2(\text{g})} + 2\text{H}^+ + 2\text{e}^- \rightleftharpoons \text{H}_2\text{O}_{2(\text{aq})}$	+0.70
Al	$\text{Al}^{3+} + 3\text{e}^- \rightleftharpoons \text{Al}_{(\text{s})}$	-1.66	Fe	$\text{Fe}_2\text{O}_3_{(\text{s})} + 6\text{H}^+ + 2\text{e}^- \rightleftharpoons 2\text{Fe}^{2+} + 3\text{H}_2\text{O}$	+0.728
Ti	$\text{Ti}^{2+} + 2\text{e}^- \rightleftharpoons \text{Ti}_{(\text{s})}$	-1.63	Fe	$\text{Fe}^{3+} + \text{e}^- \rightleftharpoons \text{Fe}^{2+}$	+0.77
Ti	$\text{Ti}^{3+} + 3\text{e}^- \rightleftharpoons \text{Ti}_{(\text{s})}$	-1.37	Ag	$\text{Ag}^+ + \text{e}^- \rightleftharpoons \text{Ag}_{(\text{s})}$	+0.80
Ti	$\text{TiO}_{(\text{s})} + 2\text{H}^+ + 2\text{e}^- \rightleftharpoons \text{Ti}_{(\text{s})} + \text{H}_2\text{O}$	-1.31	Hg	$\text{Hg}_2^{2+} + 2\text{e}^- \rightleftharpoons 2\text{Hg}_{(\text{l})}$	+0.80
Mn	$\text{Mn}^{2+} + 2\text{e}^- \rightleftharpoons \text{Mn}_{(\text{s})}$	-1.18	N	$\text{NO}_3^-_{(\text{aq})} + 2\text{H}^+ + \text{e}^- \rightleftharpoons \text{NO}_{2(\text{g})} + \text{H}_2\text{O}$	+0.80
V	$\text{V}^{2+} + 2\text{e}^- \rightleftharpoons \text{V}_{(\text{s})}$	-1.13	Fe	$2\text{FeO}_2^{2-} + 5\text{H}_2\text{O} + 6\text{e}^- \rightleftharpoons \text{Fe}_2\text{O}_3_{(\text{s})} + 10\text{OH}^-$	+0.81
Ti	$\text{TiO}^{2+} + 2\text{H}^+ + 4\text{e}^- \rightleftharpoons \text{Ti}_{(\text{s})} + \text{H}_2\text{O}$	-0.93	Hg	$\text{Hg}^{2+} + 2\text{e}^- \rightleftharpoons \text{Hg}_{(\text{l})}$	+0.85
Si	$\text{SiO}_{2(\text{s})} + 4\text{H}^+ + 4\text{e}^- \rightleftharpoons \text{Si}_{(\text{s})} + 2\text{H}_2\text{O}$	-0.91	Mn	$\text{MnO}_4^- + \text{H}^+ + \text{e}^- \rightleftharpoons \text{HMnO}_4^-$	+0.90
Fe	$\text{Fe}_2\text{O}_3_{(\text{s})} + 3\text{H}_2\text{O} + 2\text{e}^- \rightleftharpoons 2\text{Fe}(\text{OH})_{2(\text{s})} + 2\text{OH}^-$	-0.86	Hg	$2\text{Hg}^{2+} + 2\text{e}^- \rightleftharpoons \text{Hg}_2^{2+}$	+0.91
H	$2\text{H}_2\text{O} + 2\text{e}^- \rightleftharpoons \text{H}_{2(\text{g})} + 2\text{OH}^-$	-0.828	Pd	$\text{Pd}^{2+} + 2\text{e}^- \rightleftharpoons \text{Pd}_{(\text{s})}$	+0.915
Zn	$\text{Zn}^{2+} + 2\text{e}^- \rightleftharpoons \text{Zn}_{(\text{s})}$	-0.762	N	$\text{NO}_3^-_{(\text{aq})} + 4\text{H}^+ + 3\text{e}^- \rightleftharpoons \text{NO}_{(\text{g})} + 2\text{H}_2\text{O}_{(\text{l})}$	+0.96
Cr	$\text{Cr}^{3+} + 3\text{e}^- \rightleftharpoons \text{Cr}_{(\text{s})}$	-0.74	Fe	$\text{Fe}_3\text{O}_4_{(\text{s})} + 8\text{H}^+ + 2\text{e}^- \rightleftharpoons 3\text{Fe}^{2+} + 4\text{H}_2\text{O}$	+0.98
Ni	$\text{Ni}(\text{OH})_{2(\text{s})} + 2\text{e}^- \rightleftharpoons \text{Ni}_{(\text{s})} + 2\text{OH}^-$	-0.72	Br	$\text{Br}_{2(\text{aq})} + 2\text{e}^- \rightleftharpoons 2\text{Br}^-$	+1.09
Ag	$\text{Ag}_2\text{S}_{(\text{s})} + 2\text{e}^- \rightleftharpoons 2\text{Ag}_{(\text{s})} + \text{S}_{2-}^2_{(\text{aq})}$	-0.69	Ag	$\text{Ag}_2\text{O}_{(\text{s})} + 2\text{H}^+ + 2\text{e}^- \rightleftharpoons 2\text{Ag}_{(\text{s})} + \text{H}_2\text{O}$	+1.17
Pb	$\text{PbO}_{(\text{s})} + \text{H}_2\text{O} + 2\text{e}^- \rightleftharpoons \text{Pb}_{(\text{s})} + 2\text{OH}^-$	-0.58	Pt	$\text{Pt}^{2+} + 2\text{e}^- \rightleftharpoons \text{Pt}_{(\text{s})}$	+1.188
Fe	$\text{Fe}^{2+} + 2\text{e}^- \rightleftharpoons \text{Fe}_{(\text{s})}$	-0.44	Cl	$\text{ClO}^- + 2\text{H}^+ + 2\text{e}^- \rightleftharpoons \text{ClO}_3^- + \text{H}_2\text{O}$	+1.20
Cr	$\text{Cr}^{3+} + \text{e}^- \rightleftharpoons \text{Cr}_2^+$	-0.42	O	$\text{O}_{2(\text{g})} + 4\text{H}^+ + 4\text{e}^- \rightleftharpoons 2\text{H}_2\text{O}$	+1.229
Cd	$\text{Cd}^{2+} + 2\text{e}^- \rightleftharpoons \text{Cd}_{(\text{s})}$	-0.40	Cl	$\text{Cl}_{2(\text{g})} + 2\text{e}^- \rightleftharpoons 2\text{Cl}^-$	+1.36
Cu	$\text{Cu}_2\text{O}_{(\text{s})} + \text{H}_2\text{O} + 2\text{e}^- \rightleftharpoons 2\text{Cu}_{(\text{s})} + 2\text{OH}^-$	-0.36	Br	$\text{BrO}_3^- + 5\text{H}^+ + 4\text{e}^- \rightleftharpoons \text{HBrO}_{(\text{aq})} + 2\text{H}_2\text{O}$	+1.45
Pb	$\text{PbSO}_4_{(\text{s})} + 2\text{e}^- \rightleftharpoons \text{Pb}_{(\text{s})} + \text{SO}_2 - 4$	-0.36	Br	$2\text{BrO}_3^- + 12\text{H}^+ + 10\text{e}^- \rightleftharpoons \text{Br}_2_{(\text{l})} + 6\text{H}_2\text{O}$	+1.48
Pb	$\text{PbSO}_4_{(\text{s})} + 2\text{e}^- \rightleftharpoons \text{Pb}(\text{Hg}) + \text{SO}_2 - 4$	-0.35	Cl	$2\text{ClO}_3^- + 12\text{H}^+ + 10\text{e}^- \rightleftharpoons \text{Cl}_{2(\text{g})} + 6\text{H}_2\text{O}$	+1.49
Co	$\text{Co}^{2+} + 2\text{e}^- \rightleftharpoons \text{Co}_{(\text{s})}$	-0.28	Mn	$\text{MnO}_4^- + 8\text{H}^+ + 5\text{e}^- \rightleftharpoons \text{Mn}_2^+ + 4\text{H}_2\text{O}$	+1.51
Ni	$\text{Ni}^{2+} + 2\text{e}^- \rightleftharpoons \text{Ni}_{(\text{s})}$	-0.25	Au	$\text{Au}^{3+} + 3\text{e}^- \rightleftharpoons \text{Au}_{(\text{s})}$	+1.52
As	$\text{As}_{(\text{s})} + 3\text{H}^+ + 3\text{e}^- \rightleftharpoons \text{AsH}_3_{(\text{g})}$	-0.23	Pb	$\text{Pb}^{4+} + 2\text{e}^- \rightleftharpoons \text{Pb}^{2+}$	+1.69
Ag	$\text{AgI}_{(\text{s})} + \text{e}^- \rightleftharpoons \text{Ag}_{(\text{s})} + \text{I}^-$	-0.15	Mn	$\text{MnO}_4^- + 4\text{H}^+ + 3\text{e}^- \rightleftharpoons \text{MnO}_{2(\text{s})} + 2\text{H}_2\text{O}$	+1.70
Sn	$\text{Sn}^{2+} + 2\text{e}^- \rightleftharpoons \text{Sn}_{(\text{s})}$	-0.13	Ag	$\text{AgO}_{(\text{s})} + 2\text{H}^+ + \text{e}^- \rightleftharpoons \text{Ag}^+ + \text{H}_2\text{O}$	+1.77
Pb	$\text{Pb}^{2+} + 2\text{e}^- \rightleftharpoons \text{Pb}_{(\text{s})}$	-0.126	O	$\text{H}_2\text{O}_{2(\text{aq})} + 2\text{H}^+ + 2\text{e}^- \rightleftharpoons 2\text{H}_2\text{O}$	+1.78
C	$\text{CO}_{2(\text{g})} + 2\text{H}^+ + 2\text{e}^- \rightleftharpoons \text{HCOOH}_{(\text{aq})}$	-0.11	Au	$\text{Au}^+ + \text{e}^- \rightleftharpoons \text{Au}_{(\text{s})}$	+1.83
C	$\text{CO}_{2(\text{g})} + 2\text{H}^+ + 2\text{e}^- \rightleftharpoons \text{CO}_{(\text{g})} + \text{H}_2\text{O}$	-0.11	Ag	$\text{Ag}^{2+} + \text{e}^- \rightleftharpoons \text{Ag}^+$	+1.98
Fe	$\text{Fe}_3\text{O}_4_{(\text{s})} + 8\text{H}^+ + 8\text{e}^- \rightleftharpoons 3\text{Fe}_{(\text{s})} + 4\text{H}_2\text{O}$	-0.08	Mn	$\text{HMnO}_4^- + 3\text{H}^+ + 2\text{e}^- \rightleftharpoons \text{MnO}_{2(\text{s})} + 2\text{H}_2\text{O}$	+2.09
Fe	$\text{Fe}^{3+} + 3\text{e}^- \rightleftharpoons \text{Fe}_{(\text{s})}$	-0.04	Fe	$\text{FeO}_2^{4-} + 8\text{H}^+ + 3\text{e}^- \rightleftharpoons \text{Fe}_3^+ + 4\text{H}_2\text{O}$	+2.20
H	$2\text{H}^+ + 2\text{e}^- \rightleftharpoons \text{H}_{2(\text{g})}$	0.00	F	$\text{F}_{2(\text{g})} + 2\text{H}^+ + 2\text{e}^- \rightleftharpoons 2\text{HF}_{(\text{aq})}$	+3.05


Increasing oxidizing strength (decreasing reducing strength)

Procedure

Part A: Daniell galvanic cell

- ☐ *Step 1:* – Obtain a copper and zinc electrodes, a salt bridge as well as a potentiometer.
- ☐ *Step 2:* – We use a BK Precision potentiometer: make sure the black wire is connected to the COM connection whereas the red one is connected to the Hz V Ω connection. The potentiometer has to be turned towards the right to start measuring DC. Make sure you select the scale that gives Volts.
- ☐ *Step 3:* – Place 0.1M $\text{Cu}(\text{NO}_3)_2$ (approximately 50 mL) in one compartment of the cell and 0.1M $\text{Zn}(\text{NO}_3)_2$ in the other.
- ☐ *Step 4:* – Obtain a copper and zinc electrodes.
- ☐ *Step 5:* – Clean the electrodes by dipping them (no more than 3 seconds) into 40mL of 6M HNO_3 . You will notice that the surface of the electrodes will become bright. Rinse the electrodes with distilled water.

Good Lab Practice


 Concentrated acid can cause severe burns and ruin your clothes.

- ☐ *Step 6:* – Obtain a salt bridge as well as a potentiometer.
- ☐ *Step 7:* – Place the electrodes in separate 150-mL beakers, placing the salt bridge between both beakers. Connect each electrode to a terminal of the potentiometer.
- ☐ *Step 8:* – To find out whether the electrodes are properly connected to the potentiometer, read the voltage from the potentiometer. If the number is negative you need to switch the electrode connections.
- ☐ *Step 9:* – Record the cell potential in the results section, as well as the nature of both electrodes (Zn and Cu).

Part B: A series of galvanic cell

- ☐ *Step 1:* – Obtain a platinum (or graphite), copper, iron, lead and zinc electrodes.
- ☐ *Step 2:* – Clean the electrodes by dipping them (no more than 3 seconds) into 40mL of 6M HNO_3 . You will notice that the surface of the electrodes will become bright. Rinse the electrodes with distilled water.

Good Lab Practice

 Concentrated acid can cause severe burns and ruin your clothes.

 Use gloves to handle lead solutions. Lead solutions are poisonous.

- ☐ *Step 3:* – Gather the following solutions: 0.1M $\text{Cu}(\text{NO}_3)_2$, $\text{Zn}(\text{NO}_3)_2$, $\text{Pb}(\text{NO}_3)_2$, FeSO_4 , and an iron solution made of 0.05M FeSO_4 and $\text{Fe}_2(\text{SO}_4)_3$ in 1M sulfuric acid (this last solution was made by dissolving $\text{Fe}_2(\text{SO}_4)_3$ in sulfuric acid as solvent).
- ☐ *Step 4:* – Set up the following galvanic cells, writing down in the results section its voltage and the nature of the anode and cathode. Mind you can replace Pt by graphite.

- ☐ $\text{Cu} \mid \text{Cu}^{2+}(0.1 \text{ M})$ and $\text{Pb} \mid \text{Pb}^{2+}(0.1 \text{ M})$
- ☐ $\text{Cu} \mid \text{Cu}^{2+}(0.1 \text{ M})$ and $\text{Fe} \mid \text{Fe}^{2+}(0.1 \text{ M})$
- ☐ $\text{Cu} \mid \text{Cu}^{2+}(0.1 \text{ M})$ and $\text{Pt} \mid \text{Fe}^{2+}, \text{Fe}^{3+}(0.1 \text{ M})$
- ☐ $\text{Zn} \mid \text{Zn}^{2+}(0.1 \text{ M})$ and $\text{Pb} \mid \text{Pb}^{2+}(0.1 \text{ M})$


- ☐ $\text{Zn} \mid \text{Zn}^{2+}(0.1 \text{ M})$ and $\text{Fe} \mid \text{Fe}^{2+}(0.1 \text{ M})$
- ☐ $\text{Zn} \mid \text{Zn}^{2+}(0.1 \text{ M})$ and $\text{Pt(or graphite)} \mid \text{Fe}^{2+}, \text{Fe}^{3+}(0.1 \text{ M})$
- ☐ $\text{Fe} \mid \text{Fe}^{2+}(0.1 \text{ M})$ and $\text{Pt} \mid \text{Fe}^{2+}, \text{Fe}^{3+}(0.1 \text{ M})$

- ☐ *Step 5:* – For each cell, write down the measured voltage in the Results section. Every time you change electrodes, make sure you rinse the salt bridge with water.

Part C: Concentration cell

- ☐ *Step 1:* – Obtain two copper electrodes, a salt bridge as well as a potentiometer.
- ☐ *Step 2:* – Clean the electrodes by dipping them (no more than 3 seconds) into 40mL of 6M HNO_3 . You will notice that the surface of the electrodes will become bright. Rinse the electrodes with distilled water.

Good Lab Practice

 Concentrated acid can cause severe burns and ruin your clothes.

- ☐ *Step 3:* – Prepare a copper diluted solution by adding 1mL of a 0.1M concentrated $\text{Cu}(\text{NO}_3)_2$ solution into a 100mL graduated cylinder. Stir well the mixture. Fill the cylinder with distilled water until you reach the 100mL mark.
- ☐ *Step 4:* – Place some of the concentrated solution in a 150-mL beaker and some of the diluted solution in the other one. Insert one copper electrode in each beaker, placing the salt bridge between both beakers. Connect each electrode to a terminal of the potentiometer.
- ☐ *Step 5:* – Record the cell potential in the results section, as well as the nature of both electrodes ($\text{Cu} \mid \text{Cu}^{2+}(0.1 \text{ M})$ and $\text{Cu} \mid \text{Cu}^{2+}(0.01 \text{ M})$) while identifying the anode and cathode.
- ☐ *Step 6:* – Add 10 drops of the concentrated 0.1M $\text{Cu}(\text{NO}_3)_2$ solution to the most diluted electrode of the cell you built up. Record the cell potential in the results section.

Part D: Electrolysis

- ☐ *Step 1:* – Obtain two graphite electrodes and a power supply (6V battery).
- ☐ *Step 2:* – Weigh 2g of KI and transfer the solid into a 250mL Erlenmeyer flask. Add 100mL of distilled water and swirl until the solid dissolves.
- ☐ *Step 3:* – Using a PH meter, measure the PH of the KI solution. Write down the value in the results section.
- ☐ *Step 4:* – Place two graphite electrodes in a 250mL beaker. Clamp the electrodes so that they don't touch each other. Use a pair of upside-down test tubes to secure the electrodes. Attach the wires to the electrodes but do not connect them yet to the battery.
- ☐ *Step 5:* – Add 30mL of the KI solution to the beaker.
- ☐ *Step 6:* – Attach the wires to the battery and record the time in the Results section. Observe the solution.
- ☐ *Step 7:* – Allow the electrolysis to run for 15 minutes.
- ☐ *Step 8:* – After the time has elapsed, remove the electrodes from the battery. Remove the electrodes and take a good look at them.
- ☐ *Step 9:* – Stir the KI solution in the beaker and using a PH meter measure the PH. Write down the results in the results section.

Calculations

- ① Write down here the measured voltage.
- ② Write down here the cathodic reaction.
- ③ Write down here the cathodic standard potential.
- ④ Write down here the anodic reaction (revert the reaction from the table of voltages).
- ⑤ Write down here the anodic standard potential.
- ⑥ Write down here the overall reaction: add the cathodic reaction and the anodic reaction.
- ⑦ Write down here the overall standard potential.
- ⑧ Elapsed time in seconds.
- ⑨ Measured PH.
- ⑩ Calculated hydroxyl concentration:

$$[\text{OH}^-] = 1 \times 10^{-14} / 10^{-PH}$$

- ⑪ Change on hydroxyl concentration:

$$\Delta [\text{OH}^-] = [\text{OH}^-]_f - [\text{OH}^-]_0$$

- ⑫ The intensity in Amperes:

$$\textcircled{11} \times 96485 / \textcircled{8}$$

STUDENT INFO


Name:

Date:


Pre-lab Questions

Electrochemistry: galvanic cells

1. For the galvanic cell below, indicate: (a) the direction of flow of electrons (b) the direction of flow of cations (c) the direction of flow of anions




2. Calculate the standard cell potential for the galvanic cell below:



3. For the unbalanced reaction below, calculate \mathcal{E}° and indicate whether the reaction is spontaneous under standard conditions. Balance the reaction.



4. Calculate the standard cell potential for the galvanic cell below:



5. For the concentration cell below:



- (a) Identify the anode and the cathode (b) calculate the cell potential at 298.15K

STUDENT INFO

Name:

Date:

Results
EXPERIMENTElectrochemistry:
galvanic cells

Part A						
Electrodes	$\Delta \mathcal{E}_{exp}$	Cathodic reaction	$\mathcal{E}_{cath}^{\circ}$	Anodic reaction	$\mathcal{E}_{anod}^{\circ}$	Total reaction
	①	②	③	④	⑤	⑥
Cu—Zn						⑦
Part B						
Electrodes	$\Delta \mathcal{E}_{exp}$	Cathodic reaction	$\mathcal{E}_{cath}^{\circ}$	Anodic reaction	$\mathcal{E}_{anod}^{\circ}$	Total reaction
Cu—Pb						
Cu—Fe						
Zn—Pb						
Cu—Fe ³⁺ +Fe ²⁺						

Part B						
Electrodes	$\Delta \mathcal{E}_{exp}$	Cathodic reaction	$\mathcal{E}^{\circ}_{cath}$	Anodic reaction	$\mathcal{E}^{\circ}_{anod}$	Total reaction
	①	②	③	④	⑤	⑥
Zn—Fe ³⁺ +Fe ²⁺						
Fe, Fe ²⁺ —Fe ³⁺ +Fe ²⁺						
Cu ²⁺ (0.01 M)—Cu ²⁺ (0.1 M)						
Part C						
Cu ²⁺ (0.01 M)*—Cu ²⁺ (0.1 M)						

* After adding 10 drops of Cu²⁺ (0.1 M)

** Use the temperature value of 298K

Part D					
pH_0	pH_f	[OH ⁻] ₀	[OH ⁻] _f	Δ[OH ⁻]	I
⑨	⑨	⑩	⑩	⑪	⑫

STUDENT INFO

Name:

Date:

Post-lab Questions

Electrochemistry: galvanic cells

1. Calculate the percent error between the measured voltage and the calculated one for all electrodes you measured using the formula below:

$$\%Error = \frac{|\Delta\mathcal{E}_{calc} - \Delta\mathcal{E}_{exp}|}{\Delta\mathcal{E}_{calc}} \times 100$$

Electrodes	$\Delta\mathcal{E}_{exp}$ (1)	$\Delta\mathcal{E}^\circ$ (7)	Expression for $-\frac{R \cdot T}{n_e \cdot F} \ln Q$ (8)	$\Delta\mathcal{E}_{calc}$ (7)+(8)	% Error
Cu—Zn					
Cu—Pb					
Cu—Fe					
Zn—Pb					
Cu—Fe ³⁺ +Fe ²⁺					
Zn—Fe ³⁺ +Fe ²⁺					
Fe, Fe ²⁺ —Fe ³⁺ +Fe ²⁺					
Cu ²⁺ (0.01 M)—Cu ²⁺ (0.1 M)					
Cu ²⁺ (0.01 M)*—Cu ²⁺ (0.1 M)					

* After adding 10 drops of Cu²⁺ (0.1 M)

2. Using line notation, write down the cell with that gives the largest and the cell that gives the smallest voltage.

14 Radioactivity: half-life of K-40

EXPERIMENT

Radioactivity: half-life of K-40

Goal

The goal of this laboratory experiment is to calculate the half-life of potassium-40 by means of the measurement of its rate of disintegration.

Materials

- ☐ A Geiger counter
- ☐ Potassium chloride

- ☐ A 100mL beaker
- ☐ A stand and a clamp

Background

Light elements have normally stable nuclei. Differently, heavier elements with atomic numbers larger than 20 tend to often have several isotopes—remember these are atoms of the element with a different number of neutrons—that have unstable nuclei. For these unstable isotopes, the forces that keep the nucleus together are not strong enough to stabilize the nuclei. An unstable nucleus is radioactive, which means that it will spontaneously emit radiation in the form of small particles. Not all radioactivity is the same and there exist different types of radiation, which we will address in the following. Table 2 reports common nuclear symbols.

alpha radiation

Alpha radiation—referred to as α —is a type of radiation that contains alpha particles. These particles are indeed helium nuclei, with 2 protons, 2 neutrons, and a (2+) positive charge. Alpha particles are often represented as α or ${}^4_2\text{He}$.

beta radiation

Beta radiation—referred to as β —is a type of radiation that contains beta particles. These particles are indeed high-energy electrons with (−) negative charge. Beta particles are often represented as β or ${}^0_{-1}\text{e}$.

gamma radiation

Gamma radiation—referred to as γ —is a type of radiation that contains high-energy photons. These particles are indeed photons with no mass or charge. Gamma particles are often represented as γ or ${}^0_0\gamma$.

protons

Protons in this chapter are often referred to as p or ${}^1_1\text{H}^+$. These are positive charges.

positrons

Positrons are the electron antiparticle, often referred to as β^+ or ${}^0_{+1}\text{e}$. They do have a positive charge.

neutrons

Neutrons are nuclear particles with no charge, often referred to as n or ${}_0^1n$.

Table 2 Nuclear symbols

Particle Name		Symbol	Charge	Identity	Penetrating power	Discovery
Alpha	(α)	${}^4_2\text{He}$	2+	Helium nucleus	Minimal	1899
Beta	(β)	${}^0_{-1}\text{e}$	-1	Electrons	Short	1899
Gamma	(γ)	${}^0_0\gamma$	0	Electromagnetic radiation	Deep	1900
Neutrons	(n)	${}^1_0\text{n}$	0	nuclear particle	Maximal	1932
Proton	(p)	${}^1_1\text{H}^+$	+1	nuclear particle		1919
Positrons	(β^+)	${}^0_{+1}\text{e}$	+1	antiparticle		1932


Isotopes—called emitters—spontaneously decompose producing new isotopes in a process called radioactive decay. In this decay, radiation is also emitted.



In the following, we will discuss the most important type of radioactive decay.


alpha decay

Some isotopes produce alpha radiation, that is, they produce α particles on its decay. A nuclear reaction that produces an α particle (${}_2^4\text{He}$) is called alpha decay. In alpha decay, the emitter decreases its mass number A four units and its atomic number Z two units.




beta decay

Other isotopes produce beta radiation, that is, they produce β particles on its decay. A nuclear reaction that produces a β particle (${}_{-1}^0\text{e}$) is called beta decay. In beta decay, the emitter has the same mass number A as the product isotope. However, its atomic number Z decreases by one unit.



positron emission

Certain isotopes decay by producing a positron, that is, they produce ${}_{+1}^0\text{e}$ particles on its decay. A nuclear reaction that produces ${}_{+1}^0\text{e}$ is called positron emission. In a positron emission, the emitter has the same mass number A as the product isotope. However, its atomic number Z increases by one unit.



gamma decay

Some other isotopes produce gamma radiation in the form of γ particles on its decay. A nuclear reaction that produces a γ particle (${}_0^0\gamma$) is called gamma decay. In this type of decay, no new isotope is produced. Gamma emitters are normally excited, that is they have higher energy than normal; we denote this with a $*$ symbol. Excited particles tend to lose energy to become more stable. In gamma decay, the emitter and the product isotope, both have the same mass and atomic number.




Table 11 Half-life for various isotopes and chemicals

Americium-241	432.2 years	Lutetium-177	6.71 days	Hydrogen-3	12.35 years
Barium-133	10.74 years	Molybdenum-99	66 hours	Technetium-99	213,000 years
Bismuth-212	60.55 minutes	Nickel-63	96 years	Indium-111	2.83 days
Cadmium-109	464 days	Phosphorus-32	14.29 days	Technetium-99m	6.02 hours
Calcium-45	163 days	Potassium-40	1.28×10^9 years	Indium-113m	1.658 hours
Carbon-14	5730 years	Plutonium-239	24,065 years	Tin-113 115.1	days
Cesium-137	30 years	Polonium-210	138.38 days	Iodine-123	13.2 hours
Chlorine-36	301,000 years	Radium-226	1600 years	Tungsten-188	69.4 days
Chromium-51	27.704 days	Radon-222	3.8235 days	Iodine-125	60.14 days
Cobalt-57	270.9 days	Rhenium-188	16.98 hours	Uranium-235	703,800,000 years
Cobalt-58	70.8 days	Rubidium-81	4.58 hours	Iodine-129	15,700,000 years
Cobalt-60	5.271 years	Selenium-75	119.8 days	Uranium-238	4,468,000,000 years
Copper-62	9.74 minutes	Sodium-22	2.602 years	Iodine-131	8.04 days
Copper-64	12.701 hours	Sodium-24	15 hours	Xenon-127	6.41 days
Copper-67	61.86 hours	Strontium-85	64.84 days	Iron-55	2.7 years
Gallium-67	78.26 hours	Strontium-89	50.5 days	Xenon-133	5.245 days
Gold-195	183 days	Sulfur-35	87.44 days	Iron-59	44.529 days
Ondansetron	360 min	Capecitabine	2400s	Carmustine	0.25h

Radioisotopes— isotopes that decay producing radiation—are unstable and with time they eventually disappear given a more stable isotope. Some radioisotopes decay very quickly, such as the ones used in nuclear medicine to fight cancer. Other radioisotopes take longer to disappear.

The concept of half-live

The half-life of an isotope represented as $t_{1/2}$ is the time it takes for an isotope to disappear reducing the sample mass to half the initial value. For example, $t_{1/2}$ for chromium-51 is 28 days and that means that after 28 days a sample of 1 gram of the radioisotope will indeed weigh 0.5 g. Table 11 reports half-lives of numerous isotopes. Samples of radioisotopes weigh less and less with time as they decompose producing more stable isotopes. Similarly, $t_{1/2}$ for strontium-90 is 38 years which means that a one-gram sample will take 38 years to reduce its mass to 0.5g. We can use the concept of half-life to compare the speed of decomposition of different radioisotopes. For example $t_{1/2}$ for strontium-90 is 38 years whereas $t_{1/2}$ for chromium-51 is 28 years. Hence, strontium-90 will exist longer than chromium-51. The activity of an isotope is indeed its rate of the decomposition r which depends on the amount of radioactive isotope you have in the sample n ,

$$r = kn$$

where k is the rate constant for the decomposition. At the same time this rate constant is related to half-life, as decomposition is a first order reaction:

$$t_{1/2} = \frac{0.693}{k}$$

Quantifying half-live

The formula that related the amount of radioisotope with $t_{1/2}$ is:

$$N(t) = N_o \cdot 0.5^{\left(\frac{t}{t_{1/2}}\right)} \quad (8)$$

where $N(t)$ is the amount of isotope at a given time t , N_o is the initial amount of isotope, t is the time and $t_{1/2}$ is the half-life. $N(t)$ is often referred to as the activity of the radioisotope at a given time t . At the same time, while the radioisotope disappears, a new isotope—this time more stable than the radioisotope—starts forming. The amount of product formed $F(t)$

at a given time is:

$$F(t) = N_o \cdot \left[1 - 0.5^{\left(\frac{t}{t_{1/2}}\right)} \right] \quad (9)$$

After several half-lives

So if the half-life is the time it takes for a radioisotope to decompose in half, what would happen after several half-lives? For example, imagine we have 20 grams of iridium-131 with a half-life of 8 days. When we prepare or hypothetically unseal the sample, we will have 20 grams of ^{131}Ir . After one half-life (8 days) we'll have 10 grams of ^{131}Ir . After two half-lives (16 days), we'll have 5 grams of ^{131}Ir . Similarly, after three half-lives (22 days), we'll have 2.5 grams.

Example

$^{131}_{53}\text{I}$ has a half-life of 8 days. How many milligrams of a 50mg sample will remain after 10 days.

Answer: We will following steps to solve this problem:

1 **Step one:** list of the given variables.

2 **Step two:** use the half-life formula $N(t) = N_o \cdot 0.5^{\frac{t}{t_{1/2}}}$ to obtain the mass remaining of the radioisotope

$$N(t) = N_o \cdot 0.5^{\frac{t}{t_{1/2}}}$$

50mg 8d 10d

3 **Step three:** solve for $N = 50 \cdot 0.5^{\frac{10}{8}} = 50 \cdot 0.5^{1.25} = 21\text{mg}$. The result means that after 10 days from the 50mg sample of radioisotope, only 21mg will remain due to radioactive decay.

Procedure

Part A: Background measurement

- ☐ **Step 1:** – Obtain a Geiger counter. Turn it on and let it warm up for five minutes.
- ☐ **Step 2:** – Start the Geiger counter. Set up the measurement time to 60 seconds and the measuring voltage according to your professor's instructions. Mind to select a voltage of 900V for all measurements (Press Display/High Voltage/Up/Down until you reach 900V). Press measure (press Display until the light cursor is next to count; then press Count until the stop button lights up.) and write down the background radioactivity in counts per minute in the table below.
- ☐ **Step 3:** – Measure the background radiation by reading the meter 10 times.
- ☐ **Step 4:** – Compute the average activity and standard deviation.

Part B: Calculation of $t_{1/2}$ for a set of mass measurements

- ☐ Step 1: – Weight between 10 and 11 grams of KCl in a scale. Write down your measurement. Place the sample in a 100mL beaker.
- ☐ Step 2: – With the help of a stand, clamp the counter as close a possible to the sample but without touching the sample. Be very careful with the membrane at the end of the counter as it is very delicate and tears easily.
- ☐ Step 3: – Set up the measurement time to 60 seconds and the measuring voltage according to your professor's instructions. Mind to select a voltage of 900V for all measurements (Press Display/High Voltage/Up/Down until you reach 900V). Press measure (press Display until the light cursor is next to count; then press Count until the stop button lights up.) and write down the activity in counts per minute in the Results section.
- ☐ Step 4: – Measure radiation by reading the meter 10 times.
- ☐ Step 5: – Repeat Part A for a set of K masses between: 14-15g, 19-20g, and 24-25g.
- ☐ Step 6: – For each sample compute the average activity \bar{A} in cps and the number of 40-K atoms in the sample, N .
- ☐ Step 7: – Plot \bar{A} in cps in the vertical axis versus N in the horizontal axis.
- ☐ Step 8: – Compute the half-life for each mass measurement with its average.
- ☐ Step 9: – From the slope of your plot calculate the half-life of the isotope.
- ☐ Step 10: – Compare the calculated half-lives with the experimental value.

Calculations

① The mass of KCl your weighted, m .

② The number of 40-K atoms in you sample, n , given its abundance is 0.012%:

$$N = \frac{m \cdot N_o \cdot f}{AW} = \textcircled{0} \cdot 6.02 \times 10^{23} \cdot 0.00012 \cdot 39.0983 / 74.5$$

③ The average activity in cpm.

④ The average activity in cps:

$$A(cps) = \textcircled{2} / 60$$

⑤ The standard deviation in cps.

⑥ These are the activity values in cpm for the background (without K).

⑦ These are the activity values in cpm for the sample.

⑧ These are the activity values in cpm for the sample without the background. You can susbtract the average value and the standard deviation.

$$\textcircled{6} - \textcircled{5}$$

- 8 Calculated half-life of the isotope in years (the one measured in the experiment). In the formula below you need to use the Activity in cps.

$$t_{1/2}^{calc.} = \frac{0.693 \cdot \textcircled{1}}{\textcircled{7} \cdot 31536000}$$

- 9 The sum of activity times mass and mass squared

- 10 The average half-life

- 11 Half-life calculated by means of a linear regression

$$t_{1/2}^{calc.}(LR) = \frac{0.693}{31536000} \times \frac{\sum N^2}{\sum A \cdot N}$$

- 12 Experimental half-life of the isotope in years (the one obtained from the tables)

- 13 The percent error

STUDENT INFO

Name:

Date:

Pre-lab Questions**Radioactivity: half-life of K-40**

1. Research the half-life of the following isotopes: (a) Chlorine-36 (b) Cadmium-109 (c) Copper-64 (d) Chromium-51 (e) Carbon-14 (f) Gold-195 (g) Calcium-45

2. Classify the following nuclear reactions as: (a) α decay (b) β decay (c) γ decay (d) positron emission (e) electron capture
(i) ${}^{14}_6\text{C} \longrightarrow {}^{14}_7\text{N} + {}^0_{-1}\beta$ (ii) ${}^{11}_6\text{C} \longrightarrow {}^{11}_5\text{B} + {}^0_{+1}\beta^+$ (iii) ${}^{55}_{26}\text{Fe} + {}^0_{-1}\beta \longrightarrow {}^{55}_{25}\text{Mn} + \text{X-ray}$ (iv) ${}^{234}_{88}\text{Th}^* \longrightarrow {}^{234}_{88}\text{Th} + {}^0_0\gamma$ (v) ${}^{226}_{88}\text{Ra} \longrightarrow {}^{222}_{86}\text{Rn} + {}^4_2\alpha$

3. The half-life of bromine-74 is 25 min. How much of a 100 mg sample is still active after 100 min?

4. Identify the unknown radioactive particle involved in the following nuclear equations: (a) ${}^9_4\text{Be} + {}^A_Z\text{X} \longrightarrow {}^{12}_6\text{C} + {}^1_0\text{n}$ (b) ${}^{31}_{15}\text{P} + {}^1_1\text{H} \longrightarrow {}^{31}_{16}\text{S} + {}^A_Z\text{X}$ (c) ${}^3_1\text{H} + {}^2_1\text{H} \longrightarrow {}^A_Z\text{X} + {}^1_0\text{n}$ (d) ${}^{14}_6\text{C} \longrightarrow {}^A_Z\text{X} + {}^0_{-1}\beta$

5. Indicate the nuclear symbol for (a) Oxygen-18 (b) Magnesium-24 (c) Lithium-7

Results

Name:

Date:

Radioactivity: half-life of K-40

0 m (K) in g=_____

1 N (⁴⁰K) in atoms=_____

[illegible]

0 m (K) in g=_____

1 $N(^{40}\text{K})$ in atoms=_____

[illegible]

$\textcircled{0} \quad m \text{ (K) in g} = \underline{\hspace{2cm}}$
 $\textcircled{1} \quad N \text{ (}^{40}\text{K) in atoms} = \underline{\hspace{2cm}}$

Activity, A (cpm)	1	2	3	4	5	6	7	8	9	10	\overline{A} (cpm)	\overline{A} (cps)	s (cps)
Sample + Background $\textcircled{6}$													
Sample $\textcircled{7}$													

$\textcircled{0} \quad m \text{ (K) in g} = \underline{\hspace{2cm}}$
 $\textcircled{1} \quad N \text{ (}^{40}\text{K) in atoms} = \underline{\hspace{2cm}}$

Activity, A (cpm)	1	2	3	4	5	6	7	8	9	10	\overline{A} (cpm)	\overline{A} (cps)	s (cps)
Sample + Background $\textcircled{6}$													
Sample $\textcircled{7}$													

N (atoms)	$t_{1/2}$ (s)	A (cpm)	$A \cdot N$	N^2
$\textcircled{1}$	$\textcircled{8}$	$\textcircled{2}$	$\textcircled{1} \times \textcircled{2}$	$\textcircled{1} \times \textcircled{1}$
$\textcircled{9} \quad \text{Sum } (\Sigma) =$				

$\textcircled{10} \quad \overline{t_{1/2}}^{calc.} \text{ (years)} = \underline{\hspace{2cm}}$
 $\textcircled{12} \quad t_{1/2}^{theory} \text{ (years)} = \underline{\hspace{2cm}}$
 $\textcircled{13} \quad \% \text{ error} = \underline{\hspace{2cm}}$

$\textcircled{11} \quad t_{1/2}^{calc.} (LR) \text{ (years)} = \underline{\hspace{2cm}}$
 $\textcircled{12} \quad t_{1/2}^{theory} \text{ (years)} = \underline{\hspace{2cm}}$
 $\textcircled{13} \quad \% \text{ error} = \underline{\hspace{2cm}}$




Figure 12: \bar{A} cps (Y axis) vs. N (X axis)

STUDENT INFO

Name:

Date:

Post-lab Questions**Radioactivity: half-life of K-40**

-
- Classify the following nuclear reactions as: (a) α decay (b) β decay (c) γ decay (d) positron emission (e) electron capture
(i) ${}_{92}^{238}\text{U} \longrightarrow {}_{90}^{234}\text{Th} + {}_2^4\text{He}$ (ii) ${}_{19}^{42}\text{K} \longrightarrow {}_{20}^{42}\text{Ca} + {}_{-1}^0\text{e}$ (iii) ${}_{8}^{15}\text{O} \longrightarrow {}_{7}^{15}\text{N} + {}_{+1}^0\text{e}$ (iv) ${}_{88}^{228}\text{Ra} \longrightarrow {}_{89}^{228}\text{Ac} + {}_{-1}^0\text{e}$ (v) ${}_{6}^{13}\text{C} + {}_{+1}^1\text{H} \longrightarrow {}_{7}^{14}\text{N} + {}_{0}^0\gamma$
 - Indicate the name of the following nuclear symbols: (a) ${}_{+1}^0\text{e}^{+}$ (b) ${}_{0}^0\gamma$ (c) ${}_{1}^1\text{H}$
 - Research the half-life of the following isotopes: (a) Potassium-40 (b) Cesium-137 (c) Cobalt-57 (d) Bismuth-212 (e) Gallium-67 (f) Americium-241
 - The half-life of bromine-74 is 25 min. 20mg of the isotopes remain after 10 minutes of preparing the sample. Calculate the initial mass of the bromine-74 sample.

Periodic Table of Chemical Elements

18 VIIIA

1 IA

1	1.0079	2	4.0025
1	H Hydrogen	13 IIIA	He Helium
2	3	14 IVA	17 VIIA
2	4	15 VA	18
3	5	16 VIA	19
3	6	17 VIIA	20
3	7	18	21
4	8	19	22
4	9	20	23
5	10	21	24
5	11	22	25
6	12	23	26
6	13	24	27
6	14	25	28
7	15	26	29
7	16	27	30
7	17	28	31
7	18	29	32
7	19	30	33
7	20	31	34
7	21	32	35
7	22	33	36
7	23	34	37
7	24	35	38
7	25	36	39
7	26	37	40
7	27	38	41
7	28	39	42
7	29	40	43
7	30	41	44
7	31	42	45
7	32	43	46
7	33	44	47
7	34	45	48
7	35	46	49
7	36	47	50
7	37	48	51
7	38	49	52
7	39	50	53
7	40	51	54
7	41	52	55
7	42	53	56
7	43	54	57
7	44	55	58
7	45	56	59
7	46	57	60
7	47	58	61
7	48	59	62
7	49	60	63
7	50	61	64
7	51	62	65
7	52	63	66
7	53	64	67
7	54	65	68
7	55	66	69
7	56	67	70
7	57	68	71
7	58	69	72
7	59	70	73
7	60	71	74
7	61	72	75
7	62	73	76
7	63	74	77
7	64	75	78
7	65	76	79
7	66	77	80
7	67	78	81
7	68	79	82
7	69	80	83
7	70	81	84
7	71	82	85
7	72	83	86
7	73	84	87
7	74	85	88
7	75	86	89
7	76	87	90
7	77	88	91
7	78	89	92
7	79	90	93
7	80	91	94
7	81	92	95
7	82	93	96
7	83	94	97
7	84	95	98
7	85	96	99
7	86	97	100
7	87	98	101
7	88	99	102
7	89	100	103
7	90	101	104
7	91	102	105
7	92	103	106
7	93	104	107
7	94	105	108
7	95	106	109
7	96	107	110
7	97	108	111
7	98	109	112
7	99	110	113
7	100	111	114
7	101	112	115
7	102	113	116
7	103	114	117
7	104	115	118
7	105	116	119
7	106	117	120
7	107	118	121
7	108	119	122
7	109	120	123
7	110	121	124
7	111	122	125
7	112	123	126
7	113	124	127
7	114	125	128
7	115	126	129
7	116	127	130
7	117	128	131
7	118	129	132
7	119	130	133
7	120	131	134
7	121	132	135
7	122	133	136
7	123	134	137
7	124	135	138
7	125	136	139
7	126	137	140
7	127	138	141
7	128	139	142
7	129	140	143
7	130	141	144
7	131	142	145
7	132	143	146
7	133	144	147
7	134	145	148
7	135	146	149
7	136	147	150
7	137	148	151
7	138	149	152
7	139	150	153
7	140	151	154
7	141	152	155
7	142	153	156
7	143	154	157
7	144	155	158
7	145	156	159
7	146	157	160
7	147	158	161
7	148	159	162
7	149	160	163
7	150	161	164
7	151	162	165
7	152	163	166
7	153	164	167
7	154	165	168
7	155	166	169
7	156	167	170
7	157	168	171
7	158	169	172
7	159	170	173
7	160	171	174
7	161	172	175
7	162	173	176
7	163	174	177
7	164	175	178
7	165	176	179
7	166	177	180
7	167	178	181
7	168	179	182
7	169	180	183
7	170	181	184
7	171	182	185
7	172	183	186
7	173	184	187
7	174	185	188
7	175	186	189
7	176	187	190
7	177	188	191
7	178	189	192
7	179	190	193
7	180	191	194
7	181	192	195
7	182	193	196
7	183	194	197
7	184	195	198
7	185	196	199
7	186	197	200
7	187	198	201
7	188	199	202
7	189	200	203
7	190	201	204
7	191	202	205
7	192	203	206
7	193	204	207
7	194	205	208
7	195	206	209
7	196	207	210
7	197	208	211
7	198	209	212
7	199	210	213
7	200	211	214
7	201	212	215
7	202	213	216
7	203	214	217
7	204	215	218
7	205	216	219
7	206	217	220
7	207	218	221
7	208	219	222
7	209	220	223
7	210	221	224
7	211	222	225
7	212	223	226
7	213	224	227
7	214	225	228
7	215	226	229
7	216	227	230
7	217	228	231
7	218	229	232
7	219	230	233
7	220	231	234
7	221	232	235
7	222	233	236
7	223	234	237
7	224	235	238
7	225	236	239
7	226	237	240
7	227	238	241
7	228	239	242
7	229	240	243
7	230	241	244
7	231	242	245
7	232	243	246
7	233	244	247
7	234	245	248
7	235	246	249
7	236	247	250
7	237	248	251
7	238	249	252
7	239	250	253
7	240	251	254
7	241	252	255
7	242	253	256
7	243	254	257
7	244	255	258
7	245	256	259
7	246	257	260
7	247	258	261
7	248	259	262
7	249	260	263
7	250	261	264
7	251	262	265
7	252	263	266
7	253	264	267
7	254	265	268
7	255	266	269
7	256	267	270
7	257	268	271
7	258	269	272
7	259	270	273
7	260	271	274
7	261	272	275
7	262	273	276
7	263	274	277
7	264	275	278
7	265	276	279
7	266	277	280
7	267	278	281
7	268	279	282
7	269	280	283
7	270	281	284
7	271	282	285
7	272	283	286
7	273	284	287
7	274	285	288
7	275	286	289
7	276	287	290
7	277	288	291
7	278	289	292
7	279	290	293
7	280	291	294
7	281	292	295
7	282	293	296
7	283	294	297
7	284	295	298
7	285	296	299
7	286	297	300
7	287	298	301
7	288	299	302
7	289	300	303
7	290	301	304
7	291	302	305
7	292	303	306
7	293	304	307
7	294	305	308
7	295	306	309
7	296	307	310
7	297	308	311
7	298	309	312
7	299	310	313
7	300	311	314
7	301	312	315
7	302	313	316
7	303	314	317
7	304	315	318
7	305	316	319
7	306	317	320
7	307	318	321
7	308	319	322
7	309	320	323
7	310	321	324
7	311	322	325
7	312	323	326
7	313	324	327
7	314	325	328
7	315	326	329
7	316	327	330
7	317	328	331
7	318	329	332
7	319	330	333
7	320	331	334
7	321	332	335
7	322	333	336
7	323	334	337
7	324	335	338
7	325	336	339
7	326	337	340
7	327	338	341
7	328	339	342
7	329	340	343
7	330	341	344
7	331	342	345
7	332	343	346
7	333	344	347
7	334	345	348
7	335	346	349
7	336	347	350
7	337	348	351
7	338	349	352
7	339	350	353
7	340	351	354
7	341	352	355
7	342	353	356
7	343	354	357
7	344	355	358
7	345	356	359
7	346	357	360
7	347	358	361
7	348	359	362
7	349	360	363
7	350	361	364
7	351	362	365
7	352	363	366
7	353	364	367
7	354	365	368
7	355	366	369
7	356	367	370
7	357	368	371
7	358	369	372
7	359	370	373
7	360	371	374
7	361	372	375
7	362	373	376
7	363	374	377
7	364	375	378
7	365	376	379
7	366	377	380
7	367	378	381
7	368	379	

